



May 24, 2023

## FIP200-eGFP expression and purification

DOI

[dx.doi.org/10.17504/protocols.io.dm6gpbkq5lzp/v1](https://dx.doi.org/10.17504/protocols.io.dm6gpbkq5lzp/v1)

Eleonora Turco<sup>1</sup>, Justyna Sawa-Makarska<sup>1</sup>

<sup>1</sup>Max Perutz Labs, University of Vienna



Dorotea Fracchiolla

Sascha Martens lab, University of Vienna, Max Perutz Labs - ...

### Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.dm6gpbkq5lzp/v1>

**Protocol Citation:** Eleonora Turco, Justyna Sawa-Makarska 2023. FIP200-eGFP expression and purification. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.dm6gpbkq5lzp/v1>

**License:** This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working

We use this protocol and it's working

**Created:** July 13, 2022



**Last Modified:** May 31, 2024

**Protocol Integer ID:** 66613

**Keywords:** FIP200 purification, FIP200-eGFP purification, ASAPCRN, egfp expression, human fip200, egfp, purification this protocol, purification

## Abstract

This protocol describes how to express and purify human FIP200 tagged C-terminally with eGFP.

## Attachments



protocols.io FIP200....

192KB



## Materials

### Expression:

pGB-GST-3C-FIP200-GFP (Addgene ID: 187832)

Sf9 insect cells

SF921 medium with antibiotics 100 IU/ml Penicillin and 100 µg/ml Streptomycin

### Lysis Buffer:

50 mM HEPES pH 7.5

300 mM NaCl

1 mM MgCl<sub>2</sub>

10% glycerol

0.5% CHAPS

5 U/ml Benzonase (Sigma)

1 mM DTT

CIP protease inhibitor (Sigma)

cOmplete EDTA-free protease inhibitor cocktail (Roche)

### Wash Buffer:

50 mM HEPES pH 7.5

200 mM NaCl

1 mM MgCl<sub>2</sub>

1 mM DTT

### Gel-filtration buffer:

25 mM HEPES pH 7.5

200 mM NaCl

1 mM DTT

### Columns/Resin:

Glutathione Sepharose 4B (Cytiva)

Superose 6 increase 10/300 column (Cytiva)

## Troubleshooting

## Expression

- 1 To generate FIP200-GFP constructs the insect codon optimized FIP200 gene was purchased from GenScript and cloned with respective tags into pGB-02-03 (pGB-GST-3C-FIP200-GFP - Addgene ID: 187830). Generated construct was used for expression in Sf9 insect cells using the Bac-to-Bac system (ThermoFischer Scientific) .
- 2 Transfect 2.5 µg of bacmid DNA into Sf9 insect cells in 6-well plate using FuGene transfection reagent (Promega).
- 3 About 7 days after transfection the V0 virus should be ready for harvesting. Use the V0 to produce a V1 virus stock by infecting 30 ml of Sf9 cells (1 million/ml). Collect V1 about 4-5 days later. Monitor viability of the cells and green fluorescence to decide when to collect V1.
- 4 Infect 1L culture of Sf9 cells at 1-1.5 million/ml cells/volume at 99-100% viability in log phase with 1 ml of Virus 1 (V1).
- 5 After infection monitor cells for viability and fluorescence. Harvest by centrifugation when the viability drops to 80–95% and clear green fluorescence is present.
- 6 To harvest spin down the cells at 2000 rpm, for 15 min at RT (Sorvall RC6+ centrifuge, Thermo Scientific). Gently wash the cell pellets with PBS, flash-freeze in liquid nitrogen, and store at –80 °C until purification.

## Purification

- 7 Thaw a cell pellet corresponding to 1L culture by re-suspending it in 40 ml lysis buffer (50 mM HEPES pH 7.5, 300 mM NaCl, 1 mM MgCl<sub>2</sub>, 10% glycerol, 0.5% CHAPS, 5 U/ml Benzonase (Sigma), 1 mM DTT, CIP protease inhibitor (Sigma), cOmplete EDTA-free protease inhibitor cocktail (Roche)) and rolling or stirring in the cold room.
- 8 Additionally disrupt the cells with a Dounce homogenizer followed by 1 min sonication at 50% cycles and 50–60% power.
- 9 Clear the lysate by centrifugation at 72,000 × *g* for 45 min at 4 °C (Beckman Ti45 rotor)
- 10 Incubate the cleared supernatant with 5 ml of Glutathione Sepharose 4B beads slurry (Cytiva) overnight at 4°C rolling gently. The GSH slurry should be washed with water and then pre-equilibrated with lysis buffer before incubating with the lysate.



- 11 Wash the beads seven times with wash buffer (50 mM HEPES pH 7.5, 200 mM NaCl, 1 mM  $\text{MgCl}_2$ , 1 mM DTT).
- 12 Incubate the beads overnight with preCission 3C protease at 4°C (40  $\mu\text{l}$  of 6 mg/ml home-made protease).
- 13 Spin down the beads (4000 rpm, 3 min, 4°C) and collect the supernatant containing cleaved FIP200-eGFP.
- 14 Filter the supernatant through a 0.45  $\mu\text{m}$  syringe filter to remove any residual beads.
- 15 Concentrate the protein down to 0.5 ml using a 100 kDa cut-off Amicon filter and apply onto a Superose 6 Increase 10/300 column (Cytiva) pre-equilibrated with a SEC buffer containing 25 mM HEPES pH 7.5, 200 mM NaCl, 1 mM DTT. Pool fractions containing pure proteins (see attached pdf), concentrate, snap freeze in liquid nitrogen, and store at -80°C.