


Oct 16, 2019

# First strand cDNA synthesis (ThermoScientific RevertAid)

 Forked from [First strand cDNA synthesis \(ThermoScientific RevertAid\)](#)

DOI

[dx.doi.org/10.17504/protocols.io.8bmhsk6](https://dx.doi.org/10.17504/protocols.io.8bmhsk6)

Ben Kuipers<sup>1</sup>

<sup>1</sup>Wageningen University

iGEM Wageningen 2019



Ben Kuipers

OPEN  ACCESS



DOI: [dx.doi.org/10.17504/protocols.io.8bmhsk6](https://dx.doi.org/10.17504/protocols.io.8bmhsk6)

**Protocol Citation:** Ben Kuipers 2019. First strand cDNA synthesis (ThermoScientific RevertAid). **protocols.io**  
<https://dx.doi.org/10.17504/protocols.io.8bmhsk6>

**License:** This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working

**We use this protocol and it's working**

**Created:** October 16, 2019

**Last Modified:** October 16, 2019

**Protocol Integer ID:** 28749

## Abstract


The following protocol is optimized to generate first-strand cDNA for use in (q)PCR


## Materials

### MATERIALS

 5X RT Buffer **Thermo Fisher Scientific Catalog ##B91**

 dNTP Mix 10 mM each **Thermo Fisher Scientific Catalog ##R0191**

 Water, nuclease free

 RiboLock RNase Inhibitor **Thermo Fisher Scientific Catalog ##EO0381**

 RevertAid Reverse Transcriptase (200 U/μL) **Thermo Fisher Catalog #EP0442**

 Oligo(dT)18 Primer **Thermo Fisher Catalog #SO131**

### Before start

Mix and briefly centrifuge all reagents after thawing, keep on ice.

- 1 Add reaction components into sterile, nuclease-free tube on ice in the indicated order:

Template RNA	100 ng ( 1pg - 5 µg)
Oligo(dT)18	1 µl (100 pmol)
<b>Water, nuclease-free</b>	<b>to 12 µl</b>

- 2 **Optional:** If the RNA template is GC-rich or is known to contain secondary structures, mix gently, centrifuge briefly and incubate at 65 °C for 5min. Chill on ice, briefly centrifuge again and place on ice.

- 3

5X RT Buffer	4 µl
RiboLock RNase Inhibitor	1 µl (20 U)
RevertAid RT (200 U/µL)	2 µl (400 U)
10 mM dNTP Mix	1 µl
<b>Total volume</b>	<b>20 µl</b>

Mix gently and centrifuge briefly.

- 4

5 min	25 °C
60 min	42 °C (For GC-rich RNA, the reaction temperature can be increased to 45 °C)
5 min	70 °C

- 5 The cDNA product is now ready for downstream applications