

May 20, 2025

## FISH CODETECTION

DOI

<https://dx.doi.org/10.17504/protocols.io.bp2l6y631vqe/v1>

Sylvie Dumas<sup>1</sup>, Åsa allén-Mackenzie<sup>1</sup>

<sup>1</sup>Lund University



amanda schneeweis

Northwestern

### Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.bp2l6y631vqe/v1>

**Protocol Citation:** Sylvie Dumas, Åsa allén-Mackenzie 2025. FISH CODETECTION. **protocols.io**  
<https://dx.doi.org/10.17504/protocols.io.bp2l6y631vqe/v1>

**License:** This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working

We use this protocol and it's working

**Created:** May 19, 2025

Last Modified: May 20, 2025

Protocol Integer ID: 218572

**Keywords:** ASAPCRN, riboprobes fluorescent in situ, riboprobes fluorescent, fluorescein, fish codetection dig, brain section

**Funders Acknowledgements:**


ASAP

Grant ID: ASAP-025188

## Abstract

DIG and Fluorescein labeled  
riboprobes  
Fluorescent *in situ*  
hybridization  
Frozen brain sections

## Image Attribution



## Troubleshooting

## PCR production by plasmid of interest

### 1 PCR reagents/plasmid

- In ice, place in a PCR tube :

H<sub>2</sub>O: 21µl

Mix primers S et AS (10 µM each) : 3 µl

DNA plasmid (10ng/µl): 6 µl

REDEXTRACT PCR Mix : 30µl

Dispatch 2 times 19,8uL of PCR mix in two other PCR tubes

(Note that the optimal vol for this PCR is 20 µl).

### 2 PCR Cycles

	94°C	5 min
	94°C	45sec
x 35	55°C	45 sec
	72°C	1 min
	72°C	5 min
	4°C	...

### 3 PCR gel control:

Make a 1% agarose gel

-1 µl 1kb+DNA ladder+ 1µl of 10X Blue juice Gel loading buffer + 8 µl H<sub>2</sub>O

-5 µl of PCR product/ mRNA

### 4 PCR purification, (QIAquick PCR Purification kit)

o Pool the 3 PCR tubes in an 1,5 ml tube

o Add 5 volume PB buffer to 1 volume of the PCR reaction (roughly 250 µl de PB) and mix

o Place a QIAquick column in a provided 2 ml collection tube

o To bind DNA, apply the sample to the QIAquick column.

o Centrifuge

1 min at 13 000 rpm

o Discard flow-through

o Place the QIAquick column back in the same tube

o To wash, add 750 ul of PE Buffer to the QIAquick column



- o Centrifuge 1 min à 13 000 rpm .
- o Discard flow-through
- o Place the QIAquick column back in the same tube
- o Centrifuge 1 min à 13 000 rpm to remove residual wash buffer.
- o Discard flow-through
- o Place the QIAquick column in a clean 1,5 ml RNase free tube
- o To elute DNA, add 50 µl Buffer EB to the center of the QIAquick membrane and wait for 3 min.
- o Centrifuge 1 min à 13 000 rpm.

## 5 Dosage : Nano drop

## Riboprobe synthesis

### 6 Transcription reagents/riboprobe

In ice, place in a 1,5 ml tube :

-		
5X transcription buffer		4µl
10X Dig-RNA labelling mix		
<b>or</b> 10X Fluorescein-RNA labelling mix		2µl
100mM DTT		1µl
Rnasin ( 20U/µl)		0.5µl
T3/T7 RNA		
polymerase (20U/µl)	2.0µl	
PCR-product		150 ng
H2O		up to 20 µl

### 7 Transcription reaction

Incubate 2h at **37°C**

### 8 Riboprobe purification (illustra ProbeQuant G-50 Micro Columns)

#### Column preparation

- Re-suspend the resin in the column by vortexing
- Loosen the cap one-quarter turn and twist off the bottom closure



- Place the G-50 column in the supplied Collection tube
- Centrifuge 1 minute 735 × g
- Place the G-50 column in a clean 1,5 ml RNase free tube

#### Sample application

- Add to the riboprobe sample 30 µl Probe Buffer type 1.
- apply the sample in the middle of the column.
- 

Centrifuge 2 min à 750 g at 4°C

-

Place tube in ice

#### 9 Dosage (Nanodrop).

#### 10 Storage

- Add 100 µl Hyb Buffer (50% « Hybe-F » + 50% Formamide).
- Aliquot in 100 µl tube (50 µl/tube)
- Keep at -80°C.

## Tissue Preparation

- 11 Adult mice are sacrificed by cervical dislocation, and brains are rapidly dissected out and snap frozen in cold (-30°C to -35°C) 2-methylbutane (≥99%, Honeywell, M32631) on dry ice. Brains are sealed with parafilm and aluminum foil and stored at -80°C until use. Embryonic brains, which are less rich in myelin and fat, require more gentle freezing in 2-methylbutane (-20°C to -25°C) immediately after extraction from the euthanized and sacrificed mother.

## Hybridization

- 12
- 1- Thaw slides (around 5 min)
  - 2- Fix slides in 4% PFA in PBS, 10 min at RT.
  - 3- Wash 3×10 min in PBS
  - 4- Acetylation treatment:
    - Equilibrate slides in TEA 5 min
    - Prepare extemporaneously the acetylated buffer in the following ratio:

acetic anhydride 25 µl in 10 ml TEA

- Acetylation treatment, 10 min at RT
- Wash 3×5 min in PBS

-5- Hybridization:

- Make the Hyb Buffer by mixing Hyb Buffer

without Formamide + Formamide (50/50). Keep it on ice.

- For each riboprobe. Add to 20 µl Hyb Buffer the 2 riboprobes. Amount in accordance with the number of slides that will be hybridized:

- o Dig-probe (50-75ng/100µl)
- o Fluoresceinprobe (75-100 ng/100 µl)

- Denature 10 min at 85 °C
- Put in ice immediately 1 min
- Briefly centrifuge tubes at 2000g
- Put in ice
- Add XXX µl of cold Hyb Buffer adapted to the number of slides that will be hybridized.

Keep tubes in ice

- Vortex
- Take slide one by one from the PBS, remove the PBS by shaking the slide, dry the bottom of the slide with a klenex. Place slide in 3 slides hybridization box.
- Put 150 µl of (Hyb Buffer+ denaturated riboprobe) /slide
- Cover with a RNase free coverslip (Hybrislip)
- Place the 3 slides hybridization box in a humidified chamber (5xSSC)

Note: no need to put Formamide here!

- Hybridize at 65°C between 16h and 18h.
- Prewarm X5 SSC and 0,2X SSC wash solutions at 65°C

## Day 2: washes and immunological staining

- 13
- Place slide one by one in prewarmed X5 SSC and remove the coverslip
  - Transfer slides in prewarmed X5 SSC and wash 5 min at 65°C
  - Wash 3×20 min in prewarmed 0.2X SSC at 65°C
  - Wash 10 min in 0.2X SSC at RT
  - Wash 3 × 10min in MABT (RT)

14 **From now, all steps at RT**

**Fluorescein-riboprobe revelation.**

- Place slides in humidified chamber (H<sub>2</sub>O)
- Add 500 µl of BR/slide for 20 min
- Replace solution by 1/5000 anti-Fluorescein-POD in BR (500µl/lame) for 1h30
- Wash 3× 10 min in MABT



- Wash 3 × 10 min PBST

During the last wash, prepare a H<sub>2</sub>O<sub>2</sub> 0,001% solution from a 30% H<sub>2</sub>O<sub>2</sub> by cascade dilution:

- 10µL H<sub>2</sub>O<sub>2</sub> 30% + 90µL PBST ◊ H<sub>2</sub>O<sub>2</sub> 3%
  - 10µL H<sub>2</sub>O<sub>2</sub> 3% + 590µL PBST ◊ H<sub>2</sub>O<sub>2</sub> 0.05%
  - 10µL H<sub>2</sub>O<sub>2</sub> 0.05% + 490µL PBST ◊ H<sub>2</sub>O<sub>2</sub> 0.001%
- From now everything must be done away from light.
  - Place slides in humidified chamber (H<sub>2</sub>O)
  - Place the solution 1/250 Fluorescein-tyramide in H<sub>2</sub>O<sub>2</sub> 0.001% (500 µl/slide)
  - When all slides are covered, gently move the slide for a uniform distribution of the solution at the slide surface
  - Incubation 1h at RT
  - Wash 3 × 10 min PBST

## 15 **POD- inhibition:**

- 
- Wash 20 min in Glycine buffer 0,1M pH= 2.1
- Quick wash in PBS
- Wash 20 min in 3% H<sub>2</sub>O<sub>2</sub> in PBS
- Wash 3 ×10 min in PBST
- Keep the slides in the last wash at 4°C overnight

## Day3: washes and immunological staining

### 16 **Dig-riboprobe revelation:**

- Wash 1 ×10 min in PBST
  - Wash 3 × 10 min MABT
  - Place slides in humidified chamber (H<sub>2</sub>O)
  - Add 450 µl of BR/slide for 20 min
  - Replace solution by 1/1500 anti-DIG-POD in BR (450µl/lame)
- for 1h
- Wash 3× 10 min in MABT
  - Wash 3 × 10 min PBST

During the last wash, prepare a H<sub>2</sub>O<sub>2</sub> 0,001% solution

- Prepare 0.001% H<sub>2</sub>O<sub>2</sub> from 30% H<sub>2</sub>O<sub>2</sub> by cascade dilution:

- 10µL H<sub>2</sub>O<sub>2</sub> 30% + 90µL PBST ◊ H<sub>2</sub>O<sub>2</sub> 3%
  - 10µL H<sub>2</sub>O<sub>2</sub> 3% + 590µL PBST ◊ H<sub>2</sub>O<sub>2</sub> 0.05%
  - 10µL H<sub>2</sub>O<sub>2</sub> 0.05% + 490µL PBST ◊ H<sub>2</sub>O<sub>2</sub> 0.001%
  
  - Place slides in humidified chamber (H<sub>2</sub>O)
  - Place the solution 1/250 Tamra-tyramide in H<sub>2</sub>O<sub>2</sub> 0.001% (500 µl/slide)
  - When all slides are covered, gently move the slide for a uniform distribution of the solution at the slide surface
  -
- Incubation 1h at RT
- Wash 5 × 10 min PBST
  - Place slides in humidified chamber (H<sub>2</sub>O)
  - Add 450 µl/slide of DAPI for 20 min
  - Wash 3× 5min PBST
  - Mount in fluoromount.
  - Quick microscope slide analysis
- Give slide to be scanned