

Feb 14, 2025

Fixation & Tissue Processing Protocol

DOI

<https://dx.doi.org/10.17504/protocols.io.bp2l6dzy1vqe/v1>

Roy Ambli Dalmo¹

¹The Arctic University of Norway



Celine Richard

The Arctic University of Norway

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.bp2l6dzy1vqe/v1>

Protocol Citation: Roy Ambli Dalmo 2025. Fixation & Tissue Processing Protocol. protocols.io
<https://dx.doi.org/10.17504/protocols.io.bp2l6dzy1vqe/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: February 13, 2025



Last Modified: February 14, 2025

Protocol Integer ID: 120172

Keywords: tissue processing protocol fixation, fixation, processing, tissue, protocol

Abstract

Fixation & Tissue Processing Protocol



Fixative

1

Note

- * Freshly made RNase-free 4% paraformaldehyde (4% PF) used specifically for *In Situ* Hybridization
- * 4% paraformaldehyde or 10% buffered formalin used for routine histology (i.e., H&E staining)
- * 2% paraformaldehyde (bgal, GFP, RFP), 4% paraformaldehyde or 10% buffered formalin used for Immunocytochemistry

Fixation procedure

2

Note

The following protocols are for tissue sample processing for subsequent light microscopy with H&E staining, immunocytochemistry or in situ hybridization. RNase-free PBS and fixative should be used for in situ hybridization with RNA probes.

* Use the 20 ml scintillation vials with 15-20 ml solution volumes or 60 ml yellow-capped containers with 30-40 ml solution volumes to fix, wash and dehydrate the samples.

Well-fixed samples can be stored in the same fixative at 4°C for several days or even longer for H&E staining. In situ hybridization samples must be dehydrated and stored at -20 °C with 100% ethanol after 24-48 hours fixation. Immunohistochemistry samples must be dehydrated and stored at -20 °C with 100% ethanol after overnight fixation.

Tissue fixation

- 3 Isolate the tissue into **cold PBS** as soon as possible (can be omitted).
- 4 Wash tissue with **PBS** to remove all blood (can be omitted).
- 5 Place tissue in **fixative** for 10-15 minutes to one hour (can be omitted)





- 6 Cut tissue to proper size. The size can be 2X2 mm to 1X2 cm but thickness should be 3 mm for better fixation. The cutting surface of the tissue should be flat and smooth.
- 7 Transfer tissue to **fixative** and swirl the container to ensure all tissues are completely immersed in fixative. The volume of fixative must be 20-30 times the tissue volume.
- 8 Fix tissue at 4°C for overnight.
- 9 Check sample the next day to ensure proper fixation.



For best morphology and staining results take care to:

10

Note

1. Use high-quality fixative and sufficient amounts of fixative.
2. Use sharp razor blade and fine instruments such as scissors and forceps to handle all specimens. A piece of dental wax can be used as a cutting board. Treat tissues as gently as possible to minimize morphology artifacts.
3. Don't let samples dry during preparation. Samples should always be covered by solutions.