

Jul 06, 2023

FixNCut v1.0

DOI

dx.doi.org/10.17504/protocols.io.14egn3xjql5d/v1

FIX

Laura Jiménez-Gracia¹, Domenica Marchese¹, Holger Heyn¹, Luciano G Martelotto²

¹Centre for Genomic Regulation; ²University of Adelaide

Human Cell Atlas Metho...

Single Cell Core, Harvar...

1 more workspace



Luciano G Martelotto

University of Adelaide

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.14egn3xjql5d/v1>

Protocol Citation: Laura Jiménez-Gracia, Domenica Marchese, Holger Heyn, Luciano G Martelotto 2023. FixNCut v1.0. protocols.io <https://dx.doi.org/10.17504/protocols.io.14egn3xjql5d/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: June 22, 2023

Last Modified: July 06, 2023

Protocol Integer ID: 83875

Keywords: Fixation for cells, Fixation for tissues, protocol details reversible fixation for cell, reversible fixation, fixncut, omics compatibility, fixation, multiome, snrna, snatac, tissues for subsequent use

Disclaimer

DISCLAIMER – FOR INFORMATIONAL PURPOSES ONLY; USE AT YOUR OWN RISK

The protocol content here is for informational purposes only and does not constitute legal, medical, clinical, or safety advice, or otherwise; content added to **protocols.io** is not peer reviewed and may not have undergone a formal approval of any kind. Information presented in this protocol should not substitute for independent professional judgment, advice, diagnosis, or treatment. Any action you take or refrain from taking using or relying upon the information presented here is strictly at your own risk. You agree that neither the Company nor any of the authors, contributors, administrators, or anyone else associated with **protocols.io**, can be held responsible for your use of the information contained in or linked to this protocol or any of our Sites/Apps and Services.

Abstract

This protocol details reversible fixation for cells and tissues for subsequent use in sc/snRNA, sc/snATAC or Multiome. Spatial-Omics compatibility is being validated. For more information check this preprint:

<https://www.biorxiv.org/content/10.1101/2023.06.16.545221v2>

Attachments



[748-1909.docx](#)

16KB



Guidelines




Notes:

All washes and centrifugations need to be done at  4 °C unless otherwise specified.

IMPORTANT: Washing volumes may change accordingly to your needs. If you want to change the protocol, let's discuss just in case. Time, temperatures and concentrations must be maintained.

DSP has been used before for fixing cells and prep RNA, it works fine. For single cell or tissue following dissociation is what we have been studying and works well. It's still work in progress, but the key is to keep the stock fixative away from water because it neutralizes the NHS-esters quickly.

Points to take into account:

- **CRITICAL:** Prepare working solution (1x) right before fixation (no more than 5 minutes). For larger pieces replace with fresh 1x fixative a couple of times.
- make single use aliquots (20-50 uL) for 2 to 5 fixations. What's not used do not re-freeze (it is fine to re freeze let's not give them the option)
- keep aliquots at  -80 °C in a bag (with silica if possible).
- bring tubes at  Room temperature and prepare the fixation a few minutes (no longer than 10 min) before fixing. This will ensure that the NHS-ester isn't in contact with aqueous solution for too long.
- evaluate small precipitation during fixation. Too much ppt: bad. You should see a small ppt on the walls of the tube, like in the attached photo. You will notice that the first 2 drops of PBS will be generate the precipitate but will precipitate as you add more.
- do not prepare aliquots larger than  500 µL .
- do not store 1x solutions.
- We have noticed some performance variability from vial to vial purchased from Sigma.
- Viability is not a good measure because PI or Trypan Blue don't enter after fixation and the cells look alive.
- So far, the best test has been the small shift in the LMO-FAM or LMO-Cy5 on cells.

Materials

 Liberase™ TM Research Grade **Merck MilliporeSigma (Sigma-Aldrich) Catalog #5401127001** ,

 Nuclei Isolation Kit: Nuclei EZ Prep **Merck MilliporeSigma (Sigma-Aldrich) Catalog #NUC101-1KT**

Troubleshooting




Preparation of DSP (Oz Soup) stock and working solutions




1



30m

Note






DSP (dithiobis(succinimidyl propionate)) also known as Lomant's Reagent and can be purchased from Thermo: <https://www.thermofisher.com/order/catalog/product/22585>.

 DSP (dithiobis(succinimidyl propionate)) Lomants Reagent **Thermo Fisher Scientific Catalog #22585**

Equilibrate DSP vial at  Room temperature for  00:30:00 and then prepare 50x stock solution of DSP ( 50 mg/mL) in molecular biology grade dimethyl sulfoxide (Sigma, cat. no. D8418-50ML).

2 Dispense the stock into  100 μ L aliquots and store at  -80 °C .






3 Immediately before use prepare  500 μ L of  1 mg/mL DSP working solution (DSP 1x is also known as OzSoup) in molecular biology grade 1x PBS as follows: aliquot  10 μ L of stock DSP in a 1.5mL eppi tube and while vortexing (VERY IMPORTANT) add  490 μ L of PBS ( Room temperature) dropwise using a P200.



4 Filter DSP working solution using a 40- μ m Flowmi strainer (Sigma, cat. no. BAH136800040-50EA).




Fixation

1h 30m

5 Submerge a ~3×3 mm (the smaller the better) piece of tissue (or organoids) in  500 μ L of the OzSoup and incubate for  00:45:00 at  Room temperature .

45m








5.1 For cells in suspensions, wash cells in cold PBS at least twice (no media + FBS should be present). Pellet cells and resuspend (up to 2 million) cells in  500 μ L of the Oz soup and incubate for  00:30:00 at  Room temperature .

30m






6 At the 45 min mark, add  10 μL of  1 Molarity (M) Tris-HCl  7.5 , mix well by vortexing for 2-3" and sit at  Room temperature for at least  00:15:00 .


15m




7 Pellet the pieces of tissue at  500 rcf, 00:20:00 or a 5-10" in minispinner, and remove supernatant.

20m




7.1 For cells, mix by vortexing 2-3min and pellet cells  500 rcf, Room temperature, 00:05:00 .

5m

8 Add  1 mL of PBS, mix by vortexing 2-3 min, and pellet pieces (or cells) as above, remove supernatant.



8.1 For cells, mix by vortexing 2-3 min and pellet cells  500 rcf, Room temperature, 00:05:00 .

5m








Fixation: For cells only

1h 30m




9 Repeat 8 once more for a total of 2 washes. Continue on step 18-21 below. If sorting or shipping samples follow step 17.

Fixation: For tissues only

50m

10 Add  1 mL of  200 $\mu\text{g}/\text{mL}$ Liberase (Liberase™ Research Grade Sigma-Aldrich 5401127001, <https://www.sigmaaldrich.com/US/en/product/roche/libtmro>) in PBS ( 80 μL of  2.5 mg/mL Liberase +  920 μL PBS).





11 Incubate at  37 °C for  00:30:00 with agitation  800 rpm . At 15 min mark pipette up and down 5 times.

30m



12 After digestion, filter the digestion reaction through a  70 μm mesh.

13 Add  10 mL * ice-cold PBS and pellet cells for  500 rcf, 4°C, 00:05:00 (swinging bucket rotor). Pre-Wash.



5m



**Note**

For pre-wash and Washes 1-3 10 mL is a starting point. One can use less as per preference.



*This volume can be adjusted for small tissues sizes and/or pellets. Optimization may be required.

- 14 Remove supernatant and add  10 mL * of cold PBS+1% BSA and resuspend the pellet before pelleting again ( 500 rcf, 4°C, 00:05:00). Wash 1.

5m

**Note**



*This volume can be adjusted for small tissues sizes and/or pellets. Optimization may be required.

- 15 Remove supernatant and add  10 mL * of cold PBS+1% BSA and resuspend the pellet before pelleting again ( 500 rcf, 4°C, 00:05:00). Wash 2.

5m

**Note**



*This volume can be adjusted for small tissues sizes and/or pellets. Optimization may be required.

- 16 Remove supernatant and add  10 mL * of cold PBS+1% BSA and resuspend the pellet before pelleting again ( 500 rcf, 4°C, 00:05:00). Wash 3.

5m


**Note**

*This volume can be adjusted for small tissues sizes and/or pellets. Optimization may be required.

- 17 **Optional:** If storing for later processing or shipping samples, freeze cells using a cryopreservation strategy as if the main goal was to keep the cells alive. For example, use CryoStor10 and Mr Frosty for slow freezing. Include 1-2 U/uL of RNase inhibitor per sample for storage. Store  -80 °C until use. After storage, thaw in water bath at  37 °C and wash cells twice with PBS+0.5-1%BSA.





- 18 Remove supernatant and resuspend cells in  0.5-1 mL of PBS+1% BSA (optionally add +0.5-1 U/uL RNase Inhibitor).
- 19 Filter cells through Flowmi 40 um.
- 20 Count cells and bring concentration to 1000-1500 cells/uL.
- 21 Load Chromium as per manual.

Note

For ATAC or Multiome kits prpware nuclei using EzLysis Buffer (Sigma-Aldrich, Cat: NUC101-1KT), SaltyEz10/50 protocols ([dx.doi.org/10.17504/protocols.io.bx64prgw](https://doi.org/10.17504/protocols.io.bx64prgw)) or alternatives you are familiar with.