

Jun 26, 2025

Version 2

FLAG-tag mediated co-immunoprecipitation of phosphorylated Rab GTPases with cellular interacting proteins V.2

DOI

<https://dx.doi.org/10.17504/protocols.io.5qpvwjn9l4o/v2>

Eve Napier¹, Pawel Lis², Dario Alessi², Amir Khan¹

¹School of Biochemistry and Immunology, Trinity College Dublin;

²MRC PPU, School of Life Sciences, University of Dundee

AKhanLab



Amir Khan

Trinity College Dublin

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.5qpvwjn9l4o/v2>

Protocol Citation: Eve Napier, Pawel Lis, Dario Alessi, Amir Khan 2025. FLAG-tag mediated co-immunoprecipitation of phosphorylated Rab GTPases with cellular interacting proteins. **protocols.io**

<https://dx.doi.org/10.17504/protocols.io.5qpvwjn9l4o/v2> Version created by [Amir Khan](#)



License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: June 25, 2025

Last Modified: June 26, 2025

Protocol Integer ID: 221081

Keywords: phosphorylated rab gtpase, tagged rab protein, rab protein, interacting proteins leucine rich repeat kinase, subset of rab gtpase, proteins leucine rich repeat kinase, rab gtpase, interactions between rab, phosphorylated flag, interacting protein, rab, protein

Funders Acknowledgements:

Research Ireland

Grant ID: 20/FFP-A/8446

The Michael J. Fox Foundation for Parkinson's Research

Grant ID: ASAP-000463

UK Medical Research Council

Grant ID: MC_UU_00018/1

Abstract

Leucine Rich Repeat Kinase 2 phosphorylates a subset of Rab GTPases at a conserved Ser/Thr residue in their switch 2 region. This modification disrupts interactions between Rabs and effectors/regulators while enabling binding to novel phospho-specific effectors. In this protocol, we use purified phosphorylated and non-phosphorylated FLAG tagged Rab proteins as 'bait' to immunoprecipitate interacting proteins. This is followed by mass spectrometry to identify proteins that bind preferentially to either phospho- or non-phospho Rabs.

Guidelines

Note on lysate preparation: Interacting proteins are immunoprecipitated from both A549 cell lysate and whole mouse brain lysate treated with the LRRK2 inhibitor MLi2. Lysates are prepared with cold lysis buffer freshly supplemented with cOmplete Mini (EDTA-free) protease inhibitor (Roche), PhosSTOP phosphatase inhibitor (Roche) and clarified by centrifuging at 17,000 x g for 15 min at 4°C. Protein concentrations are determined by Bradford Assay and the lysates stored at -80°C.



Materials

Buffers

- PBS-T: PBS supplemented with 0.01% Tween-20, 5mM MgCl₂
- Lysis buffer: 50mM Tris-HCl pH 7.4, 150mM NaCl, 10% (w/v) glycerol, 10mM sodium-glycerophosphate, 10mM sodium pyrophosphate, 0.5% (v/v) NP40 - supplemented with cOmplete Mini (EDTA-free) protease inhibitor (Roche), PhosSTOP phosphatase inhibitor (Roche) and 5mM MgCl₂
- High salt lysis buffer: lysis buffer supplemented with 500mM NaCl

Reagents

- Phosphorylated FLAG-Rab
- FLAG-Rab
- PBS (GIBCO brand)
- Pierce™ Anti-DYKDDDDK Magnetic Agarose (supplied as 25% slurry in PBS, 0.01% Tween-20 detergent, 0.02% sodium azide, pH 7.2 with a confirmed binding capacity of ~3.2mg FLAG Rab (~22kDa)/ml settled beads
- cOmplete Mini (EDTA-free) protease inhibitor (Roche)
- PhosSTOP phosphatase inhibitor (Roche)

Equipment

- Disc rotator
- Magnetic separation stand

Immunoprecipitation of Interacting Proteins

- 1 6 replicates are performed for each phospho-FLAG Rab and FLAG Rab, following the steps below for each replicate.
- 2 Aliquot 5µl settled volume of Pierce™ Anti-DYKDDDDK Magnetic Agarose in low binding Eppendorf tube and wash with 500µl PBS-T on a rotator for 5 minutes.
- 3 Place tube in magnetic separation stand and remove supernatant. Repeat washing step 3 times.
- 4 Add 20µg of FLAG protein to the washed beads and make up the volume to 500µl with PBS-T. Incubate at 4°C on a rotator for 1 hour.
- 5 Place tube in magnetic separation stand and remove supernatant without disturbing the bead pellet.
- 6 Wash the beads with 500µl PBS-T twice and remove supernatant completely.
- 7 Add cell lysate (~400µg A549 lysate or ~750µg brain tissue lysate) to the beads on ice and make up the volume to 500µl with cold lysis buffer. Incubate on a rotator for 1 hour at 4°C.

Note

Note on lysate preparation: Interacting proteins are immunoprecipitated from both A549 cell lysate and whole mouse brain lysate treated with the LRRK2 inhibitor MLI2. Lysates are prepared with cold lysis buffer freshly supplemented with cOmplete Mini (EDTA-free) protease inhibitor (Roche), PhosSTOP phosphatase inhibitor (Roche) and clarified by centrifuging at 17,000 x g for 15 min at 4°C. Protein concentrations are determined by Bradford Assay and the lysates stored at -80°C until required for the IP.

- 8 On a magnetic stand, remove supernatant and wash twice with 500µl cold high salt lysis buffer followed by twice with 500µl cold lysis buffer.

Note

Ensure supernatant is removed completely after the final wash.



- 9 Store samples at -80°C for mass spectrometry analysis or immunoblotting.