



Jun 04, 2025

🌐 Flag-TBK1 N455S and Flag-TBK1 M559R cloning (NEB Site directed mutagenesis kit, #E0554S)

DOI

<https://dx.doi.org/10.17504/protocols.io.81wgbzynogpk/v1>

Jailyn Izu¹

¹NIMH



Christopher Bartley

Translational Immunopsychiatry Unit, NIMH

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.81wgbzynogpk/v1>

Protocol Citation: Jailyn Izu 2025. Flag-TBK1 N455S and Flag-TBK1 M559R cloning (NEB Site directed mutagenesis kit, #E0554S). **protocols.io** <https://dx.doi.org/10.17504/protocols.io.81wgbzynogpk/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: May 15, 2024

Last Modified: March 26, 2026

Protocol Integer ID: 99881

Keywords: mutagenesis of tbk1, tbk1 m559r cloning, mutagenesis kit, directed mutagenesi, tbk1 n455, cloning, tbk1

Funders Acknowledgements:

National Institute of Mental Health

Grant ID: ZIA MH002982

National Institute of Mental Health

Grant ID: R01MH122471

National Institute of Neurological Disorders and Stroke

Grant ID: U01NS120836

Foundation for the National Institutes of Health

Grant ID: Deeda Blair Research Initiative for Disorders of the Brain

Abstract

site-directed mutagenesis of TBK1 variant vectors



Materials

- Template DNA: TBK1 expression plasmids (OmicsLink™ Expression Clone EX-U0617-M09/EX-U0617-M11; GeneCopoeia, Inc., Rockville, MD)
- Flag-tagged EGFP expression plasmid (OmicsLink™ Expression Clone EX-EGFP-M11; GeneCopoeia, Inc., Rockville, MD), will not go through mutagenesis
- HA-tagged EGFP expression plasmid (OmicsLink™ Expression Clone EX-EGFP-M09; GeneCopoeia, Inc., Rockville, MD), will not go through mutagenesis
- Q5® Site-Directed Mutagenesis Kit (NEB, #E0554S)
- Q5 Hot Start High-Fidelity 2X Master Mix (New England BioLabs, #M0494S)
- Forward and reverse primers:

TBK1-N455S mutation:

Forward Primer: GATGATTACAgTGAAACTGTTCAC (Invitrogen)

Reverse Primer: TTTAATTAATTCAATCAGCCATC (Invitrogen)

TBK1-M559R mutation:

Forward Primer: TTAAATTGCAgGACAGAGATTTAC (Invitrogen)

Reverse Primer: CAGGACTTGTAGTTTTTCTAC (Invitrogen)

- Thermocycler ()
- MAX Efficiency DH5α competent *E. coli* cells (Thermo Fisher Scientific, #18258012)
- SOC Medium (Thermo Fisher Scientific, #15544034)
- Miller's LB Broth (Corning, #46-050-CM)
- 14 mL falcon tube (Corning, #352059)
- carbenicillin ()
- NanoDrop Eight Spectrophotometer (ThermoFisher)



Prepare template DNA (Flag- or HA-TBK1 wild-type)

- 1 Follow the general protocol for transforming, inoculating and isolating purified GeneCopoeia ORF cDNA Clones in E. coli:

Protocol

NAME

Transforming, inoculating and isolating purified GeneCopoeia ORF cDNA Clones in E. coli

CREATED BY

Jailyn Izu

Preview

Protocol-specific parameters:

- For **Step 1**, the obtain the following plasmids from GeneCopoeia:

-Flag-tagged TBK1 expression plasmid (OmicsLinkTM Expression Clone EX-U0617-M11; GeneCopoeia, Inc., Rockville, MD)

-HA-tagged TBK1 expression plasmid (OmicsLinkTM Expression Clone EX-U0617-M09; GeneCopoeia, Inc., Rockville, MD)

-Flag-tagged EGFP expression plasmid (OmicsLinkTM Expression Clone EX-EGFP-M11; GeneCopoeia, Inc., Rockville, MD), will not go through mutagenesis

-HA-tagged EGFP expression plasmid (OmicsLinkTM Expression Clone EX-EGFP-M09; GeneCopoeia, Inc., Rockville, MD), will not go through mutagenesis

- For **Step 2**

-Add  1 μL of the 465 ng/ μL purified stock Flag-TBK1 plasmid to  999 μL 10 Mm Tris pH 7.4 to obtain 1 mL of Flag-TBK1 (465 pg/mL)

-For the HA-TBK1 plasmid, add  1 μL of the 262 ng/ μL purified stock plasmid to  499 μL 10 Mm Tris pH 7.4 to obtain 0.5 mL of HA-TBK1 (524 pg/mL).



-Dilute EGFP purified stocks in 10 Mm Tris pH 7.4 to obtain a concentration of 1-10 ng per 5 µl.

Cloning Flag- or HA-TBK1 N455S and M559R

- 2 Follow the general protocol for cloning using the NEB Site-Directed Mutagenesis Kit

Protocol

NAME

Cloning using the NEB Q5® Site-Directed Mutagenesis Kit

CREATED BY

Jailyn Izu

Preview

Protocol-specific parameters:

- For **Step 1**, the template DNA is the TBK1 expression plasmid (OmicsLink™ Expression Clone EX-U0617-M09 or EX-U0617-M11; GeneCopoeia, Inc., Rockville, MD)
- For **Step 1**, the forward and reverse primer sequences are:
 - For the TBK1-N455S mutation:
Forward Primer: GATGATTACAgTGAAACTGTTCAC (Invitrogen)
Reverse Primer: TTTAATTAATTCAATCAGCCATC (Invitrogen)
 - For the TBK1-M559R mutation:
Forward Primer: TTAAATTGCAGGACAGAGATTTAC (Invitrogen)
Reverse Primer: CAGGACTTGTAGTTTTTCTAC (Invitrogen)

Sequence-verify wild-type TBK1 and the mutagenic clones

- 3 Follow the general protocol for plasmid sequencing with GeneWiz:



Protocol

NAME

GeneWiz sequencing

CREATED BY

Jailyn Izu

[Preview](#)

Protocol-specific parameters:

- For **Step 1**, TBK1 is about 8 kB, so prepare 80 ng/μl template DNA
- For **Step 2**, primer sequences are:

For TBK1 plasmids:

For beginning:

- CMV: CGCAAATGGCGGTAGGCGTG (Invitrogen)

For middle:

For **TBK1 wild-type**:

- Forward: AGGGGAAGATGGACAGTCTGTGT (Invitrogen)
- Reverse: CCACCTGGAGATAATCTGCTGTCTG (Invitrogen)

For **TBK1 N455S**:

- Forward: GATGATTACAGTGAAACTGTTCAC (Invitrogen)
- Reverse: TTTAATTAATTCAATCAGCCATC (Invitrogen)

For **TBK1 M559R**:

- Forward: TTAAATTGCAGGACAGAGATTTAC (Invitrogen)
- Reverse: CAGGACTTGTAAGTTTTTCTAC (Invitrogen)

For end:

- T7: TAATACGACTCACTATAGGG (Invitrogen)