

Sep 25, 2024

FLAM-seq (with Kapa mRNA enrichment)_CMS_edit-2024-09-25

Forked from a private protocol

DOI

<https://dx.doi.org/10.17504/protocols.io.j8nlk8q46l5r/v1>

Corneliu Sologon¹

¹University of Miami



Corneliu Sologon

University of Miami

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.j8nlk8q46l5r/v1>

Protocol Citation: Corneliu Sologon 2024. FLAM-seq (with Kapa mRNA enrichment)_CMS_edit-2024-09-25. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.j8nlk8q46l5r/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: In development

We are still developing and optimizing this protocol



Created: September 25, 2024

Last Modified: September 25, 2024

Protocol Integer ID: 108373

Keywords: kapa mrna, flam, seq

Abstract

fork__CMS_edit-2024-09-25

Materials

Kapa mRNA HyperPre Kit Illumina (cat. KK8440, Roche; Replaces TruSeq mRNA)

- mRNA Capture Beads
- mRNA Bead Binding Buffer (BBB)
- mRNA Bead Wash Buffer(BWB)

Poly(A) Tail Length Assay Kit (cat. 764551KT, Thermo Fisher)

- 5X tail buffer mix
- 10X tail enzyme mix
- Universal Reverse (Univ. RV Primer; used in Sect 4: cDNA Amplification)

SMARTer PCR cDNA cDNA Synthesis Kit (cat. 634926, Takara)

- 5X First strand buffer
- DTT 20mM
- dNTP mix 10 mM
- RNase Inhibitor
- SMARTScribe RT
- 5' PCR Primer II A (used in Sect 4: cDNA Amplification)

Advantage 2 DNA polymerase mix (cat 639207, Takara)

- 10X Advantage 2SA PCR buffer
- dNTP Mix (10 mM each)

RNAClean XP Beads (cat. A63987, Beckman Coulter)

XP DNA beads (cat. A63881, Beckman Coulter)

Custom Oligonucleotides

isoTSO - Template switch oligo (where "i" indicates stereoisomers of dC and dG as in the associated publication):

iCiGiCAAGCAGTGGTATCAACGCAGAGTACATrGrGrG

RT primer 1:

GGTAATACGACTCACTATAGCGAGANNNNNNNNNNNNCCCCCCCCCTTT

PCR primer 1 (to use in combination with RT primer 1; not used; replaced by 5' PCR Primer II A):

GGTAATACGACTCACTATAGCGAG

RT primer 2 (to use in alternative to 1):

TGAGTCGGCAGAGAACTGGCGAANNNNNNNNNNNNCCCCCCCCCTTT,

PCR primer 2 (to use in combination with RT primer 2; not used; replaced by 5' PCR Primer II A):

TGAGTCGGCAGAGAACTGGCGAA



Poly(A)+ RNA preparation (using Kapa mRNA Beads)

10m 10s

1 Prepare mRNA Beads

1.1  52.5 μ L mRNA Beads can be scaled for multiple samples

 500 rpm, 00:01:00

1.2  Room temperature  00:05:00 magnet

Discard supernatant

5m

1.3 Remove from magnet

 52.5 μ L BBB

 500 rpm, 00:01:00

1.4  Room temperature  00:05:00 magnet

Discard supernatant

1.5 Remove from magnet

 52.5 μ L BBB

 500 rpm, 00:01:00

2  2-10 μ g  ds DNA in 0.2 mL PCR tube

 50 μ L QS; NAF

3  50 μ L re-suspended mRNA Beads (from 1.5)

 1000 rpm, 20°C, 00:01:00

4  65 °C  00:02:00

 20 °C  00:05:00

 Undetermined, 00:00:10 , pulse

7m 10s

5  Room temperature  00:05:00 magnet

Discard supernatant

5m



- 6 Remove from magnet
🧪 200 μ L BWB
🌀 1000 rpm, 20°C, 00:01:00
- 7 🌡 Room temperature ⌚ 00:05:00 magnet 5m
Discard supernatant
- 8 🧪 50 μ L NAF
🌀 1000 rpm, 20°C, 00:01:00
- 9 🌡 70 °C ⌚ 00:02:00 7m 10s
🌡 20 °C ⌚ 00:05:00
🔴 Undetermined, 00:00:10 , pulse
- 10 🧪 50 μ L BBB 5m
🌀 1000 rpm, 20°C, 00:01:00
🌡 20 °C ⌚ 00:05:00
- 11 🌡 Room temperature ⌚ 00:05:00 magnet 5m
Discard supernatant
- 12 Remove from magnet
🧪 200 μ L BWB
🌀 1000 rpm, 20°C, 00:01:00
- 13 🌡 Room temperature ⌚ 00:05:00 magnet 5m
Discard supernatant
- 14 Remove from magnet
🧪 16 μ L NAF
🌀 1000 rpm, 20°C, 00:01:00
- 15 🌡 70 °C ⌚ 00:02:00 2m
- 16 🌡 Room temperature ⌚ 00:05:00 magnet 5m



17  16 μ L  ds DNA to new 0.2 mL PCR tube

GI Tailing (using USB poly(A) tail length assay, Thermo Fisher).

18 Prepare on ice (20 μ L total):

 14 μ L RNA - poly(A) +

 4 μ L 5X tail buffer mix

 2 μ L 10X tail enzyme mix

19  37 $^{\circ}$ C  01:00:00

1h

20  1.5 μ L Stop Solution

2m

 4 $^{\circ}$ C  00:02:00

21  38.7 μ L XP RNA Beads (1.8X)

5m

 Room temperature  00:05:00

22  Room temperature  00:03:00 magnet

3m

Discard supernatant

23  50 μ L 80% EtOH  00:00:30

30s

Discard supernatant

24  go to step #23 x1

25  Undetermined, 00:00:10 , pulse

10s

Discard supernatant

26 Remove from magnet

 17 μ L NAF

 1000 rpm, 20 $^{\circ}$ C, 00:01:00

27  16 μ L  ds DNA to new 0.2 mL PCR tube



cDNA synthesis (using SMARTscribe reverse transcriptase kit, Clontech).

28 Prepare RT master mix at **room temperature** (22 uL total; 1.1X per sample):

🧴 8 µL 5X First STrand buffer

🧴 1.5 µL DTT (20 mM)

🧴 4 µL dNTP mix (10 mM)

🧴 2 µL RNase Inhibitor

🧴 2 µL isoTSO (12 uM)

🧴 2 µL SMARTScribe RT

🧴 2.5 µL NAF

29 In 0.2 mL PCR tube from **step 27**

🧴 16 µL 🧬 ds DNA

🧴 2 µL RT Primer 1 (10 uM); custom 🧬 ds DNA

30 🔄 1000 rpm, Room temperature , 00:01:00

10s

⚙️ Undetermined, 00:00:10 , pulse

31 🌡️ 72 °C ⌚ 00:03:00

1h 13m

🌡️ 42 °C **HOLD** - add 🧴 22 µL RT master mix

🌡️ 42 °C ⌚ 01:00:00

🌡️ 70 °C ⌚ 00:10:00

🌡️ 4 °C **HOLD**

32 🧴 24 µL XP DNA Beads (0.6X)

5m

🌡️ Room temperature ⌚ 00:05:00

33 🌡️ Room temperature ⌚ 00:03:00 magnet

Discard supernatant

34 🧴 50 µL 80% EtOH ⌚ 00:00:30

Discard supernatant

35 ➡️ go to step #23 x1



- 36  Undetermined, 00:00:10 , pulse
Discard supernatant
- 37 Remove from magnet
 42 μL NAF
 1000 rpm, 20°C, 00:01:00
- 38  40 μL  ds DNA to new 0.2 mL PCR tube

cDNA library amplification (using Advantage 2 PCR enzyme system, Clontech).

11m 25s

- 39 Prepare the cDNA library amplification rxn:
 40 μL  ds DNA
 10 μL 10X Advantage 2SA PCR buffer
 2 μL dNTP (10 mM)
 2 μL 5' PCR Primer II A (SMARTer PCR kit)
 2 μL Univ. RV Primer (Poly(A) Tail Kit)
 42 μL NAF
 2 μL Advantage 2 Polymerase Mix
- 40 Perform PCR
-  98 °C **HOLD** - add cDNA rxn
 -  98 °C  00:01:00
 -  98 °C  00:00:10
 -  63 °C  00:00:15
 -  68 °C  00:03:00 **repeat 40.3 -40.5 24X**
 -  68 °C  00:07:00
 -  4 °C **HOLD**
- 41  60 μL XP DNA Beads (0.6X)
 Room temperature  00:05:00
- 42  Room temperature  00:03:00 magnet

11m 25s



Discard supernatant

43



200 μ L 80% EtOH



00:00:30

Discard supernatant

44



go to step #23 x1

45



Undetermined, 00:00:10 , pulse

Discard supernatant

46

Remove from magnet



42 μ L NAF



1000 rpm, 20°C, 00:01:00

47



40 μ L



ds DNA

to new 0.2 mL PCR tube

Sequencing library preparation (using the SMRTbell™ Template Prep Kit)

48

Performed in PacBio sequencing core