



Sep 21, 2022

## Flow cytometry

DOI

[dx.doi.org/10.17504/protocols.io.kqdg3956eg25/v1](https://dx.doi.org/10.17504/protocols.io.kqdg3956eg25/v1)

gurvir.virdi<sup>1</sup>

<sup>1</sup>UCL Institute of Neurology



gurvir.virdi

### Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.kqdg3956eg25/v1>

**Protocol Citation:** gurvirdi 2022. Flow cytometry. **protocols.io**  
<https://dx.doi.org/10.17504/protocols.io.kqdg3956eg25/v1>

**License:** This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working

**We use this protocol and it's working**

**Created:** September 21, 2022



**Last Modified:** May 31, 2024

**Protocol Integer ID:** 70344

**Keywords:** ASAPCRN, midbrain dopaminergic neurons for flow cytometry analysis, downstream flow cytometry acquisition, flow cytometry, flow cytometry analysis, flow cytometry this protocol, midbrain dopaminergic neuron, dopaminergic neuron, neuron, fixed midbrain, midbrain







## Abstract

This protocol is for immunolabelling fixed midbrain dopaminergic neurons for flow cytometry analysis and downstream flow cytometry acquisition.







## Troubleshooting



## Sample preparation

- 1 Cells were washed once with PBS
- 2 The cells were incubated with Accutase (Gibco) to generate a single-cell suspension for  00:05:00 at  37 °C 5m
- 3 A cell suspension of 500k/ml was prepared in media
- 4 Cells were then spun down at  200 x g, 00:05:00 , and the supernatant was removed. 5m
- 5 Cell pellet was resuspended gently in 4ml of 4% paraformaldehyde and briefly vortexed at a low speed before being spun on a rotation spinner for  00:10:00 at room temperature. 10m
- 6 After fixation, samples were spun down and supernatant removed
- 7 Cells were resuspended in 2ml of 0.1% BSA (Sigma) in PBS.
- 8 After resuspension, cells were filtered through a 70µm strainer (Miltenyi Biotec) to filter out any cell clumps
- 9 Cells were then centrifuged  200 x g, 00:05:00 , and the supernatant was removed. 5m
- 10 Cell pellets were then resuspended in 1ml of permeabilization/blocking buffer (0.1% Triton X-100, 1% BSA, 10% normal goat serum (Sigma) in PBS), and incubated on a rotation spinner for  00:30:00 at room temperature. 30m



- 11 After permeabilization/blocking, cells were centrifuged  200 x g, 00:05:00 and the supernatant was removed. 5m
- 12 Cells were then resuspended in the primary antibodies (1:200) made up in 0.1% BSA in PBS, and incubated on the rotation spinner for  01:00:00 at room temperature. 1h
- 13 After primary antibody incubation, cells were centrifuged  200 x g, 00:05:00 , supernatant removed and washed once in 0.1% BSA in PBS. 5m
- 14 Cells were then resuspended in the species-specific secondary antibodies (AlexaFluor 488, 647) at a dilution of 1:500 made up in 0.1% BSA in PBS and incubated in the dark on a rotation spinner for  00:30:00 . 30m
- 15 After incubation, cells were centrifuged  200 x g, 00:05:00 , supernatant removed and washed once in PBS, followed by incubation with DAPI made up in PBS for  00:05:00 . 10m
- 16 The DAPI + PBS was then removed, followed by one wash in PBS, before being analysed on the flow cytometer.

## Cell sorting and analysis

- 17 The samples were run on the LSRii (BD) cell sorter. Scattering was initially used to discard debris as well as cell doublets and larger clumps.
- 18 The single-cell population was then gated to include DAPI positive only cells (negative control).
- 19 The gating threshold for measured channels was determined using the control lacking the antibody of interest (Fluorescence minus one (FMO) control), for both channels being recorded.
- 20 Once the parameters had been set, 10,000 cell events were recorded, and data were processed and analysed on the FlowJo software.