

Apr 20, 2024

Version 3

Flow Cytometry ICS Nuclear Antigens V.3

DOI

<https://dx.doi.org/10.17504/protocols.io.rm7vz3qdxgx1/v3>



Michael Betts¹, Gregory Golden¹

¹University of Pennsylvania

Human Islet Research N...



Sandy Beshir

City of Hope

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.rm7vz3qdxgx1/v3>

Protocol Citation: Michael Betts, Gregory Golden 2024. Flow Cytometry ICS Nuclear Antigens. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.rm7vz3qdxgx1/v3> Version created by **Sandy Beshir**

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: April 20, 2024

Last Modified: April 20, 2024

Protocol Integer ID: 98530

Keywords: extracellular epitope staining (ECS), intra-cellular epitope staining (ICS), HIRN, flow cytometry ics nuclear antigen, standardized protocol for cryopreserved sample recovery, cryopreserved sample recovery, flow cytometry, combined analysis of extracellular epitope, parameter flow cytometry, ics procedures for cell, extracellular epitope, preserved cell suspension, proper cryopreservation, cell suspension, cryopreservation, cell population, labeled antibody, discrete cell population, biological sample, sample recovery, cell, ics procedure

Funders Acknowledgements:

NIH

Grant ID: U01 DK112217

Abstract

High-parameter flow cytometry enables identification and characterization of a wide range of cell populations within a biological sample. A combined analysis of extracellular epitope staining (ECS) and intra-cellular epitope staining (ICS) using a collection of fluorophore-labeled antibodies sufficiently identifies discrete cell populations and their respective phenotypes. Importantly, ECS/ICS can be applied to cryo-preserved cell suspensions recovered in tissue culture media, enabling samples to be conveniently analyzed after collection and proper cryopreservation. However, consistent cryopreserved sample recovery and ECS procedures are critical to data comparison across multiple experiments. Herein, we describe a standardized protocol for cryopreserved sample recovery and ECS/ICS procedures for cell-surface and intra-cellular epitope labeling.



Materials

Materials Required

1.  1x phosphate saline buffer (PBS) **Corning Catalog #21-031-CM**
2.  Bovine serum albumin (BSA) **Gemini Bio-Products Catalog #700-101P**
3.  Sodium Azide **Fisher Scientific Catalog #S2271-500**
4. FACS buffer (1xPBS, 10g/L BSA, 1 g/L sodium azide)
5.  RPMI **Corning Catalog #10-040-CM**
6.  Fetal calf serum (FCS) **Gemini Bio-Products Catalog #900-108**
7.  Penicillin/streptomycin 10000 U/mL penicillin 10000 µg/mL streptomycin **Lonza Catalog #17-602E**
Lonza ([RRID:SCR_000377](https://pubchem.ncbi.nlm.nih.gov/compound/SCR_000377))
8. "R10" medium (RPMI, 10% FCS, 1% penicillin/streptomycin)
9.  DNase I **Roche Catalog #04716728001**
10.  Live/Dead fixable Aqua Dead Cell Stain kit **Invitrogen - Thermo Fisher Catalog #L34966**
diluted 1:60 in 1x PBS (prepare fresh each day from DMSO stock)
11.
 e-Bioscience Foxp3 / Transcription Factor Staining Buffer Set **Invitrogen - Thermo Fisher Catalog # 00-5523-00**
12.  Paraformaldehyde (PFA) **Electron Microscopy Sciences Catalog #15712-S** diluted to 1% in PBS
13. Fluorophore-labelled antibodies of choice
14. Hemocytometer

Procedure

1 Thawing and Resting

a. Pre-warm R10 media in a  37 °C water bath.

b. Thaw samples in-vial using a  37 °C water bath.

c. Add thawed cells to  14 mL of R10, then spin cells at 500 xg for 5 min.

d. Resuspend cell pellet in  3 mL of R10 and count cells.

e. Rest cells at least 3 hours (up to overnight) at 2×10^6 cells/mL in R10 medium + 1 μ L/mL DNase I at  37 °C , 5% CO₂.

Note: during resting, prepare antibody cocktail master mix. Adjust volume of ECS for  50 μ L per test with FACS buffer.

f. After resting, add PBS up to  15 mL or  50 mL (whichever is closer, rounding up) to cells and transfer to  15 mL or  50 mL conical tube.

g. Spin cells at 500 xg for 5 minutes at room temperature (RT).

h. Resuspend cells in PBS to 1×10^7 cells/mL and count. If cells are too dilute, re-spin cells and resuspend at 1×10^7 cells/mL. Transfer  200 mL of cells (2×10^6 cells) into each well of a V-bottom 96 well plate.

2 Viability and extracellular staining (ECS)

a. Spin plate at 500 xg for 5 minutes at RT

b. Using a multichannel pipette, carefully remove the supernatant.

c. Add  5 μ L of 1:60 Aqua viability dye directly to cell pellet and resuspend cells.

d. Incubate for 10 minutes at RT in the dark.



- e. Add  50 μL of ECS antibody cocktail to cells and incubate for 20 minutes at RT in the dark (prepare 1% PFA fixation buffer in the meantime).
- f. Add  100 μL of FACS buffer to each well and spin plate at 500 xg for 5 minutes. Remove supernatant.

3 Fixation, Permeabilization, and ICS Staining

- a. Fix cells with  100 μL of 1xFixation/permeabilization buffer (1 part concentrate + 3 parts diluent) for 30 minutes at RT in the dark.
Note: Make ICS abs here
- b. Add  100 μL of 1x perm/wash buffer and centrifuge 800 xg for 5 min.
- c. Remove supernatant and add  100 μL ICS cocktail (made in perm/wash buffer, made with diH₂O) to the cells.
- d. Incubate for 1h at RT in the dark.
- e. Add  100 μL perm/wash buffer. Centrifuge at 800 xg for 5 min.
- f. Discard supernatant and resuspend pellet in  200 μL 1% PFA and transfer to FACS tubes.
- g. Wrap in aluminum foil and store at  4 °C until flow cytometry.