

Jul 18, 2018

## Fluorescence analysis using CF imager

DOI

[dx.doi.org/10.17504/protocols.io.q4hdyt6](https://dx.doi.org/10.17504/protocols.io.q4hdyt6)

Steven J Burgess<sup>1</sup>

<sup>1</sup>University of Illinois at Urbana-Champaign



Steven J Burgess

University of Illinois at Urbana-Champaign

---

OPEN  ACCESS



DOI: [dx.doi.org/10.17504/protocols.io.q4hdyt6](https://dx.doi.org/10.17504/protocols.io.q4hdyt6)

**Protocol Citation:** Steven J Burgess 2018. Fluorescence analysis using CF imager. **protocols.io**  
<https://dx.doi.org/10.17504/protocols.io.q4hdyt6>

**License:** This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** In development

**We are still developing and optimizing this protocol**

**Created:** June 18, 2018

**Last Modified:** August 14, 2018

**Protocol Integer ID:** 13161



## Preparing plants

- 1 Dark adapt plants for at least 20 minutes prior to taking measurements.

### Note

This is done to ensure the photosynthetic electron transport chain is fully oxidized and reaction centres are open. In an ideal situation plants are allowed to dark adapt overnight prior to measurement.

A properly adapted, healthy plant should give a Fv/Fm value of ~0.8. This has been shown to be highly stable between species. Significant deviation from this (e.g. <0.7) either suggests incomplete adaptation or stressed plant material.

- 2 Turn on the cf imager and open the FluorImager software.






3

Start the FluorImager software



4

The surface of the leaf should be 140mm from the base of the imaging chamber, and can be adjusted by lowering or raising the plant under analysis. Position plant/leaf in the chamber




## Set focus

- 5 Set the focus by adjusting the dial above the chamber, and lock in position by turning the screw on the side.


## Set exposure

- 6 Manually adjust the aperture as shown on the right to allow an optimal amount of light into the imager so as not to overexpose measurements




7 Isolate the plant or leaf of interest


right click  
on image




After applying isolation background areas will be masked out in blue as shown below.



- 8 Remove residual noise from image
- Sometimes the software picks up background noise as real signals, shown below as black dots on the blue background. It is advisable to mask these, otherwise they will be counted as a separate 'colony' during analysis and measurements will be recorded for each of these spots in the final data sheet. Noise can be masked by moving the cursor over the dot, pressing CTRL+left click simultaneously.




## 9 Load protocol



10

11





- 12 Once you are happy with the scheduled program click on the protocol icon in the toolbar (black arrow below) to start the run



- 13 Export data

