

Apr 26, 2024

Version 1

Fluorescence assay for Enterovirus D68 3C protease activity measurement V.1

DOI

<https://dx.doi.org/10.17504/protocols.io.261ge54jyg47/v1>

Charline Giroud¹, Oleg Fedorov¹

¹Center for Medicines Discovery, University of Oxford

ASAP Discovery



Charline Giroud

Center for Medicines Discovery

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.261ge54jyg47/v1>

Protocol Citation: Charline Giroud, Oleg Fedorov 2024. Fluorescence assay for Enterovirus D68 3C protease activity measurement. [protocols.io https://dx.doi.org/10.17504/protocols.io.261ge54jyg47/v1](https://dx.doi.org/10.17504/protocols.io.261ge54jyg47/v1)

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: April 26, 2024

Last Modified: April 26, 2024

Protocol Integer ID: 98854

Keywords: FRET, enzymatic assay, protease assay, Picornaviridae , Enterovirus, enterovirus d68 3c protease activity measurement, enterovirus d68 3c protease, activity of viral protease, enterovirus d68, viral protease, 3c protease activity measurement, protease activity measurement, detection of the fluorescence, cleavage of the peptide, sequence specific to the tested protease, tested protease, fluorescence emission, fluorescence, 3c protease, peptide, assay, protease, labelled peptide, proximity of dabcyl, fluorophore edan, d68, terminal by the fluorophore edan, quencher dabcyl, d68 3c

Funders Acknowledgements:

National Institutes of Health/National Institute Of Allergy and Infectious Diseases (NIH/NIAID)

Grant ID: U19AI171399

Disclaimer

The content is solely the responsibility of the authors and does not necessarily represent the official views of the National Institutes of Health.

Abstract

This protocol details the fluorescence assay for Enterovirus D68 (EV-D68) 3C protease activity measurement. This method is intended to measure the activity of viral proteases by using a specific labelled-peptide that allows the detection of the cleaved product. The substrate contains the cleavage-sequence specific to the tested protease and is labeled in C-terminal by the fluorophore Edans (ex 336 nm; em: 455 nm) and in N-terminal by the quencher Dabcyl (ex 472 nm). In the case of a non-cleaved substrate, the proximity of Dabcyl to Edans prevents the emission and the detection of the fluorescence at 455 nm. The cleavage of the peptide by the protease allows Edans' fluorescence emission and detection.

Attachments



[nyjwb5id7.docx](#)







203KB

Materials

Reagents

- **Assay buffer:**


	A	B
	Tris pH 7.0	50 mM
	NaCl	150 mM
	Glycerol	10%
	DTT (optional)	0.5 mM

- **Incubation:**  01:00:00 at  Room temperature .
- **EV-D68 3C:** Protein stocks were stored at  -80 °C and used as 2x solution ( 1 micromolar (μM) ,  0.5 micromolar (μM) final assay concentration) in assay buffer.
- **Substrate:** Dabcyl-KEALFQGPPQFE-Edans (LifeTein, USA) prepared as a stock solution at  5 millimolar (mM) in DMSO and used at 2x solution ( 40 micromolar (μM) ,  20 micromolar (μM) final concentration assay concentration) in assay buffer.
- **Positive control:** GC376 (Pubchem CID 71481119),  50 micromolar (μM) top final assay concentration.
- **Plates:** ProxiPlate-384 Plus, white, Greiner cat# 6008280.
- **Liquid handler:** Echo® acoustic liquid handler (Beckman Coulter, USA).
- **Plate reader:** Pherastar FS, BMG Labtech (Germany), 350-460 FI optic module, the plate is read every 30 s for 2 hours and shacked during 5 s before the first reading.

EV-D68 3C Pro Assay

3h

- 1 This method is intended to measure the activity of viral proteases by using a specific labelled-peptide that allows the detection of the cleaved product. The substrate contains the cleavage-sequence specific to the tested protease and is labeled in C-terminal by the fluorophore Edans (ex 336 nm; em: 455 nm) and in N-terminal by the quencher Dabcyl (ex 472 nm). In the case of a non-cleaved substrate, the proximity of Dabcyl to Edans prevents the emission and the detection of the fluorescence at 455 nm. The cleavage of the peptide by the protease allows Edans' fluorescence emission and detection.



Asset URL:

Reagents and equipment

- 2 **Assay buffer:** [1M] 50 millimolar (mM) Tris pH 7.0, [1M] 150 millimolar (mM) NaCl, [1M] 10 % (v/v) glycerol and [1M] 1 millimolar (mM) TCEP (optional).
Incubation: ⌚ 01:00:00 at 🌡 Room temperature .
EV-D68 3C: protein stocks were stored at -80C and used as 2x solution ([1M] 1 micromolar (μM) , [1M] 0.5 micromolar (μM) final assay concentration) in assay buffer.
Substrate: Dabcyl-KEALFQGPPQFE-Edans (LifeTein, USA) prepared as a stock solution

1h



at 10 millimolar (mM) in DMSO and used at 2 Mass Percent solution (40 micromolar (μ M) , 20 micromolar (μ M) final concentration assay concentration) in assay buffer.

Plates: Non-binding, black 384-plate, small volume Greiner for the assay. Non-binding, 96-well polypropylene plates for the serial dilutions of the protease and substrate.


Plate reader: Pherastar FS, BMG Labtech (Germany), 350-460 FI optic module, the plate is read every 30 seconds for 2 hours and shaken 5 s before each reading.

EV-D68 3C Pro assay development and optimization

3h

- 3 This step aims to minimize the amount of costly reagents used for the assay and to find the best protease/substrate ratio allowing a high robustness determined by Z' calculation. A serial dilution of the protease is added to a substrate serial dilution and monitored for fluorescence intensity.
- 4 Prepare a serial dilution of the protease in assay buffer starting from 1 micromolar (μ M) 2 Mass Percent , at least 1000 μ L , and perform a cascade dilution by transferring 500 μ L of the first 2x solution to the next well which contains 500 μ L of buffer. Transfer 10uL of 2x solution onto 384 well plate by repeat-pipetting, from the left to the right to cover several columns of the plate.
- 5 Prepare a serial dilution of the substrate starting from a solution of 40 micromolar (μ M) 2 Mass Percent , 1000 μ L , by transferring in cascade 500 μ L to the next well containing the same amount of buffer. Transfer 10uL of the 2x solution onto the protease solution by repeat-pipetting from the bottom to the top of the plate.
- 6 Read the fluorescence intensity at 350/460 nm every 00:00:30 for 02:00:00 in kinetic mode with a shaking step of the plate before each measurement. 2h 0m 30s
- 7 Plot the fluorescent intensity as a function of time for each substrate concentration and determine a reaction time for the analysis which shows a complete (or almost complete) reaction. In the example below, the reaction time for the analysis is set to 120 min.

8



Asset URL:

- 9 Calculate the mean (μ) and the standard deviation (σ) of fluorescence intensity and then calculate the signal-to-background ratio and the Z' or Zfactor. (s: signal; c: control).

$$\text{Estimated Z-factor} = 1 - \frac{3(\hat{\sigma}_s + \hat{\sigma}_c)}{|\hat{\mu}_s - \hat{\mu}_c|}$$

Asset URL:

10

11

		Substrate (μM)															
		20				10				5				2.5			
		average fluo	stdev fluo	S/BG	Z'	average fluo	stdev fluo	S/BG	Z'	average fluo	stdev fluo	S/BG	Z'	average fluo	stdev fluo	S/BG	Z'
[EV-D68 3C] (μM)	0.5	29967	3247	10.340	0.633	16832	605	8.491	0.861	10259	688	6.537	0.751	5350	287	4.271	0.744
	0.25	29931	2641	10.327	0.700	16922	842	8.537	0.814	10147	95	6.466	0.955	5883	328	4.696	0.746
	0.125	28805	2736	9.939	0.676	15629	1136	7.884	0.732	10234	356	6.521	0.866	5415	323	4.323	0.721
	0.0625	17527	1026	6.048	0.776	10653	453	5.374	0.814	6877	127	4.382	0.910	4163	371	3.323	0.552
	0	2898	67	ND	ND	1982	84	ND	ND	1569	32	ND	ND	1253	63	ND	ND

Asset URL:

Choose the substrate and protease concentrations that allow the best signal-to-background and the highest Z', in this example: 5 μM substrate and 0.25 μM protease.

EV-D68 3C Pro IC50 Measurement

3h

- 12 Add  10 μL of 2x protein [M] 0.5 micromolar (μM) solution to each well containing the compounds to be tested previously dispensed onto the plate. 
- 13 Incubate the mix for  01:00:00 at  Room temperature and initiate the enzymatic reaction by the addition of  10 μL of 2x ([M] 10 micromolar (μM)) substrate solution using the plate reader injector. 
- 14 Read the fluorescence intensity at 350/460 nm every 30 seconds for  02:00:00 in kinetic mode, which includes a shaking step of the plate before each measurement. 
- 15 Calculate the IC50 by plotting the initial velocity against various concentrations of tested inhibitor by using a four-parameter dose-response curve in Prism (v8.0) software.

