

Jun 20, 2025

Version 2

Fluorescence assay for MERS-CoV Mpro protease activity measurement for screening against antiviral molecules V.2

DOI

<https://dx.doi.org/10.17504/protocols.io.8epv5rzm4g1b/v2>

Charline Giroud^{1,2,3}, Oleg Fedorov^{1,2,3}

¹Center for Medicines Discovery; ²University of Oxford; ³ASAP Discovery Consortium

ASAP Discovery



Charline Giroud

Center for Medicines Discovery

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.8epv5rzm4g1b/v2>

Protocol Citation: Charline Giroud, Oleg Fedorov 2025. Fluorescence assay for MERS-CoV Mpro protease activity measurement for screening against antiviral molecules. **protocols.io**

<https://dx.doi.org/10.17504/protocols.io.8epv5rzm4g1b/v2> Version created by [Mary-Ann Xavier](#)

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: February 09, 2025

Last Modified: June 20, 2025

Protocol Integer ID: 119854

Keywords: FRET assay, enzymatic assay, viral proteases, Coronavirus, MERS-CoV , MPro, IC50, activity of viral protease, viral protease, cov mpro protease activity measurement, protease activity measurement, sequence specific to the tested protease, protocol for protein expresssion, antiviral molecules this protocol detail, cleavage of the peptide, tested protease, protein production protocol, detection of the fluorescence, protein expresssion, antiviral molecule, protease, fluorescence, protein production, fluorescence emission, labelled peptide, peptide, assay, proximity of dabctl, purification, protein

Funders Acknowledgements:

National Institutes of Health/National Institute Of Allergy and Infectious Diseases (NIH/NIAID)

Grant ID: U19AI171399

Disclaimer

The content is solely the responsibility of the authors and does not necessarily represent the official views of the National Institutes of Health.

Abstract

This protocol details the fluorescence assay for MERS-CoV Mpro protease activity measurement. This method is intended to measure the activity of viral proteases by using a specific labelled peptide that allows the detection of the cleaved product. The substrate contains the cleavage-sequence specific to the tested protease and is labeled in C-terminal by the fluorophore Edans (ex 336 nm; em: 455 nm) and in N-ternimal by the quencher Dabctl (ex 472 nm). In the case of a non-cleaved substrate, the proximity of Dabctl to Edans prevents the emission and the detection of the fluorescence at 455 nm. The cleavage of the peptide by the protease allows Edans' fluorescence emission and detection. In this current version we added the protein production protocol and the plasmid used for this assay. In this new version, we added the protocol for protein expresssion and purification and the Addgene id.

Attachments



[nyjwb5id7.docx](#)

203KB



Materials

Reagents

- **Assay buffer:** [M] 25 millimolar (mM) HEPES pH 7.5, [M] 150 millimolar (mM) NaCl, [M] 10 millimolar (mM) glycerol, [M] 7 millimolar (mM) DTT, [M] 5 millimolar (mM) EDTA, [M] 1 millimolar (mM) TCEP (optional) and [M] 0.05 % volume Tween-20 .
- **Incubation:** ⌚ 01:00:00 at 🌡 Room temperature .
- **MERS-CoV Mpro:** Protein stocks were stored at 🌡 -80 °C and used as 2x solution ([M] 10 micromolar (μM) , [M] 5 micromolar (μM) final assay concentration) in assay buffer.
- **Substrate:** Edans-GVLQSGLV-LysDabcyI-K (LifeTein, USA) prepared as a stock solution at [M] 10 millimolar (mM) in DMSO and used at 2x solution ([M] 20 micromolar (μM) , [M] 10 micromolar (μM) final concentration) in assay buffer.
- **Positive control:** GC376 (Pubchem CID 71481119), [M] 50 micromolar (μM) top final assay concentration.
- **Plates:** Greiner, 384-well plate, black, small volume, non-binding #7840900
- **Liquid handler:** Echo® acoustic liquid handler (Beckman Coulter, USA).
- **Plate reader:** Pherastar FS, BMG Labtech (Germany), 350-460 FI optic module, the plate is read every 30 s for 2 hours and shaken during 5 s before the first reading.

Protocol materials

🔗 MERS-CoV Mpro [addgene Catalog #228646](#)



Protein production

2d

- 1 The protein used in this assay was expressed and purified following the protocol below with plasmid  MERS-CoV Mpro **addgene Catalog #228646**

Protocol

NAME

MERS-CoV Mpro large scale purification protocol

CREATED BY


korvus.wang

[Preview](#)

MERS-CoV Mpro Assay

3h

- 2 This method is intended to measure the activity of viral proteases by using a specific labelled-peptide that allows the detection of the cleaved product. The substrate contains the cleavage-sequence specific to the tested protease and is labeled in C-terminal by the fluorophore Edans (ex 336 nm; em: 455 nm) and in N-terminal by the quencher Dabcyl (ex 472 nm). In the case of a non-cleaved substrate, the proximity of Dabcyl to Edans prevents the emission and the detection of the fluorescence at 455 nm. The cleavage of the peptide by the protease allows Edans' fluorescence emission and detection.



Asset URL:

Reagents and equipment

3h 0m 35s

3 **Assay buffer:** [1] 25 millimolar (mM) HEPES pH 7.5, [1] 150 millimolar (mM) NaCl, [1] 10 millimolar (mM) glycerol, [1] 7 millimolar (mM) DTT, [1] 5 millimolar (mM) EDTA, [1] 1 millimolar (mM) TCEP (optional) and [1] 0.05 % volume Tween-20 .

3h 0m 35s

Incubation: ⌚ 01:00:00 at room temperature.

MERS-MPro: protein stocks were stored at -80C and used as [1] 2 x solution ([1] 10 micromolar (μM) , [1] 5 micromolar (μM) final assay concentration) in assay buffer.

Positive control: GC376 (Pubchem CID [71481119](https://pubchem.ncbi.nlm.nih.gov/compound/GC376)), [1] 25 micromolar (μM) top final assay concentration.

Substrate: Edans-GVLQSGLV-LysDabcyl-K (LifeTein, USA) prepared as a stock solution at 10 mM in DMSO and used at 2x solution ([1] 20 micromolar (μM) , [1] 10 micromolar (μM) final concentration assay concentration) in assay buffer.

Liquid handler: Echo[®] acoustic liquid handler (Beckman Coulter, USA).

Plate reader: Pherastar FS, BMG Labtech (Germany), 350-460 FI optic module, the plate is read every ⌚ 00:00:30 for ⌚ 02:00:00 and shaken ⌚ 00:00:05 before each reading.



MERS-MPro IC50 Measurement


3h 0m 30s

- 4 Add  10 μL of 2x protein [M] 10 micromolar (μM) solution to each well containing the compounds to be tested previously dispensed onto the plate. 
- 5 Incubate the mix for  01:00:00 at  Room temperature to initiate the enzymatic reaction by the addition of  10 μL of 2x ([M] 20 micromolar (μM)) substrate solution using the plate reader injector.  1h
- 6 Read the fluorescence intensity at 350/460 nm every  00:00:30 for  02:00:00 in kinetic mode, which includes a shaking step of the plate before each measurement. 2h 0m 30s
- 7 Calculate the IC50 by plotting the initial velocity against various concentrations of tested inhibitors by using a four-parameter dose-response curve in Prism (v8.0) software.
- 8 Calculate the mean (μ) and the standard deviation (σ) of fluorescence intensity and then calculate the signal-to-background ratio and the Z' or Z-factor. (s: signal; c: control).

$$\text{Estimated Z-factor} = 1 - \frac{3(\hat{\sigma}_s + \hat{\sigma}_c)}{|\hat{\mu}_s - \hat{\mu}_c|}$$

Asset URL:

9



Asset URL: