



Dec 09, 2024

Fluorescence immunohistochemistry

DOI

dx.doi.org/10.17504/protocols.io.n92ldrmd8g5b/v1

Ashley Seifert¹

¹University of Kentucky



Ashley Seifert

University of Kentucky

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.n92ldrmd8g5b/v1>

Protocol Citation: Ashley Seifert 2024. Fluorescence immunohistochemistry. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.n92ldrmd8g5b/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: December 06, 2024



Last Modified: December 09, 2024

Protocol Integer ID: 114535

Keywords: ASAPCRN, fluorescence immunohistochemistry this protocol, fluorescence immunohistochemistry, procedure for immunohistochemistry, immunohistochemistry, seifert lab, fluorescence, embedded tissue section, tissue section, paraffin, lab

Disclaimer

Note that any protocol involving animals should be reviewed and approved by your Institutional Animal Care and Use Committee (IACUC) before use.

Abstract

This protocol describes the procedure for immunohistochemistry detected by fluorescence used in the Seifert Lab. It assumes a starting sample of a paraffin-embedded tissue section.

Protocol materials

⊗ Xylene **Fisher Scientific Catalog #X3P-1GAL**

⊗ 100% Ethanol

⊗ 100% Ethanol

⊗ Tris-buffered saline (TBS), 1x solution **Fisher Scientific Catalog #BP24721**

⊗ Hoechst 33342 **Cell Signaling Technology Catalog #4082**


⊗ Prolong Gold **Thermo Fisher Scientific Catalog #P36930**

Troubleshooting



Day 1




1 Deparaffinize slides to  100% Ethanol by a series of:

- 2 washes of  00:05:00 each with 2X

 Xylene **Fisher Scientific Catalog #X3P-1GAL**

- 2 washes of  00:02:00 each with  100% Ethanol

2 Rehydrate to ddH₂O by a series of:

-  00:03:00 wash with 90% Ethanol
-  00:01:00 wash with 70% Ethanol
- 2 washes of  00:01:00 each with ddH₂O


3 *(SKIP IF FROZEN)*




Perform retrieval, optimized for antigen and treatment

Note

For a first run, can try heat + citrate buffer pH 6, heat + DAKO high pH 9, proteinase K
*After heat, transfer directly to ddH₂O to cool before proceeding to first wash with TBS
**For matrix proteins, try 2' proteinase K

4 Wash for  00:05:00 with

 Tris-buffered saline (TBS), 1x solution **Fisher Scientific Catalog #BP24721**

5 Block for at least  00:30:00 at  Room temperature in  15 µL/mL appropriate serum in TBS

Note

The serum is typically whatever the secondary antibody was raised in; for multiple labeling, use a 50:50 mixture of the two sera

6 *(SKIP IF NOT AMPLIFYING SECONDARY ANTIBODY)*



Avidin/Biotin Block

- 6.1 Block with avidin for 00:15:00 at Room temperature
- 6.2 Wash for 00:05:00 in TBS
- 6.3 Block with biotin for 00:15:00 at Room temperature
- 6.4 Wash for 00:05:00 in TBS
- 7 Incubate with primary antibody or control in 15 $\mu\text{L}/\text{mL}$ serum in TBS Overnight at 4 $^{\circ}\text{C}$ in a humidity box

Note

The preferred control is IgG at the same concentration as the primary antibody, but otherwise just leave out the primary antibody










Day 2

- 8 Wash for 00:05:00 in TBS
- 9 Incubate with secondary antibody (biotinylated secondary to the species of the primary) at 1:400 in 15 $\mu\text{L}/\text{mL}$ serum in TBS for 00:30:00 at Room temperature

Note

1:400 is usually a good dilution, but this can be reduced if background is too high



- 10 Wash for  00:05:00 in TBS
- 11 Incubate with desired AlexaFluor-conjugated secondary antibody at 1:400 in TBS (without serum) for  00:30:00 at  Room temperature
- 12 Wash for  00:05:00 in TBS
- 13 Wash for  00:05:00 in ddH₂O
- 14 Incubate with  1 µg/mL
 Hoechst 33342 **Cell Signaling Technology Catalog #4082** for  00:05:00 at  Room temperature
- 15 Wash with ddH₂O
- 16 Stand slide at an angle and allow to drip-dry for a few minutes
- 17 Cover slip

Note

If using  Prolong Gold **Thermo Fisher Scientific Catalog #P36930** can cover slip when a bit wet