

Oct 01, 2019

FLUORESCENCE LOSS ASSAY

DOI

<https://dx.doi.org/10.17504/protocols.io.7mfhk3n>

Alexander Niederau¹, Despoina Trasanidou¹

¹Wageningen University

iGEM Wageningen 2019



Alexander Niederau

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.7mfhk3n>

Protocol Citation: Alexander Niederau, Despoina Trasanidou 2019. FLUORESCENCE LOSS ASSAY. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.7mfhk3n>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: In development

We are still developing and optimizing this protocol

Created: September 23, 2019

Last Modified: October 01, 2019

Protocol Integer ID: 28039

Keywords: fluorescence, loss



Day 1 (=Transformation → ONC in 96-well plate)

- 1 50ul of chemically competent *E. coli* DH10B_gfp cells were transformed (heat-shock: 42°C, 30sec) with 3ng plasmid (Cplasmid*2ul=3ng/ul*x → x=...ul=2ul undiluted plasmid + rest ul MQ) and recovered in 450 ul LB for 1h at 37°C.
(PC = *E. coli* DH10B_gfp + pACYC184)
(NC = *E. coli* DH10B + pACYC184)

Day 1 (=Transformation → ONC in 96-well plate)

- 2 2ul of recovered cells were inoculated in 198ul M9TG+Cam15. [1 Masterblock] (B = Blank, just medium)

Day 1 (=Transformation → ONC in 96-well plate)

- 3 Incubate at 37°C under shaking (900rpm) for 21-22 hours.

Day 2 (=Dilution 1/100 NO INDUCERS → Dilution 1/100 WITH INDUCERS → ONC in triplicates 96-well plates)

- 4 2ul of each preculture were inoculated in 198ul M9TG+Cam15. [1 Masterblock]
- 5 Shake the plate on the thermoblock (900rpm)
- 6 2ul of each dilution were inoculated in 198ul M9TG+Cam15+inducer. [4 Masterblocks=triplicates + control plate] Plate 1,2,3 (triplicates)
- 7 Incubate at 37°C under shaking (900rpm) for 21-22 hours.

Day 3 (=SPECTROPHOTOMETRY: 1/5 DILUTION → 100ul)

- 8 Mix at thermoblock the preculture → 40ul of each sample were inoculated in 160ul of 1xPBS [4 normal 96-well plates]
- 9 Transfer 100ul in 4 black plates with transparent bottom

**10** Procedure Details:

	Plate Type	96 WELL PLATE
	Eject plate on completion	
	Set Temperature	Setpoint 25°C
		Preheat before moving to next step
	Shake	Fast, 0:10 (MM:SS)
	Read	Cell Density
		Absorba nce Endpoint
		Full Plate
		Wavelen gths: 60 0
		Read Speed: Normal, Delay: 100 msec, Measure ments/D ata Point: 8
	Read	Fluoresc ence Endpoint
		Full Plate
		Filter Set 1



		Excitation: 395/20, 0, Emission: 508/20, 0
		Optics : Top, Gain: 50
		Filter Set 2
		Excitation: 395/20, 0, Emission: 508/20, 0
		Optics : Top, Gain: 75
		Filter Set 3
		Excitation: 395/20, 0, Emission: 508/20, 0
		Optics : Top, Gain: 100
		Read Speed: Normal, Delay: 100 msec, Measurements/Data Point: 10
		Read Height: 8 mm

