Oct 18, 2019

Fluorescence measurement

DOI

<https://dx.doi.org/10.17504/protocols.io.8f9htr6>

iGEM Dusseldorf¹

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Igem Dusseldorf

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS

DOI: <https://dx.doi.org/10.17504/protocols.io.8f9htr6>

Protocol Citation: iGEM Dusseldorf 2019. Fluorescence measurement. **protocols.io**

<https://dx.doi.org/10.17504/protocols.io.8f9htr6>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: October 18, 2019

Last Modified: October 18, 2019

Protocol Integer ID: 28897

Keywords: fluorescence, measurement

Materials

LB-medium

antibiotic resistance

Tergitol

24 well plate

96 well plate

- 1 Prepare 5 ml LB medium with antibiotic resistance and 0.5% ml Tergitol
- 2 Inoculate one colony in the prepared LB medium
- 3 Add 0.01-1 mM of the desired fatty acid to the culture medium
- 4 Inoculate the culture for at least 16 hours at 37 °C
- 5 Measure the absorption with a photometer at a wavelength 600 nm
- 6 Dilute the culture to an OD₆₀₀ of 0.05 with LB medium
- 7 Fill a 24 well plate with 800 µl of 3 biological with 3 technical replicats respectively like shown in the table below

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- 9 Inoculate the 24 well plate at 37°C for at least 16 hours

10 Transfer the samples from the 24 well plate on a 96 well plate like shown in the table below

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12 Measure the desired fluorescence at the appropriate wavelength and the OD₆₀₀ with a plate reader