

Sep 14, 2019 Version 1

Fluorescence microscope operation procedure V.1

DOI

dx.doi.org/10.17504/protocols.io.7cmhiu6

宏亮 董¹

¹Northeast Forest University

2019 iGEM NEFU_China

Tech. support email: shengyiyanwork@gmail.com



宏亮 董

OPEN  ACCESS



DOI: dx.doi.org/10.17504/protocols.io.7cmhiu6

Protocol Citation: 宏亮 董 2019. Fluorescence microscope operation procedure. **protocols.io**

<https://dx.doi.org/10.17504/protocols.io.7cmhiu6>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: In development

We are still developing and optimizing this protocol

Created: September 14, 2019

Last Modified: September 14, 2019

Protocol Integer ID: 27757

Keywords: Fluorescence microscope



- 1 Production of temporary fitting
 - 1.1 Add 15ul Bacterial solution to microslide at the clean workbench.
 - 1.2 Add cover slide on it.
- 2 Use of fluorescence microscopy
 - 2.1 Turn off the light in the room and turn on the microscope mercury lamp.
 - 2.2 Select the appropriate filter according to the luciferin labeled with the sample.
 - 2.3 Place the sample and find the right view. If you want to take a picture, please make sure that the camera is loaded with color film.
 - 2.4 Turn on the selfie device and select the manual file, usually within 0.5 to 10 seconds.
 - 2.5 End of use, turn off all power and record use.