

Aug 26, 2023

Fluorescence recovery after photobleaching (FRAP)

DOI

<https://dx.doi.org/10.17504/protocols.io.5qpobj29l4o/v1>

Xinbo Wang^{1,2}, Pietro De Camilli^{1,2}

¹1. Departments of Neuroscience and of Cell Biology, Howard Hughes Medical Institute, Program in Cellular Neuroscience, Neurodegeneration and Repair, Yale University School of Medicine, New Haven, Connecticut 06510, USA;

²2. Aligning Science Across Parkinson's (ASAP) Collaborative Research Network, Chevy Chase, MD, 20815



xinbo.wang

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.5qpobj29l4o/v1>

Protocol Citation: Xinbo Wang, Pietro De Camilli 2023. Fluorescence recovery after photobleaching (FRAP). **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.5qpobj29l4o/v1>

License: This is an open access protocol distributed under the terms of the **[Creative Commons Attribution License](#)**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited



Protocol status: Working

We use this protocol and it's working

Created: August 19, 2022

Last Modified: May 31, 2024

Protocol Integer ID: 68888

Keywords: Photobleaching, Fluorescence, LRRK2, ASAPCRN, fluorescence recovery, induced liposome tubule, liposome tubule, liposome, analysis of Irrk2

Abstract

This protocol details methods of the FRAP analysis of LRRK2-induced liposome tubules *in vitro*

Attachments



[juufbvk9p.docx](#)

18KB



Fluorescence recovery after photobleaching (FRAP)

40m 30s

1 Prepare LRRK2-liposome mixtures in a PCR tube with [M] 300 nanomolar (nM) GFP-LRKK2, [M] 20 micromolar (μ M) liposomes (labeled with trace amounts of rhodamine-PE) and [M] 1 millimolar (mM) GMPPNP.



2 Immediately deposit 6 μ L - 10 μ L samples of step 1 on a 35 mm glass bottom dish and incubate at 37 °C for 00:30:00 .

30m



Note

Note: Drop some buffer in the dish to prevent samples from drying out due to evaporation during incubation.

3 Perform FRAP experiments with a Spinning disk confocal (SDC) microscopy at Room temperature on a Nikon Ti-E inverted microscope using the Improvion UltraVIEW VoX system, with the settings as:



3.1 Acquire the time-lapse images at every 00:00:15 .

15s

3.2 Acquire three images before bleaching.

3.3 Bleach three ROIs with a 488 nm laser for 500 ms.

3.4 Acquire post-bleach images up to 00:10:00 at 00:00:15 intervals.

10m 15s