



Jan 10, 2026

Fluorescence recovery after photobleaching (FRAP) in neuro

DOI

<https://dx.doi.org/10.17504/protocols.io.kxygx4do4l8j/v1>

Akio Mori¹, Robert Edwards¹

¹University of California San Francisco



Akio Mori

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.kxygx4do4l8j/v1>

Protocol Citation: Akio Mori, Robert Edwards 2026. Fluorescence recovery after photobleaching (FRAP) in neuro. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.kxygx4do4l8j/v1>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: November 20, 2025



Last Modified: January 10, 2026

Protocol Integer ID: 233008

Keywords: FRAP, α -synuclein, synuclein in synaptic bouton, synuclein, fluorescence recovery, synaptic bouton, neuro this protocol, frap, cultured hippocampal neuron, neuro, neuron, recovery

Abstract

This protocol describes a fluorescence recovery after photobleaching (FRAP) assay used to measure the mobility of α -synuclein in synaptic boutons of cultured hippocampal neurons.

Materials

Reagents

- **Tyrode's buffer:** 119 mM NaCl, 2.5 mM KCl, 2 mM CaCl_2 , 2 mM MgCl_2 , 30 mM glucose, 25 mM HEPES, pH 7.4
- D-AP5 (ab120271, abcam)
- CNQX (ab120044, abcam)

Consumables

- Glass-bottom dishes (P35G-1.5-10-C, MatTek)

Equipment

- Nikon CSU-W1 spinning disk confocal microscope with Andor DU-888 EMCCD camera
- Rapp programmable illumination system
- 60 $\times$ oil objective (PlanApo λ , NA 1.40)
- Temperature and CO_2 controller (37°C, 5% CO_2)

Neuron preparation and expression of α -synuclein-GFP in hippocampal neurons

- 1 Culture synuclein triple-knockout hippocampal neurons on MatTek glass-bottom dishes.
- 2 At DIV7, transduce neurons with lentivirus encoding human α -synuclein-GFP (α -syn-GFP).
- 3 Maintain the neurons for 14–18 days after transduction.

Photobleaching

- 4 Calibration using a calibration solution (0.01 mg/ml fluorescein for photobleaching calibration). This step is crucial for obtaining consistent results. Ensure you use the same type of glass-bottom dish that will be used for the actual experiment. Optimize the laser focus to ensure a sharp photobleaching spot, thereby avoiding unintended bleaching of adjacent synaptic boutons.
- 5 Replace culture medium with Tyrode's buffer with D-AP5 (2-amino-5-phosphonopentanoic acid) and CNQX before imaging.
- 6 Place the dish on the microscope stage and maintain 37°C with 5% CO₂ throughout imaging.
 - Wait for 5-10 minutes until the temperature stabilizes to avoid focal drift during image acquisition.
- 7 Identify α -syn-GFP-positive synaptic boutons using the spinning disk confocal microscope. 
- 8 Select a single-pixel ROI within a synaptic bouton for bleaching. 
- 9 Photobleach the ROI using a 405 nm laser at 30 mW for 50 ms (no repeat). Immediately begin acquiring fluorescence recovery using 488 nm excitation (5% laser power, 200 ms exposure with an EM gain of 200) for 100 frames (no interval) over 20 seconds using the Andor DU-888 EMCCD camera and Rapp illumination system. 

Data analysis



- 10 Using Fiji (NIH), select the bleached synaptic bouton as an ROI and measure the 488 nm fluorescence intensity using the Time Series Analyzer V3 plugin. 
- 11 Normalize intensity values as follows: Pre-bleach fluorescence = 100%, Post-bleach fluorescence = 0%
Generate fluorescence recovery curves and fit the first 10 s of data with a one-phase exponential model in GraphPad Prism to calculate the recovery time constant (τ).
Determine the extent of fluorescence recovery at 20 s after photobleaching to quantify the recovery rate. 

Protocol references

Fortin, D.L.*et al.* Neural activity controls the synaptic accumulation of alpha-synuclein. *J Neurosci* **25**, 10913-10921 (2005).