

Aug 28, 2024

Fluorescent Immunolabelling for Alpha-Synuclein in neuronal primary culture (Testing PFF toxicity)

DOI

<https://dx.doi.org/10.17504/protocols.io.kxygxzr6wv8j/v1>

Rong Chen^{1,2}, Ted Dawson^{3,2}, Roberta Marongiu^{4,2}, Sabina Marciano⁵

¹Johns Hopkins University;

²Aligning Science Across Parkinson's (ASAP) Collaborative Research Network, Chevy Chase, MD 20815, USA;

³Johns Hopkins University School of Medicine;

⁴Department of Neurological Surgery, Weill Cornell Medical College, New York, NY 10065;

⁵Weill Cornell medicine



Jacquelyn Haytayan

Weill Cornell Medicine

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.kxygxzr6wv8j/v1>

Protocol Citation: Rong Chen, Ted Dawson, Roberta Marongiu, Sabina Marciano 2024. Fluorescent Immunolabelling for Alpha-Synuclein in neuronal primary culture (Testing PFF toxicity). **protocols.io**

<https://dx.doi.org/10.17504/protocols.io.kxygxzr6wv8j/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: August 23, 2022

Last Modified: September 23, 2024

Protocol Integer ID: 69056

Keywords: ASAPCRN, synuclein in neuronal primary culture, fluorescent immunolabelling on neuronal primary culture, synuclein, marker of the pff toxicity, testing pff toxicity, fluorescent immunolabelling, protocol for neuronal primary culture, pff toxicity, phosphorylated alpha, pff incubation, neuronal primary culture, protocol for pff preparation, pff preparation, vivo injection

Funders Acknowledgements:

Aligning Science Across Parkinson's

Grant ID: 020608

Abstract

This protocol is designed to perform fluorescent immunolabelling on neuronal primary culture after PFF incubation.


The labelling of phosphorylated Alpha-synuclein is considered as a marker of the PFF toxicity and should be performed to test the PFF before *in vivo* injection.

See protocol for PFF preparation

<https://www.protocols.io/view/production-of-alpha-synuclein-preformed-fibrils-pf-b39qr56>

and protocol for neuronal primary culture

<https://www.protocols.io/view/primary-culture-cortical-hippocampal-neurons-e15-1-cfiztkf6>



Fluorescent immunolabelling of phosphorylated Alpha Synuclein (green) in neuronal primary culture.

Image Attribution

image by Sabina Marciano



Materials

Reagents:

⊗ Purified anti- α -Synuclein Phospho (Ser129) Antibody **BioLegend Catalog #MMS-5091**

⊗ Anti-MAP2 antibody (ab5392) **Abcam Catalog #ab5392**

⊗ Donkey Anti-Mouse IgG H&L (Alexa Fluor® 488) **Abcam Catalog #ab150105**

⊗ Donkey anti-Chicken IgY (H L) Highly Cross Adsorbed Secondary Antibody Alexa Fluor™ 594 **Invitrogen - Thermo Fisher Catalog #A78951**

Solutions:

4% (wt/vol) Paraformaldehyde/4% (wt/vol) sucrose/PBS

For 50 ml in PBS, pipette 10 ml of 20% paraformaldehyde into ~40 ml of PBS, add 2 g of sucrose and mix it until it dissolves. Bring the volume to 50 ml with PBS. This solution works best if it is used at room temperature, although others report optimal results when it is warmed to 37degrees. This solution can be stored at 4 degrees celsius, protected from light, for 1 week.

PBS/BSA0.005%/TWEEN0.05%

Dissolve 1mg BSA in 10mL PBS. Use 10uL of this solution for 20mL of PBS. Add 10uL Tween 20. Prepare this fresh on the day of use. All the washing steps are done for 10min.



- 1 Aspirate the medium from coverslips (aspirate from only a few coverslips at a time).
- 2 Add 4% (wt/vol) PFA /4% (wt/vol) sucrose (1 ml for a 12-well plate).
- 3 Incubate the plate at Room temperature for 00:15:00 . 15m
- 4 Wash 3 times 00:10:00 with PBS/BSA/TWEEN solution. 10m
- 5 Proceed or store in the fridge in PBS, protected from light, for up to 1 week.
- 6 Set up humidifier chambers.
- 7 Incubate with 0.5% Triton in PBS/BSA/TWEEN for 00:05:00 . 5m
- 8 Wash 3 times 00:10:00 with PBS/BSA/TWEEN. 10m
- 9 Block in Donkey Serum 1:20 in PBS for 00:30:00 . 30m
- 10 During the blocking, prepare the primary antibody solution.
In PBS/BSA/TWEEN, add Biolegend psyn (MMS-5091) Mouse (1:1000) and Abcam MAP2 (ab5392) chicken Millipore (1:2000).
- 11 Incubate with primary antibody for 01:00:00 at Room temperature . 1h
- 12 Wash 3 times 00:10:00 with PBS/BSA/TWEEN. 10m
- 13 During the washes, prepare the secondary antibody solution:



In PBS/BSA/TWEEN, add Alexa Fluor 488 donkey anti-mouse (1:1000) + Alexa Fluor 594 donkey anti-chicken (1:1000).

- 14 Incubate in secondary antibody solution for  01:00:00 at  Room temperature in 1h the DARK.
- 15 Wash 3 times  00:10:00 with PBS/BSA/TWEEN. 10m
- 16 Incubate with DAPI 1:10,000 in PBS  00:10:00 at  Room temperature . 10m
- 17 Wash 3 times  00:10:00 with PBS. 10m
- 18 Mount the coverslips onto glass slides with Prolong Gold mounting medium.