

Aug 28, 2024

Fluorescent Immunolabelling for endogenous mouse LRRK2

DOI

<https://dx.doi.org/10.17504/protocols.io.261genqkog47/v1>

Roberta Marongiu^{1,2}, andrew.west west^{3,2}

¹Department of Neurological Surgery, Weill Cornell Medical College, New York, NY 10065;

²Aligning Science Across Parkinson's (ASAP) Collaborative Research Network, Chevy Chase, MD 20815, USA;

³Duke University



Jacquelyn Haytayan

Weill Cornell Medicine

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.261genqkog47/v1>

Protocol Citation: Roberta Marongiu, andrew.west west 2024. Fluorescent Immunolabelling for endogenous mouse LRRK2. protocols.io <https://dx.doi.org/10.17504/protocols.io.261genqkog47/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: August 19, 2022

Last Modified: September 23, 2024

Protocol Integer ID: 68930

Keywords: ASAPCRN, LRRK2, endogenous mouse Irrk2, endogenous mouse Irrk2 this protocol, endogenous Irrk2 in mouse, fluorescent immunolabelling, endogenous Irrk2, fluorescent, fixed brain tissue, brain tissue, mouse

Funders Acknowledgements:

Aligning Science Across Parkinson's

Grant ID: 020608

Abstract

This protocol is used to label endogenous LRRK2 in mouse.

It has been edited from West et al., 2014

Protocol optimized for fresh fixed brain tissue sectioned at 40 um.

Materials

 Primary Antibody anti-LRRK2 c41-2 **Abcam Catalog #ab133474**

 DAPI **Thermo Fisher Scientific Catalog #D1306**

 Donkey serum **Merck MilliporeSigma (Sigma-Aldrich) Catalog #S30-100ML**

Solutions:

Tris Buffer Solution TBS - 2L

- 50mM Tris pH 7.5 – 100 mL of 1M stock

- 150 mM NaCl – 17.53 g

- Distilled Water – 1.9 L



- 1 Select brain slices, from fresh, unfrozen, fixed tissue - slices are 40um thick
- 2 Wash 3 times for  00:10:00 in TBS 10m
- 3 Sections were "quenched" in 100% Methanol for  00:15:00 at  4 °C on shaker. 15m
- 4 Wash 3 times for  00:10:00 in TBS 10m
- 5 Incubate the slices with Sodium Citrate – 10 mM Sodium Citrate pH 6.0 w/0.05% Tween for  00:30:00 at  37 °C on shaker. 30m
- 6 Wash 3 times for  00:10:00 in TBS 10m
- 7 During the washes, prepare the blocking solution : 5% normal donkey serum (matched to secondary AB) in TBS with 0.3% Triton
- 8 Incubate the slides with blocking solution for  01:00:00 at  4 °C on shaker. 1h
- 9 Wash 3 times for  00:10:00 in TBS 10m
- 10 During the washes, prepare the primary antibody solution : 5% normal donkey serum in TBS with NO TRITON + Primary Antibody anti-LRRK2 (1:500 or 1:1000)
- 11 Incubate at  4 °C on shaker for  24:00:00  48:00:00 3d
- 12 Wash 3 times for  00:10:00 in TBS 10m
- 13 During the washes, prepare the secondary antibody solution : 5% normal donkey serum in TBS with NO TRITON + Alexa Fluor secondary AB (1:1000)



-
- 14 Incubate at 4 °C on shaker for 18:00:00 18h
- 15 Wash 2 times for 00:05:00 in TBS 5m
- 16 During the washes, prepare the DAPI solution : DAPI (1mg/ml): 1:10,000, in TBS
- 17 Incubate the slides with the DAPI solution for 00:10:00 10m
- 18 Wash 3 times for 00:10:00 in TBS 10m
- 19 Mount sections on superfrost + slides and coverslip with Invitrogen Prolong Antifade reagent.