

Sep 08, 2020

## Fluorescent-Reporter Based Assay

DOI

[dx.doi.org/10.17504/protocols.io.bk5rky56](https://dx.doi.org/10.17504/protocols.io.bk5rky56)

Piyush Jain<sup>1</sup>, Long T Nguyen<sup>1</sup>, Santosh Rananaware<sup>1</sup>

<sup>1</sup>University of Florida



s.rananaware

OPEN  ACCESS



DOI: [dx.doi.org/10.17504/protocols.io.bk5rky56](https://dx.doi.org/10.17504/protocols.io.bk5rky56)

**Protocol Citation:** Piyush Jain, Long T Nguyen, Santosh Rananaware 2020. Fluorescent-Reporter Based Assay. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.bk5rky56>

**License:** This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** In development

**We are still developing and optimizing this protocol**

**Created:** September 08, 2020

**Last Modified:** September 08, 2020

**Protocol Integer ID:** 41873

**Keywords:** CRISPR, SARS-CoV-2, COVID-19 Diagnostic,

## Abstract

The CRISPR-Enhance SARS-CoV-2 detection kit has been designed to detect fragments of the Nucleocapsid ("N") gene and Envelope gene (E) of SARS-CoV-2. An included third target is the human RNase P POP7 gene ("RP") which serves as a control for the extraction of the clinical sample in the absence of a positive SARS-CoV-2 result. Amplification can be performed using a heat block, and CRISPR complex activation and reporter cleavage can be run in a standard microplate reader capable of fluorescence detection. The entire reaction from RT-LAMP amplification to CRISPR-based detection of the target analytes can be performed in approximately one hour.


The CRISPR-Enhance kit comprises of two steps. Step one is a reverse transcriptase loop-mediated amplification (RT-LAMP) where targeted SARS-CoV-2 genomic RNA is reverse transcribed to DNA, and this DNA is amplified by a strand-displacing DNA polymerase. Step two is the transcription of the amplified DNA to activate the collateral cleavage activity of a CRISPR complex programmed to the target RNA sequence. This collateral activity results in cleavage of nucleic acid reporters, resulting in a fluorescent readout detected by a plate reader.


## Guidelines

All procedures should be performed in a BSL2 laboratory, and specimens should be handled within a Biological Safety Cabinet. All necessary safety precautions should be taken according to the Laboratory guidelines. Precautions must also be taken to prevent cross-contamination of samples.


## Materials


### MATERIALS

 QuickExtract™ RNA Extraction Kit **Lucigen Catalog #QER090150**

 WarmStart®Colorimetric LAMP 2X Master Mix with UDG (Cat.No. M1804S) **New England Biolabs Catalog #M1804S**


### STEP MATERIALS


 QuickExtract™ RNA Extraction Kit **Lucigen Catalog #QER090150**


 WarmStart®Colorimetric LAMP 2X Master Mix with UDG (Cat.No. M1804S) **New England Biolabs Catalog #M1804S**





## Protocol materials

 QuickExtract™ RNA Extraction Kit **Lucigen Catalog #QER090150**

 WarmStart®Colorimetric LAMP 2X Master Mix with UDG (Cat.No. M1804S) **New England Biolabs Catalog #M1804S**


 QuickExtract™ RNA Extraction Kit **Lucigen Catalog #QER090150**

 WarmStart®Colorimetric LAMP 2X Master Mix with UDG (Cat.No. M1804S) **New England Biolabs Catalog #M1804S**

 WarmStart®Colorimetric LAMP 2X Master Mix with UDG (Cat.No. M1804S) **New England Biolabs Catalog #M1804S**

 QuickExtract™ RNA Extraction Kit **Lucigen Catalog #QER090150**

## Safety warnings

-  1. Handle all infectious samples with appropriate CDC approved methods
- 2. Wear appropriate PPE such as lab coats, gloves, N95 respirators, safety goggles etc when handling infectious samples
- 3. Discard all biohazard waste appropriately
- 4. Clean all work surfaces with bleach and IPA after use

## Nucleic Acid Extraction

18m

1


18m


- The CRISPR-Enhance SARS-CoV-2 detection kit uses QuickExtract™ RNA Extraction Kit

 QuickExtract™ RNA Extraction Kit **Lucigen Catalog #QER090150**

- Add  10 µL of patient sample to  10 µL of pre-aliquoted QuickExtract solution.
- Heat the above mixture at  65 °C for  00:15:00 followed by  98 °C for  00:03:00 .

## RT-LAMP Master Mix Preparation

- 2
- Label a new  1.5 mL microcentrifuge tube for each target ( N, E and RP) and prepare a RT-LAMP Master Mix consisting of the WarmStart® Colorimetric LAMP 2X Master Mix with UDG.

 WarmStart® Colorimetric LAMP 2X Master Mix with UDG (Cat.No. M1804S) **New England Biolabs Catalog #M1804S**


and the appropriate 10x Primer Mix using the recipe in Table 1 below. Make enough of each master mix for all samples to be tested and the necessary controls for each run.






Reagent Name	Volume per reaction	Total Volume
WarmStart®Colorimetric LAMP 2X Master Mix with UDG	12.5 µL	12.5 µL x (N+1)
10x Primer Mix (N, E or RP)	2.5 µL	2.5 µL x (N+1)
RNase-Free Water	5 µL	5 µL x (N+1)
<b>Total Volume</b>	<b>20 µL</b>	<b>20 µL x (N+1)</b>

**Table 1: Target Specific RT-LAMP Master Mix Recipe**

N = number of extracted samples plus number of controls. Prepare enough for 1 extra (N + 1) sample to allow for overage during reaction set-up.

## RT-LAMP Amplification

- 3
- Label a strip tube (  0.2 mL ) with the target name (e.g. N) and strip number corresponding to each sample.







- Add  20  $\mu\text{L}$  of the RT-LAMP Master Mix from the previous step into one well for each sample and control to be amplified. Repeat for the remaining 2 targets using a new strip for each target (e.g. E or RP)
- Add  5  $\mu\text{L}$  of extracted RNA in each respective strip tube containing the RT-LAMP Master mix. Vortex the strip tube for  00:00:03 and spin down for  00:00:03 in microcentrifuge with a  0.2 mL tube adaptor.

Reagent	Volume per reaction
RT-LAMP Master Mix	20 $\mu\text{L}$
RNA Sample or Controls	5 $\mu\text{L}$
<b>Total Volume</b>	<b>25 <math>\mu\text{L}</math></b>

**Table 2: RT-LAMP Assay Components and reaction volume**

- Heat the mixture at  65  $^{\circ}\text{C}$  for  00:40:00

## CRISPR-Cas Reaction Preparation


- 4
  - Preheat a fluorescence microplate reader to  37  $^{\circ}\text{C}$ .
  - For each target tested label a  1.5 mL tube with the target name (e.g. N, E or RNASE-P) and "Cas Mix". Prepare a CRISPR Cas Master Mix using the following recipe in Table 3 below, scaling as required for the number of assays to be run (one Cas assay for every RT-LAMP reaction).
  - Incubate the mixture at  37  $^{\circ}\text{C}$  for  00:15:00
  - Pulse vortex for  00:00:03 and spin down for  00:00:03 in a microcentrifuge after all components are added.

Reagent Name	Volume per Reaction	Volume Total
NEB 2.1 Buffer	1.2 $\mu\text{L}$	1.2 $\mu\text{L}$ x (N+1)
3 $\mu\text{M}$ crRNA (N or E or RP)	0.8 $\mu\text{L}$	0.8 $\mu\text{L}$ x

		(N+1)
3 $\mu$ M IbCas12a	0.4 $\mu$ L	0.4 $\mu$ L x (N+1)
RNAse-Free Water	9.6 $\mu$ L	9.6 $\mu$ L x (N+1)
Total Volume	12 $\mu$ L	12 $\mu$ L x (N+1)

Table 3: Target CRISPR Cas Master Mix Recipe

*N* = number of extracted samples plus number of controls. Prepare enough for 1 extra (*N* + 1) sample to allow for overage during reaction set-up.

- In a separate  1.5 mL microcentrifuge tube prepare the fluorescence reporter as per Table 4

Reagent Name	Volume per Reaction	Volume Total
FAM-FQ	0.2 $\mu$ L	0.2 $\mu$ L x ( <i>N</i> *+1)
RNAse-Free Water	25.8 $\mu$ L	25.8 $\mu$ L x ( <i>N</i> *+1)
Total Volume	12 $\mu$ L	12 $\mu$ L x ( <i>N</i> *+1)

Table 4: Fluorescence Reporter Mix Recipe

*N*\* = number of extracted samples multiplied by the total number of genes plus number of controls. Prepare enough for 1 extra (*N*\* + 1) sample to allow for overage during reaction set-up

## CRISPR-Cas Detection



- 5
  - Add 26  $\mu\text{L}$  of the Fluorescent reporter Mix made in step 4 to each well of a 384 well-plate corresponding to the number of samples and controls for every gene (N, E or RP).
  - Add 2  $\mu\text{L}$  of the RT-LAMP product from step 3 to each well containing the fluorescent reporter.
  - Add 12  $\mu\text{L}$  of the CRISPR-Cas master mix to the wells containing the corresponding RT-Lamp product and fluorescent reporter.
  - Seal the plate with an Optical seal
  - Open the plate reader software to create a read procedure. Set temperature to 37°C
  - Select "Kinetic" run reading with a total read time of 30 min, and data collection intervals at 2.5 mins
  - Save experiment in a designated place with an appropriate unique name
  - When plate loader extends, load plate. Ensure plate is loaded in correct orientation. Read the data.