

Dec 06, 2022

🌐 Fluorescently labeled polyamine uptake (via Flow Cytometry)

DOI

<https://dx.doi.org/10.17504/protocols.io.n92ldp8qxl5b/v1>



Marine Houdou¹, Peter Vangheluwe¹

¹KU Leuven



Marine Houdou

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.n92ldp8qxl5b/v1>

External link: <https://doi.org/10.3390/biom13020337>

Protocol Citation: Marine Houdou, Peter Vangheluwe 2022. Fluorescently labeled polyamine uptake (via Flow Cytometry). protocols.io <https://dx.doi.org/10.17504/protocols.io.n92ldp8qxl5b/v1>

**Manuscript citation:**

Houdou M, Jacobs N, Coene J, Azfar M, Vanhoutte R, Haute CVd, Eggermont J, Daniëls V, Verhelst SHL, Vangheluwe P, Novel Green Fluorescent Polyamines to Analyze ATP13A2 and ATP13A3 Activity in the Mammalian Polyamine Transport System. Biomolecules 13(2). doi: [10.3390/biom13020337](https://doi.org/10.3390/biom13020337)

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: December 02, 2022

Last Modified: May 31, 2024

Protocol Integer ID: 73483

Keywords: Fluorescently labeled polyamines, ASAPCRN, labeled polyamine uptake, polyamine uptake, labeled polyamine, acquisition via flow cytometry, flow cytometry, fluorescence, cytometry, uptake capacity

Disclaimer

Protocol Particulars Video

The video below is a supplement with extra context and tips, as part of the Aligning Science Across Parkinson's (ASAP) Protocol Particulars video interview series, featuring conversations with protocol authors.

https://www.youtube.com/embed/c0GkyegZaG8?si=ZnQ8VCKG_nkCvcng

Abstract

Assess polyamine uptake capacities of a specific cell line after incubation with fluorescently labeled polyamines and mean fluorescence intensity acquisition *via* flow cytometry.

Guidelines

Make sure to stay protect from the light when manipulating fluorescent materials.



Materials

0.25% Trypsin-EDTA: Gibco, 25200056

Albumin Fraction V: Carl Roth, 8076.4

Dulbecco's Phosphate Buffered Saline modified without calcium chloride and magnesium chloride (DPBS (-/-)): Gibco, D8537

TrypLE Express Enzyme: Gibco, 12604021

Versene Solution: Gibco, 15040

Before start

- Prepare appropriate dilution of fluorescently labeled polyamine in cell culture medium.
- Prepare Flow Cytometry (FC) buffer made of 1% Albumin Fraction V in DPBS (-/-).
- Prepare Eppendorf tubes and FC tubes by labeling them and keeping them on ice.



- 1 Seed cells in 12 well-plate that they reach 70%-80% confluency the day of the experiment.

Note

This protocol is made for a 12 well-plate format. Adaptations need to be done to scale up or down.

- 2 The day of the experiment, remove cell culture medium and add fluorescently labeled polyamines to the cells in a final volume of  500 μ L / well .

Note

Concentration and incubation time may vary depending on the cell type. It's best to try different incubation times and doses to better appreciate the kinetics.

Note

Stay protect from the light.

- 3 Incubate the desired time at 37°C, 5% CO₂, protected from the light.
- 4 Harvest cells and prepare samples for Flow Cytometry (FC).
 - 4.1 Discard medium.
 - 4.2 Wash with  500 μ L /well Versene or DPBS (-/-).
 - 4.3 Discard Versene or DPBS (-/-).



- 4.4 Add  200 μ L /well 0.25% Trypsin or TrypLE. Incubate  00:05:00 at RT protect 5m
from the light.
- 4.5 Add  500 μ L /well 1%-Albumin Fraction V in DPBS (-/-), collect the cells and transfer
in Eppendorf tubes on ice.
- 4.6 Centrifuge for  00:05:00 at 400 g and 4°C. 5m
- 4.7 Discard supernatants and resuspend cell pellets in  250-500 μ L 1%-Albumin
Fraction V in DPBS (-/-), depending on the size of the cell pellet.
- 4.8 Filter cell suspension through Nylon filter into FC tubes.
- 5 Keep on ice until acquisition at the flow cytometer. Record 10,000 live events per sample.