

Aug 07, 2023

FM1-43 endocytic uptake assay in HIPSC derived neurons

DOI

<https://dx.doi.org/10.17504/protocols.io.3byl4qj52vo5/v1>

Sakthi Kumar¹

¹Waisman Center, University of Wisconsin-Madison, Madison, WI

Daniel's workspace



Daniel El Kodsi

Brigham and Women's Hospital and Harvard Medical School

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.3byl4qj52vo5/v1>

Protocol Citation: Sakthi Kumar 2023. FM1-43 endocytic uptake assay in HIPSC derived neurons. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.3byl4qj52vo5/v1>

License: This is an open access protocol distributed under the terms of the **[Creative Commons Attribution License](#)**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited



Protocol status: Working

We use this protocol and it's working

Created: August 07, 2023

Last Modified: May 31, 2024

Protocol Integer ID: 86065

Keywords: ASAPCRN, Endocytosis, FM1-43, neurons fm1, internalized within the recycled synaptic vesicle, derived neurons fm1, recycled synaptic vesicle, increase in the fluorescence intensity, endocytic uptake assay in hipsc, lipophilic styryl dye, fluorescence intensity, endocytic uptake assay, exocytosis, decrease in the fm1, fluorescent, releasing neurotransmitter, nontoxic to the cell, neuron, plasma membrane, dye, cell, fm1, containing vesicle, membrane, derived neuron, nerve terminal

Funders Acknowledgements:

MJFF

Grant ID: 000301

Abstract

FM1-43 is a lipophilic styryl dye, nontoxic to the cells, virtually nonfluorescent in aqueous media. The dye will become fluorescent once it binds to the plasma membrane. In neurons that are actively releasing neurotransmitters, these dyes become internalized within the recycled synaptic vesicles, and the nerve terminals become brightly stained. This shows an increase in the fluorescence intensity of FM1-43. When the cells undergo another stimulation (through KCL or electrical) the FM1-43 dye containing vesicles will be released through exocytosis. Thus, resulting in a decrease in the FM1-43 fluorescence intensity.



Reagents needed:

- 1 FM1-43 dye (Thermo fisher scientific, catalog no: T3163).
This dye is available in both live and fixable versions. Also, available in different colors like FM4-64

Mature neurons- hipsc derived Dopaminergic, Cortical or GABAergic neurons:

3w

- 2
 1. mDA, CTX or GABAergic neuron progenitors were plated in 12mm coverslips, precoated with PDL+ Laminin.
 2. After  504:00:00 post plating, neurons were subjected to FM1-43 dye uptake.
Note: We have not tried in early immature neurons; this experiment was also performed in more mature neurons like 40 days post plating.
 3. Some cell types may mature differentially. So, the user must validate the protocol based on their cell type and other experimental conditions.

3w

Required buffer:

- 3
 1. ACSF
 2. ACSF- containing high KCL (50mM) for stimulation.
 3. ACSF- containing no Calcium to prevent spontaneous release.

Normal ACSF composition:

Standard External (pH 7.4 NaOH)

		milliMolar	Gms/500ml
	NaCl	140	4.090 gms
	KCl	4	0.150
	CaCl ₂	2	0.150
	MgCl ₂	2	0.204
	TES(HEPES)	10(10)	1.146(1.191)
	Glucose	5	0.440

Spontaneous FM1-43 uptake:

10m



- 4
1. Transfer the coverslips to normal ACSF.
 2. Add 10uM of FM1 43 dye to the cells and wait for  00:10:00 at  Room temperature .
 3. Remove the dye and add ACSF without calcium and wash for 3 times. At this step, the user may want to add antagonists for AMPA/ NMDA or Na⁺ channel to inhibit spontaneous release.
 4. After 3 washes transfer the coverslips to the imaging chamber (for live imaging) or add PFA for fixation.

10m

PFA fixation:

25m

- 5
1. After adding 4% PFA wait for  00:15:00 for fixation.
 2. Then wash 3 times with PBS.
 3. Add Hoechst and incubate for  00:10:00 .
 4. Then wash 3 times with PBS and mount the coverslips on slides with Fluoromount G.
 5. Image the same day or keep the slides at  4 °C if imaging the next day.

25m

Stimulation induced FM1-43 dye uptake:

11m

- 6
1. Transfer the coverslips to normal ACSF.
 2. Add ACSF containing high KCL (50mM) for  00:01:00 to stimulate the cells.
 3. Then Add 10um of FM1 43 dye to the cells and wait for  00:10:00 at  Room temperature to measure the post stimulus induced endocytosis.
 4. Remove the dye and add ACSF without calcium and wash for 3 times. At this step, the user may want to add antagonists for AMPA/ NMDA or Na⁺ channel to inhibit spontaneous release.
 5. After 3 washes transfer the coverslips to the imaging chamber (for live imaging) or add PFA for fixation.
 6. If fixing the cells, follow the same steps as mentioned above in PFA fixation.

11m

Quantification:



- 7 Use ImageJ based program to quantify FM1-43 fluorescence intensity.