



Mar 15, 2024

🌐 Forebrain neural progenitor cells (fbNPCs) differentiation from hiPSCs (dual SMAD inhibition)

DOI

<https://dx.doi.org/10.17504/protocols.io.kxygx9ozzg8j/v1>

anita.adami¹

¹Laboratory of Molecular Neurogenetics, Department of Experimental Medical Science, Wallenberg Neuroscience Center and Lund Stem Cell Center, BMC A11, Lund University, 221 84 Lund, Sweden.



Anita Adami

Lund University

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.kxygx9ozzg8j/v1>

Protocol Citation: anita.adami 2024. Forebrain neural progenitor cells (fbNPCs) differentiation from hiPSCs (dual SMAD inhibition). **protocols.io** <https://dx.doi.org/10.17504/protocols.io.kxygx9ozzg8j/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: January 31, 2023

Last Modified: March 15, 2024

Protocol Integer ID: 76132

Keywords: dual smad inhibition, neural progenitor cell, using dual smad inhibition, differentiation from hipsc, hipsc, cell

Funders Acknowledgements:

Aligning Science Across Parkinson's through the Michael J. Fox Foundation for Parkinson's Research

Grant ID: ASAP-000520

Swedish Research Council

Grant ID: 2018-02694

Swedish Brain Foundation

Grant ID: FO2019-0098

Cancerfonden

Grant ID: 190326

Barncancerfonden

Grant ID: PR2017-0053

NIHR Cambridge Biomedical Research Centre

Grant ID: NIHR203312

Swedish Society for Medical Research

Grant ID: S19-0100

National Institutes of Health

Grant ID: HG002385

Swedish Research Council

Grant ID: 2021-03494

Swedish Research Council

Grant ID: 2020-01660

Abstract

This protocol described how to differentiate hiPSCs to forebrain neural progenitor cells using dual SMAD inhibition

hiPSCs dissociation

15m

- 1 Human induced pluripotent stem cells (hiPSCs) are rinsed once with DPBS (GIBCO) and dissociated when 70-90% confluent with 0.5 mM EDTA (75 ml/cm²; GIBCO) at  37 °C for  00:07:00
- 2 Dissociated cells are then carefully washed away from the wells and collected in **wash medium** ( 9.5 mL DMEM/F-12 (31330-038; GIBCO) and  0.5 mL knockout serum replacement (GIBCO)). 
- 3 The cells are then centrifuged at  400 x g for  00:05:00 and resuspended in iPS brew medium (StemMACS iPS-Brew XF and 0.5% penicillin/streptomycin (GIBCO)). 

hiPSCs differentiation into fbNPCs via dual SMAD inhibition

2w

- 4 The protocol is based on the one described in [Nolbrant et al., 2017](#) and [Grassi et al., 2019](#).

Day 0 of differentiation. The resuspended cells are then plated on LN111-coated (1.14µg/cm²; Biolamina) coated Nunc multidishes at a density of 10000 cells/cm² and grown in **N2 medium** (1:1 DMEM/F-12 (21331020; GIBCO) and Neurobasal (21103049; GIBCO) supplemented with 1% N2 (GIBCO), 2 mM L-glutamine (GIBCO), and 0.2% penicillin/streptomycin). 10 µM of Y27632 (Rock inhibitor; Miltenyi) are also added to the cells when plating. Finally, 10 µM SB431542 (Axon) and 100 ng/ml noggin (Miltenyi) for dual SMAD inhibition are supplemented to the media.
- 5 The media, supplemented with dual SMAD inhibitors, is changed every 2-3 days until *day 9 of differentiation*. The amount of media is slightly increased at every feeding as the cells' confluency increases.
- 6 *Day 9 of differentiation.* The cells are fed with N2 media without dual SMAD inhibitors.
- 7 *Day 11 of differentiation.* The cells are dissociated, centrifuged, and resuspended as described above (hiPSCs dissociation SECTION). They are then replated on LN111-coated Nunc multidishes at a density of 800000 cells/cm² in **B27 medium** (Neurobasal supplemented with 1% B27 without vitamin A (GIBCO), 2 mM L-glutamine and 0.2% penicillin/streptomycin Y27632 (10 µM), BDNF (20 ng/ml; R&D), and L-



ascorbic acid (0.2 mM; Sigma)). 10 μ M of Y27632 are also added. The cells are kept in the same medium until *day 14 of differentiation*, when they are harvested for downstream analyses (e.g. bulk RNA sequencing).