



Nov 01, 2022

Foreskin Tissue DNA Extraction

DOI

dx.doi.org/10.17504/protocols.io.4r3l2774jg1y/v1

Brandon Maust^{1,2}

¹Center for Global Infectious Disease Research, Seattle Children's Research Institute;

²Division of Infections Disease, Department of Pediatrics, University of Washington School of Medicine



Brandon Maust

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.4r3l2774jg1y/v1>

Protocol Citation: Brandon Maust 2022. Foreskin Tissue DNA Extraction. **protocols.io**

<https://dx.doi.org/10.17504/protocols.io.4r3l2774jg1y/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: October 31, 2022

Last Modified: November 01, 2022

Protocol Integer ID: 72091

Keywords: foreskin tissue dna extraction, total dna from foreskin tissue, foreskin tissue, dna extraction, total dna, dna, extraction

Disclaimer

DISCLAIMER – FOR INFORMATIONAL PURPOSES ONLY; USE AT YOUR OWN RISK

The protocol content here is for informational purposes only and does not constitute legal, medical, clinical, or safety advice, or otherwise; content added to [protocols.io](https://www.protocols.io) is not peer reviewed and may not have undergone a formal approval of any kind. Information presented in this protocol should not substitute for independent professional judgment, advice, diagnosis, or treatment. Any action you take or refrain from taking using or relying upon the information presented here is strictly at your own risk. You agree that neither the Company nor any of the authors, contributors, administrators, or anyone else associated with [protocols.io](https://www.protocols.io), can be held responsible for your use of the information contained in or linked to this protocol or any of our Sites/Apps and Services.

Abstract

Extracts total DNA from foreskin tissue, as used in Maust et al. 2022

Materials

Supplies

2 ml screw top microcentrifuge tubes
3mm steel beads (1 per sample)
scalpels
plate sealing film
proteinase K (20 mg/ml)
Absolute ethanol
PowerSoil Pro plates (Cat no 19311)
DNeasy 96 Powersoil Pro QIAcube HT kit (Cat no 47021)

Equipment

scale
heat block
Qiagen TissueLyser II
Qiagen QIAcube HT

Troubleshooting



Preparation

- 1 Add steel bead to 1.5 ml tube
- 2 Add 800 μ L CD1 to tube
- 3 Weigh tube (with bead and solution)

Tissue

- 4 Dissect approximately 25 mg (3 mm³) piece of tissue
Fragment tissue with scalpel
Add tissue to tube (with bead and solution)
- 5 Weigh tube and calculate net tissue weight.

Pre-robot processing

- 6 Process tubes on TissueLyzer II, 2 min at 30 Hz
Rotate block 180° and repeat
- 7 Centrifuge tubes briefly to collapse foam
- 8 Spin PowerBead Pro plate to ensure beads are settled at the bottom
- 9 Add contents of each specimen tube (except steel bead) to a PowerBead plate well
- 10 Add 5 μ L of 20 mg/ml proteinase K to each well
- 11 Seal plate with film
Vortex briefly to mix



- 12 Incubate plate at 65° C for 60 min or until tissue is mostly digested
- 13 Process plate on TissueLyser II, 5 min at 25 Hz
Rotate plate 180° and repeat
- 14 Centrifuge plate at 3000 x *g* for 7.5 min
- 15 Transfer approximately 350 µL supernatant from each well to fresh S-block using well-vator
- 16 add 300 µL of solution CD2 to each well and mix thoroughly by pipetting
- 17 Seal the plate with film
Centrifuge at 3000 x *g* for 7.5 min at room temperature
- 18 Avoiding the pellet, transfer 500 µL of supernatant to fresh S-block using well-vator

QIAcube

- 19 Follow DNeasy 96 PowerSoil Pro Protocol for QIAcube HT (page 17 onward)
Elute in maximum volume (120 µL)
Include vacuum performance check
- 20 Seal plate and freeze at -20°C