

Aug 23, 2025

Free Floating DAB Staining

 In 1 collection

DOI

<https://dx.doi.org/10.17504/protocols.io.q26g7pxqqgwz/v1>

Ashley Harms¹, Jhodi Webster¹

¹University of Alabama at Birmingham



Jhodi Webster

University of Alabama at Birmingham

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.q26g7pxqqgwz/v1>

Protocol Citation: Ashley Harms, Jhodi Webster 2025. Free Floating DAB Staining. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.q26g7pxqqgwz/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working



Created: October 14, 2023

Last Modified: August 23, 2025

Protocol Integer ID: 89279

Keywords: ASAPCRN, immunohistochemical staining, floating tissue section, floating dab, antibody penetration by processing section, permanent chromogenic precipitate for the microscopic visualization, antibody penetration, antigen distribution, analysis of antigen distribution, tissue section, dab, tissue integrity, permanent chromogenic precipitate, microscopic visualization

Abstract

This protocol describes the procedure for performing 3,3'-Diaminobenzidine (DAB) immunohistochemical staining on free-floating tissue sections. The method is optimized to preserve tissue integrity and enhance antibody penetration by processing sections without being mounted to a slide. It results in a stable, permanent chromogenic precipitate for the microscopic visualization and analysis of antigen distribution.

Safety warnings

- ⚠ Be sure to do DAB developing in a fume hood and dispose of all apparatus that touched DAB solution in bucket of bleach and water.



DAY 1

- 1 Wash sections 3×5min in TBS
- 2 Quench sections with solution for 5min at RT
Quenching solution (10mL – 3% H₂O₂ in 1XTBS):
4.5mL TBS
4.5mL MeOH
1mL 30% H₂O₂
- 3 Wash 2×5min TBS
- 4 Put sections in antigen retrieval for 30mins at 37C with agitation
Antigen retrieval in TBS:
10mM sodium citrate
0.05% Tween-20
- 5 Wash 3×5min in TBS
- 6 Block in 5% serum + TBST for at least 1hr @ room temp
- 7 24-well plate (500μL/well): dilute primary Ab in TBST + 1% serum (of secondary host animal).
- 8 Incubate overnight at 4C

Day 2

- 9 Wash sections 3×10min in TBST
- 10 Incubate sections with diluted biotinylated secondary Ab [(1:1000) in TBST + 1% serum] for 2hrs at RT in the dark
- 11 Wash sections 3×10min in TBS in the dark



- 12 Dilute ABC solution (1:5 in TBS) and incubate sections for 30min at RT (using 2mL/well)
- 13 Wash 3×10mins in TBS
- 14 Follow DAB kit instructions (work in the hood, set up a bucket with bleach and water)
- 15 Prepare substrate working solution. To every 5mL of diH₂O, add:
 - 16 2 drops of Buffer Stock Solution
 - 17 4 drops of DAB Stock Solution
 - 18 2 drops of Hydrogen Peroxide Solution
 - 19 If gray-black reaction is desired, add 2 drops of Nickel Solution
- 20 Move samples to substrate solutions and develop until sufficient colour change is evident; around 3-6mins
- 21 Immediately wash 3×5min in TBS
- 22 Mount sections on slides and dry for approximately 30min to 1hr or just until sections are adhered to slide.
- 23 Dehydrate:
 - 70% ethanol 2×3min
 - 90% ethanol 2×3min
 - 95% ethanol 2×3min
 - 100% ethanol 2×3min
 - Citrisolve 2×3min



24 Using Permout, cover with coverslip and dry several hours up to 2 days @RT.