

Jul 23, 2020

Freezing cancer cell lines

DOI

dx.doi.org/10.17504/protocols.io.bgtyjwpw

Emily Souster¹, Verity Goodwin¹, Adam Jackson¹, Charlotte Beaver¹, Rizwan Ansari¹, Fiona Behan¹, Mathew Garnett¹

¹Wellcome Sanger Institute

Cellular Generation and ...



Emily Souster

OPEN  ACCESS



DOI: dx.doi.org/10.17504/protocols.io.bgtyjwpw

Protocol Citation: Emily Souster, Verity Goodwin, Adam Jackson, Charlotte Beaver, Rizwan Ansari, Fiona Behan, Mathew Garnett 2020. Freezing cancer cell lines. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.bgtyjwpw>

License: This is an open access protocol distributed under the terms of the **[Creative Commons Attribution License](#)**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: May 25, 2020


Last Modified: July 23, 2020

Protocol Integer ID: 37464

Abstract

This protocol outlines routine banking of cancer cell lines and Ca9 transduced cancer lines.

Process diagram:



Guidelines

- As a guideline, we usually bank 5 cryovials from a 70% confluent T150 flask, each containing 1ml cell suspension.



Materials

MATERIALS

✕ Falcon™ 15mL Conical Centrifuge Tubes **Fisher Scientific Catalog #14-959-53A**

✕ TrypLE®; Express Enzyme (1X), no phenol red **Thermo Fisher Catalog #12604021**

✕ Nunc®; Biobanking and Cell Culture Cryogenic Tubes, 1.8mL, 48mm, external thread, printed **Thermo Fisher Catalog #375418**

✕ DMSO **Sigma Aldrich Catalog #D2650**

✕ DPBS Invitrogen - **Thermo Fisher Catalog #14190**

Select an appropriate culture media for your cell line. Common culture medias used for cancer cell lines are serum supplemented Advanced DMEM F-12 or RPMI, in the presence of pen-strep.

Equipment

Light Microscope

Microbiological Safety Cabinet (MSC)

Pipette Boy

Stripettes

Pipettes and tips

🔥 37 °C , 5% CO₂ incubator

Centrifuge

CoolCell or appropriate freezing container


-80C freezer

Liquid Nitrogen storage

Before start



- Pre-warm complete culture media to room-temperature.
- Check the cells under the microscope and record percentage confluency. Cancer cells should be banked when ~70% confluent.



- 1 Prepare  1 mL freezing media per vial as follows: complete culture media + 10% DMSO.
- 2 Detach and collect cells from a flask, by following Steps 1-6 of the protocol: **Passaging adherent cancer cell lines.**
- 3 Aspirate the supernatant, taking care to avoid disturbing the cell pellet. Resuspend the pellet in an appropriate volume of **freezing media**- depending on the number of vials being frozen. Mix well to ensure a single cell suspension.



For example, if 5 vials are being frozen from a T150, resuspend the cell pellet in

 5 mL of freezing media.

- 4 Transfer  1 mL aliquots of the cell suspension to pre-labelled 1.8ml cryovials.
- 5 Place vials in a 'CoolCell' or appropriate freezing container and store at  -80 °C overnight.

Note

Appropriate freezing containers will ensure that the liquid freezes at a controlled rate of around

 -1 °C per minute at  -80 °C .

- 6 Transfer the vials to liquid nitrogen for long-term storage.