

Oct 19, 2019

Version 2

## Frequently sampled Insulin glucose tolerance test V.2

DOI

[dx.doi.org/10.17504/protocols.io.8e8hthw](https://dx.doi.org/10.17504/protocols.io.8e8hthw)



Timothy Nichols<sup>1</sup>, David Clemmons<sup>1</sup>

<sup>1</sup>University of North Carolina at Chapel Hill

Diabetic Complications Consortium

Tech. support email: [rmcindoe@augusta.edu](mailto:rmcindoe@augusta.edu)



Lili Liang

### Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.8e8hthw>

External link: <https://www.diacomp.org/shared/document.aspx?id=50&docType=Protocol>



**Protocol Citation:** Timothy Nichols, David Clemmons 2019. Frequently sampled Insulin glucose tolerance test. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.8e8hthw>

**License:** This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working

**We use this protocol and it's working**

**Created:** October 19, 2019

**Last Modified:** October 19, 2019

**Protocol Integer ID:** 28864

**Keywords:** glucose tolerance, insulin sensitivity, cardiovascular, nephropathy, neuropathy, retinopathy, uropathy, insulin glucose tolerance test summary, sampled insulin glucose tolerance test summary, insulin sensitivity in pig, assessment of insulin sensitivity, diabetic complication, diacomp, assay

## Abstract

### Summary:

This assay is used by the DiaComp to measure glucose tolerance and insulin sensitivity in pigs.

### Diabetic Complications:



Cardiovascular



Nephropathy



Neuropathy



Retinopathy



## Uropathy

### Reference:

1. Bergman RN, Finegood DT, Ader M: Assessment of insulin sensitivity in vivo. *Endocr Rev* 1985, 6:45-86.

### Materials

#### MATERIALS

YSI Glucose Analyzer YSI Life Sciences

ICN Insulin RIA kit

#### Reagents Quantity Required

Reagent/Material	Quantity Required
Intravenous catheter	2
ICN Insulin RIA kit	1 kit
YSI Glucose Analyzer	1

### Troubleshooting

## 1 **FSIGT or Bergman analysis**

Pigs are studied after an overnight fast. The food intake of the animals is monitored for 3 days prior to the fast to ensure adequate carbohydrate intake. Two intravenous catheters are placed, one for sampling and one for infusing glucose and insulin. A bolus of glucose (0.3 gm/kg) is administered as a 50% solution over ~5 min. Blood samples are obtained at -15, -10, -5, -1, 0, 2, 3, 4, 5, 6, 8, 10, 12, 14, 16, and 19 minutes. At 20 minutes an insulin bolus (0.03U/kg) is injected and frequent blood samples for insulin and glucose measurements are collected up to the 180 minute time point. Insulin is measured by RIA (ICN) and glucose is measured on a YSI instrument (Yellow Springs, Ohio). The data were analyzed by the Bergman method to calculate an insulin sensitivity index ( $S_I$ ) using MINMOD Millennium version 6.02.<sup>1</sup>