



Oct 23, 2024

Fresh Cell Isolation RNA Extraction

DOI

<https://dx.doi.org/10.17504/protocols.io.261gerpd7l47/v1>

Malu G Tansey¹

¹University of Florida



Janna Jernigan Posey

Indiana University

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.261gerpd7l47/v1>

Protocol Citation: Malu G Tansey 2024. Fresh Cell Isolation RNA Extraction . **protocols.io**

<https://dx.doi.org/10.17504/protocols.io.261gerpd7l47/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: October 23, 2024



Last Modified: October 23, 2024

Protocol Integer ID: 110754

Keywords: ASAPCRN, fresh cell isolation rna extraction, extracting rna, rna extraction, qiagen rneasy mini kit, isolated immune cell population, rna, immune cell population, extraction

Funders Acknowledgements:

Aligning Science Across Parkinson's (ASAP) initiative

Grant ID: ASAP-020527

Abstract

This protocol describes the process of extracting RNA from freshly isolated immune cell populations using the Qiagen RNeasy Mini Kit.

Guidelines

Always wear gloves and spray with RNaseZap and spread solution all over gloves.

Make sure to use RNA only pipettes and filter tips.

Perform all steps of the procedure at room temperature. During the procedure, work quickly.

Perform all centrifugation steps at 20–25°C in a standard microcentrifuge. Ensure that the centrifuge does not cool below 20°C.

Materials

RNeasy mini kit (74106, Qiagen)

QIAshredder (79656, Qiagen)

96-100% Ethanol

Nuclease Free Water (AM9938, Thermo Scientific)

2-mercaptoethanol (M6250, Sigma Aldrich)

RNaseZAP (AM9780, Ambion)

1.5ml tubes

Safety warnings

⚠ Wear appropriate PPE and work with concentrated β -ME and ethanol in a chemical safety hood.



Before start

Prepare BME lysis buffer

Add 20 μ l β -ME per 1 ml Buffer RLT. Dispense in a fume hood. Buffer RLT containing β -ME can be stored at room temperature for up to 1 month.

Prepare RPE buffer

Add 4 volumes of ethanol (96–100%) as indicated on the bottle to concentrate buffer RPE to obtain a working Solution.



RNA Extraction

- 1 Resuspend cell pellet in 350uL of BME+RLT lysis buffer and vortex well.
- 2 Transfer lysis solution to (purple) qiashredder tubes and then spin 21300xg (full speed) for 2min.
- 2.1 Can freeze flow-through of cell lysate at this stage and store in -80 degrees celsius for several months
- 3 Add 350uL of 70% ethanol to the homogenized lysate, and mix well by pipetting.
- 4 Transfer 700 µl of the sample, including any precipitate that may have formed, to an RNeasy spin column placed in a 2 ml collection tube. Close the lid gently, and centrifuge for 30s at 9,391xg (10,000rpm). Discard the flow-through .
- 5 Add 700 µL Buffer RW1 to the RNeasy spin column. Close the lid gently, and centrifuge for 30s at 9,391xg (10,000 rpm) to wash the spin column membrane. Discard the flow-through.
- 6 Add 500 µl Buffer RPE to the RNeasy spin column. Close the lid gently, and centrifuge for 30s at 9,391xg (10,000rpm) to wash the spin column membrane. Discard the flow-through.
- 7 Repeat step 6.
- 8 Spin for 2 minutes at full speed to ensure any residual buffers are gone.
- 9 Place the RNeasy spin column in a new 1.5 ml collection tube. Add 20µl RNase-free water directly to the spin column membrane. Close the lid gently, and centrifuge for 1 min at 9,391xg (10,000rpm) to elute the RNA.



- 10 After first spin transfer elution back to the membrane of the RNeasy spin column and spin again for 1 min at 9,391xg (10,000rpm) (re-elution results in greater RNA yield). Samples are now ready for downstream assays or should be stored at -80 degrees.