

Feb 04, 2025

Functional assays via live–cell imaging

 In 1 collection

DOI

<https://dx.doi.org/10.17504/protocols.io.3byl4wkn8vo5/v1>

Maria Jose Perez J.¹

¹Institut Imagine

Team Deleidi



Federico Bertoli

ASAP

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.3byl4wkn8vo5/v1>

Protocol Citation: Maria Jose Perez J. 2025. Functional assays via live–cell imaging. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.3byl4wkn8vo5/v1>

License: This is an open access protocol distributed under the terms of the **[Creative Commons Attribution License](#)**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited



Protocol status: Working

We use this protocol and it's working

Created: November 26, 2024

Last Modified: February 04, 2025

Protocol Integer ID: 112817

Keywords: imaging functional assay, functional assay, cell, imaging, assay

Funders Acknowledgements:

ASAP

Abstract

Functional assays via live–cell imaging

Functional assays via live–cell imaging

- 1 Incubate cells with specific fluorescent dyes tailored to each desired measurement for 30 minutes at 37°C.
- 2 Wash cells twice with PBS
- 3 Conduct imaging using a Leica TCS SP8 confocal microscope equipped with a 63×/1.4 numerical aperture oil immersion objective
- 4 Perform all measurements and analyses on data obtained from over 40 cells across three independent experiments

For mitochondrial function assessment

- 5 Use 100 nM tetramethylrhodamine methyl ester perchlorate (TMRM) and 100 nM MitoTracker Green (Invitrogen)
- 6 Use an excitation wavelength of 568 nm for TMRM and 488 nm for MitoTracker Green
- 7 Measure emitted fluorescence at wavelengths shorter than 574 nm for TMRM and between 490–530 nm for MitoTracker Green
- 8 Quantify TMRM intensity in mitochondria stained with MitoTracker Green using Fiji (version 2.7.0, RRID: SCR_002285)
- 9 Set basal TMRM fluorescence intensity at 100% and determine background fluorescence after mitochondrial depolarization with FCCP. Subtract the background value from the basal value for corrected intensity

For calcium measurements

- 10 Incubate iPSC-derived neurons and astrocytes with 1 μM Fluo-4 AM (Invitrogen)



- 11 Use an excitation wavelength of 488 nm and emission range of 490-530 nm for Fluo-4 AM
- 12 Acquire images every 25 seconds
- 13 Stimulate neurons with 0.6 mM KCl and astrocytes with 1 mM glutamate for 15 minutes following a 1-minute equilibration period to establish baseline fluorescence intensity

For oxidative stress evaluation

- 14 Incubate cells with 1 μ M DCF (Invitrogen)
- 15 Use an excitation wavelength of 488 nm and emission range of 490-530 nm
- 16 Measure fluorescence intensity using Fiji (version 2.7.0, RRID: SCR_002285)
- 17 Normalize fluorescence intensity values to baseline values.