#### Oct 26, 2020

# Fungal CTAB DNA Extraction

In 1 collection

DOI

dx.doi.org/10.17504/protocols.io.bhx8j7rw

## Derreck Carter-House<sup>1</sup>, Jason E Stajich<sup>1</sup>, Sarah Unruh<sup>1</sup>, Tania Kurbessoian<sup>1</sup>

Chytrid Fungi - Molecula...

<sup>1</sup>UC Riverside

High molecular weight D...

8

Jason E Stajich





DOI: <u>dx.doi.org/10.17504/protocols.io.bhx8j7rw</u>

**Protocol Citation:** Derreck Carter-House, Jason E Stajich, Sarah Unruh, Tania Kurbessoian 2020. Fungal CTAB DNA Extraction . **protocols.io** <u>https://dx.doi.org/10.17504/protocols.io.bhx8j7rw</u>

1 more workspace

**License:** This is an open access protocol distributed under the terms of the **<u>Creative Commons Attribution License</u>**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

#### Protocol status: Working We use this protocol and it's working

Created: June 26, 2020

Last Modified: November 02, 2020

Protocol Integer ID: 38624

Keywords: Fungi, high molecular weight DNA, CTAB, sequencing



## Abstract

This protocol is a CTAB DNA extraction method for filamentous fungi. Its purpose is to extract high molecular weight genomic DNA for genome sequencing.

## Materials

### MATERIALS

🔀 RNase A (10 mg/mL) Thermo Fisher Scientific Catalog #EN0531

- 🔀 PVP Sigma
- Sodium Acetate 3M, pH 5.2 Thermo Scientific Catalog #R1181
- X potassium acetate Sigma Aldrich Catalog #P1190
- X Proteinase K Sigma Aldrich Catalog #P2308
- BUFFER A: 0.35 M sorbitol; 0.1 M Tris-HCl, pH 9; 5 mM EDTA, pH 8
- BUFFER B: 0.2 M Tris-HCI, pH 9; 50 mM EDTA, pH 8; 2 M NaCI; 2% CTAB
- BUFFER C: 5% Sarkosyl (N-lauroylsarcosine sodium salt SIGMA L5125)
- Polyvinylpyrrolidone (PVP) 1 %
- Proteinase K (20 mg/ml)
- Potassium Acetate 5M (KAc precipitate polysaccharides) pH 7.5
- RNAse A (10 mg/ml)
- (PCI) Phenol:Chloroform:Isoamyl alcohol (25:24:1)
- (CI)Chloroform:Isoamyl alcohol (24:1)
- Sodium Acetate (NaAc) 3M
- Isopropanol 100%
- Ethanol 70%

## **Protocol materials**

- **X** RNase A (10 mg/mL) **Thermo Fisher Scientific Catalog #**EN0531
- **X** PVP Merck MilliporeSigma (Sigma-Aldrich)
- Sodium Acetate 3M, pH 5.2 Thermo Scientific Catalog #R1181
- X potassium acetate Merck MilliporeSigma (Sigma-Aldrich) Catalog #P1190
- X Proteinase K Merck MilliporeSigma (Sigma-Aldrich) Catalog #P2308
- 🔀 3M Sodium Acetate solution
- 🔀 Ethanol
- 🔀 Ethanol
- 🔀 Qubit™ dsDNA BR Assay Kit Thermo Fisher Catalog #Q32853
- X Ethidium bromide [EB, EtBr] Bio Basic Inc. Catalog #EB0195.SIZE.1g

DNA	A Extraction steps	1m
1	Each tube of lysis buffer will be split in half so prepare one tube of lysis buffer for two samples. Prepare Lysis Buffer by adding to each 2mL microcentrifuge tube $4650 \mu$ L $\mu$ L Buffer A, $4650 \mu$ L $\mu$ L Buffer B, $4260 \mu$ L $\mu$ L Buffer C, $4175 \mu$ L $\mu$ L [M] 1 % (v/v) PVP, and $410 \mu$ L Proteinase K to microcentrifuge tube, mix, and then split equally into two 2 mL tubes.	
	Note	
	65 degree preheated solutions are much less viscous and allow for better mixing and equal splitting.	
2	Place in hot plate and heat to 8 65 °C.	
3	Heat a water bath to 37 °C for use in step 12 of this protocol	
4	Take your liquid nitrogen, carefully and completely immerse the inside of your mortar.	1m
	Add your pestle inside the liquid to super cool it as well.	
	Add your fungi using a spatula cleaned with ethanol, make sure to grab only the fungi and as little agar as possible.	
	Grind fungal tissue in liquid nitrogen, add 🛛 100 mg of tissue to each tube.	
	Note	
	Take your time on this step by slowly crushing the pieces of frozen fungi. Make sure the mortar stays cold by carefully adding more liquid nitrogen. (Add a little at a time to ensure the crushed fungi remains in the mortar). Grinding can take up to a minute, make sure to put enough force to get a fine powder.	

5	Incubate 😒 00:30:00 at 🖡 65 °C mixing by inversion frequently ( 😒 00:02:00 -
	Note
	Sometimes this step can generate gas due to heat expansion. As you are mixing the inversion, make sure to open the cap a bit to relieve the pressure.
6	Add $\boxed{280 \ \mu L}$ 5M Potassium Acetate to each tube, mix by inversion, incubate on ice for $\bigcirc 00:05:00$ .
7	Add $4500 \mu\text{L}$ - $4700 \mu\text{L}$ (as much as the microcentrifuge tube can reasonably hold) Phenol:Choloform:Isoamyl alcohol, 25:24:1 (PCI), mix by inversion (>5 min) or vortex briefly (5-10 seconds) then incubate for 2 min at room temp (RT).
	Note
	DO NOT overfill the tube though. The phenol:chloroform:isoamyl alcohol can seep out of the edges and make it very slippery to hold. (Vortexing will be difficult).
	If you have over-filled, tip out a bit of your solution onto a napkin and throw it away.
8	Spin at 6000 x g, 00:10:00
9	Take supernatant, transfer to new 2 ml microcentrifuge tube add equal volume Choloform:Isoamyl alcohol, 24:1 (CI) (usually about 1000ul).
	Note
	If you can, use cut pipette tips here to ensure DNA doesn't shred due to the force of pipetting through a small opening. Cut pipette tips must be autoclaved before use.

10	Mix by inversion ( 🕥 00:05:00 at least ) then incubate at room temperature for
	00:02:00
11	<b>3</b> 6000 x g, 00:10:00
12	Take supernatant (usually $\boxed{4}$ 700 $\mu$ L) and add it to a new 2 ml tube then: a. RNAse treatment ( $\boxed{4}$ 2.5 $\mu$ L RNAse A, $\boxed{6}$ 37 °C, $\boxed{6}$ 01:30:00 - $\boxed{6}$ 02:00:00
	) b. Optional 1-2 additional CI washes
	Note
	Additional CI washes will result in cleaner DNA but will also have a lower yield.
13	Add 1/10 vol 3M Sodium Acetate, mix, add 1 vol Isopropanol
	X 3M Sodium Acetate solution
	Note
	How to tell what 1/10 volume will be? Check the side of the eppendorf tube with your liquid inside. Is the line at 750uL? Divide this by 10 and that will be the volume of Sodium Acetate you add. (in this case 75uL)
14	<b>Gently</b> mix by inversion and incubate at Room temperature 00:05:00, should
	start to see lots of DNA threads.

15

# Note This may look like a super sugary solution when you add water to it. It'll have ripples inside. Don't worry if this ends up dissolving after some time. Equipment NAME Centrifuge TYPE **Benchtop Centrifuge** BRAND Eppendorf SKU 5405000441 LINK https://online-shop.eppendorf.us/US-en/Centrifugation-44533/Centrifuges-44534/Centrifuge-5425-PF-243560.html SPECIFICATIONS Any benchtop centrifuge will suffice

**3**000 x g, 00:02:00 **pour out the supernatant.** 

Note A tip here is to align your tubes all in the same orientation. This will help with clear DNA or very low yields of DNA. The next step requires you to pour our the supernatant. It's good to have a general idea of where the DNA is as you pour out. Make sure to pour slowly and in one step (vs pouring a little and orienting the tube to stand, back and forth which essentially mixes your DNA and dislodges it from the edge of the tube).

16 Wash with 4 ImL freshly prepared, cold [M] 70 % (v/v)

🔀 Ethanol

## 17

18

Equipment	
Centrifuge	NAME
Benchtop Centrifuge	TYPE
Eppendorf	BRAND
5405000441	SKU
https://online-shop.eppendorf.us/US-en/Centrifugation-44533/Centri 44534/Centrifuge-5425-PF-243560.html	fuges- <sup>LINK</sup>
Any benchtop centrifuge will suffice	SPECIFICATIONS
Note	
🔀 Ethanol emove as much EtOH as possible before drying.	
Pour ethanol out slowly and remove excess EtOH	
erv pellet at RT for 🔊 00:10:00 - 🔊 00:15:00 and/or 📱 65 °C for	< <
3 00:02:00 to dry any leftover ethanol OR if you have a vacuum char	nber, place you
amples there for 10-20 minutes under constant vacuum.	
. Resuspend in $\boxed{100 \ \mu L}$ - $\boxed{100 \ \mu L}$ TE (adjusted to $\bigcirc$ 9) at	₿ 65 °C
DNA still has high protein:	
	add egual
optional CI wash add 600-800 TE buffer at 📲 65 °C , resuspend DNA,	

	Note	
	*For example, take 500-600 uL supernatant if added 800 uL CI	
19	Assess DNA quality on NanoDrop. DNA concentration may not be accurate but the absorbance values will determine purity of the sample.	e the
	Equipment	
	NanoDrop <sup>™</sup> 3300 Fluorospectrometer	
	Fluorospectrometer TYPE	
	NanoDrop <sup>™</sup> 3300 BRAND	
	nd-3300 SKU	
	https://www.thermofisher.com/order/catalog/product/ND-3300#/ND-3300 <sup>LINK</sup>	
	Note	
	Remember: Ranges for high quality DNA 260:280 - 1.80-2.00 260:230 - 2.00-2.20	
20	Access DNA concentration with Oubit Elucromator with DNA Broad Dance (DD) or	
	Note	say.
	* Dillute $\boxed{I}_2 \mu$ of DNA with water and use all $\boxed{I}_10 \mu$ for the Qubit Assay.	

X Qubit™ dsDNA BR Assay Kit **Thermo Fisher Catalog** #Q32853

21 Assess DNA quality on Agarose gel

a. Make 0.8-1% agarose gel with  $\Delta$  3  $\mu$ L EtBr solution per 100mL of gel

X Ethidium bromide [EB, EtBr] Bio Basic Inc. Catalog #EB0195.SIZE.1g

b. Mix equal parts loading dye and DNA sample and load into wells (leaving space for DNA ladder)

- c. Add equal amount of ladder as loading dye + DNA to the first and last well
- d. Run the gel (75-120v) for 🚫 00:45:00 to 🚫 01:00:00 or until samples have run
- 3/4 of the way through the gel
- e. Visualize bands to assess size of DNA fragments

#### Note

What you are looking for is high quality DNA.

Most should stay in the well.

What you don't want to see is shearing of DNA, which will have a large blur down your entire well.

Also a thing to look for is a general smearing at the bottom of the gel, indicating protein contamination. If you have this, follow **Part A** in **Step 18**.

For reference genome sequencing or other applications where a high volume of DNA is required, tubes of the same sample can be combined. Be sure to check the concentration of the final sample.

For applications requiring very clean DNA, it may be helpful to perform additional chloroform: isoamyl washes (Step 12b). As this reduces final yield, it may be helpful to perform multiple DNA extractions of the same organism/sample.