

Jul 10, 2023

Fura-2 imaging of ionomycin response, with and without R568; a CaSR positive modulator

 In 2 collections

DOI

<https://dx.doi.org/10.17504/protocols.io.kqdg39p4zg25/v1>

Peter Kilfeather¹

¹University of Oxford



Peter Kilfeather

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.kqdg39p4zg25/v1>

Protocol Citation: Peter Kilfeather 2023. Fura-2 imaging of ionomycin response, with and without R568; a CaSR positive modulator. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.kqdg39p4zg25/v1>

License: This is an open access protocol distributed under the terms of the **[Creative Commons Attribution License](#)**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited



Protocol status: Working

We use this protocol and it's working

Created: October 18, 2022

Last Modified: July 10, 2023

Protocol Integer ID: 71468

Keywords: translatomic profiling of dopaminergic neuron, ionomycin response, dopaminergic neuron, single cell spatial, imaging

Funders Acknowledgements:

ASAP

Medical Research Council

Abstract

Fura-2 imaging protocol to accompany Kilfeather, Khoo et al., 2023: **Single cell spatial transcriptomic and translatomic profiling of dopaminergic neurons in health, ageing and disease**

Data Capture

- 1 Prepare Fura-2 AM to 5 mM concentration in calcium-free HBSS supplemented with 20 mM HEPES.
- 2 Dilute Fura-2 AM to the final concentration of 2.5 mM in Neurobasal medium supplemented with B27 and L-glutamine, with DMSO (0.1%) or the CaSR positive modulator R568 (10 mM).
- 3 Incubate iPSC-derived dopaminergic neurons (DANs) in the solution of Fura-2 AM and drug for 1 h at 37°C, 5% CO₂ and then image on a FlexStation 3 Multi-Mode Microplate Reader (Molecular Devices) at 37°C.
- 4 Excite the dye at 340 nm and 380 nm and detect at 510 nm. Image each well every 4 s for 100 s and inject with ionomycin (final concentration 5 mM) after a baseline of 28 s.

Analysis

- 5 For the analysis, compute the 340/380 ratio, and then subtract the baseline from all the timepoints, to obtain a normalised trace for each well. Find the maximum intensity (peak amplitude) of the normalised trace and calculate the area under the curve (AUC) using the left rectangular approximation method.