


Feb 15, 2022

Gallyas-silver stain

DOI

dx.doi.org/10.17504/protocols.io.b44cqysw



Christiana.bjorkli¹, María José Lagartos Donate²

¹The Norwegian University of Science and Technology; ²University of Oslo

Sandvig lab



Christiana.bjorkli

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.b44cqysw>

Protocol Citation: Christiana.bjorkli, María José Lagartos Donate 2022. Gallyas-silver stain. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.b44cqysw>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited



Protocol status: Working

We use this protocol and it's working

Created: February 15, 2022

Last Modified: February 15, 2022

Protocol Integer ID: 58212

Keywords: pathological staining in rodent, silver staining protocol, pathological staining, silver stain, gallya, staining protocol, rodent, human tissue

Abstract

Simplified Gallyas-silver staining protocol that works for pathological staining in rodent and human tissue

Troubleshooting

Solutions that can be made in advance and stored in a fume hood

1 Make 4% of **NaOH** (sodium hydroxide) - use gloves, store in fume hood

1.1 1g **NaOH** + 25ml **H₂O**

2 Make 2 × 600ml **dH₂O** – use gloves, store in fume hood

2.1 1.5g **ammonium nitrate** + 1.2g **silver nitrate** + 3.6g of 4% **NaOH**

Safety information

Ammonium nitrate can cause an explosion with metals

Reaction between **silver nitrate** and ethanol is explosive

Reaction between **silver nitrate** and **NaOH** can form inflammable gases/vapors

3 Make **acetic acid** – use gloves, store in fume hood

3.1 0.5%: 3ml **acetic acid** in 597ml **H₂O**

0.1%: 0.6ml **acetic acid** in 600ml **H₂O**

0.05%: 50ml 0.5% **acetic acid** + 450ml **H₂O**

4 Make 0.2% **potassium ferricyanide** – can be stored up to a week, use gloves, store in fume hood

4.1 0.6g **potassium ferricyanide** in 300ml **H₂O**

5 Make 0.5% **sodium thiosulfate** - use gloves, store in fume hood

5.1 1.5g **sodium thiosulfate** in 300ml **H₂O**

6 Make 25%, 50%, 70%, 95% and 100% **ethanol**



Solutions that need to be prepared on the same day of staining

7 Make **pyridine + acetic acid** - use gloves, store in fume hood

7.1 100ml **pyridine** + 50ml **acetic acid**

Safety information

Pyridine is incompatible with rubber, plastics, and metals! Very important to use butyl-rubber gloves with this chemical

8 Make **silver nitrate** (add in this order) - use gloves, store in fume hood

8.1 100ml **dH2O**
0.192g **ammonium nitrate**
0.2g **silver nitrate**
0.6ml of 4% **NaOH**
pH must be 7.5

Safety information

Ammonium nitrate can cause an explosion with metals
Reaction between **silver nitrate** and **ethanol** is explosive
Reaction between **silver nitrate** and **NaOH** can form inflammable gases/vapors

9 Make **A + B + C** (1L) - use gloves, store in fume hood

9.1 **Solution A**
5g anhydrous sodium carbonate + 100ml dH2O

Solution B (add in this order)
100ml dH2O + 0.19g ammonium nitrate + 0.2g silver nitrate + 1g silicongulistic
(tungolistic) acid

Solution C (add in this order)



100ml dH₂O + 0.19g ammonium nitrate + 0.2g silver nitrate + 1g silicongulistic (tungolistic) acid + 0.66ml 37% formaldehyde

A + B + C - use gloves, store in fume hood

50ml **A** + 37.5ml **B** + 37.5ml **C** (25±2°C)

Safety information

Ammonium nitrate can cause an explosion with metals

Reaction between **silver nitrate** and **ethanol** is explosive

Reaction between **silver nitrate** and **NaOH** can form inflammable gases/vapors

Formaldehyde may cause cancer

Consumables for staining

- 10
- 17 x glass petri dishes
 - 20 x glass pipettes
 - 1 x rubber pump for glass pipettes
 - Butyl rubber gloves
 - Lab coat with long arms
 - Protective eyewear
 - Face mask

Note

Need to use glass petri dishes and glass pipettes for each wash

Gallyas-silver staining

11

	Step	Reagent	Minutes	Comments
	1	dH₂O	3 or more	Put in dH ₂ O while preparing other solutions
	2	Pyridine + acetic acid	60	On shaker
	3	50% ethanol	3	
	4	25% ethanol	3	



5	0.05% acetic acid	2.5		
6	0.1% acetic acid	2.5		
7	0.05% acetic acid	10	Leave longer if necessary, to prepare silver nitrate	
8	Silver nitrate	60	On shaker – make A+ B + C	
9	0.5% acetic acid	10		
10	A + B + C	6 (15min human)	Do not let the tissue become dark, time is dependent on sample	
11	0.5% acetic acid	1	New	
12	0.2% potassium ferricyanide	5	On shaker	
13	dH2O	1		
14	0.5% acetic acid	1	New, can be reused	
15	0.5% sodium thiosulfate	2	On shaker	
16	dH2O	4		
17	dH2O	4		
18	0.5% acetic acid	1	Use the one from the previous step	
19	A + B + C	3-4	On shaker	
20	0.5% acetic acid		Use the one from the previous step	
21	0.2% potassium ferricyanide	3	On shaker	
22	dH2O	1		
23	0.5% sodium thiosulfate	2	On shaker	
24	dH2O	4		



25	dH2O	4		
26	0.5% acetic acid	1	New, can be reused	
27	A + B + C	3-4	On shaker	
28	0.5% acetic acid	1	Use the one from the previous step	
29	0.2% potassium ferricyanide	10	On shaker	
30	dH2O	1		
31	0.5% sodium thiosulfate	2	On shaker	
32	dH2O	Fast wash		
33	dH2O	5		
34	Dehydrate	3 in each	50% - 70% -95% - 100% - 100% ethanol	
35	Xylene x 2	3 in each		
36	Coverslip			