

Apr 05, 2024

Gallyas Silver Staining

DOI

<https://dx.doi.org/10.17504/protocols.io.5qpvo366zv4o/v1>

madalynn.erb Erb¹

¹Van Andel Research Institute

Madalynn's Workspace



madalynn.erb Erb

Van Andel Research Institute

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.5qpvo366zv4o/v1>

Protocol Citation: madalynn.erb Erb 2024. Gallyas Silver Staining. protocols.io
<https://dx.doi.org/10.17504/protocols.io.5qpvo366zv4o/v1>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited



Protocol status: Working

We use this protocol and it's working

Created: January 29, 2024

Last Modified: May 31, 2024

Protocol Integer ID: 94364

Keywords: ASAPCRN, mouse brain sections stain, staining gallyas silver, gallyas silver, mouse brain section, mouse brain, degenerating neurite, mouse

Disclaimer

The [protocols.io](https://www.protocols.io) team notes that research involving animals and humans must be conducted according to internationally-accepted standards and should always have prior approval from an Institutional Ethics Committee or Board.

Abstract

Gallyas Silver staining in mouse brain sections stains degenerating neurites black.

Day 1

- 1 For silver staining use the FD NeuroSilver™ Kit II (FD NeuroTechnologies INC Cat# PK301)
- 2 Staining of free-floating brain sections is performed in glass staining dishes (Pyrex 36754-60) using 8-section staining nets (Ted Pella 36154-64)
 - 2.1 The sections can be transferred between wells using a paint brush.
- 3 Wash sections 35µm mouse brain sections in 0.1M Phosphate buffer (PB) 3 times (5 minutes each wash) at room temp to remove cryoprotectant solution
- 4 Incubate sections in 0.1M PB + 4% paraformaldehyde pH 7.4 for 7 days at 4°C
 - 4.1 Cover dish in parafilm
 - 4.2 Lightly shake sections during this step - *we use a shaker in the cold room*

Day 8

- 5 *Perform all silver staining steps in a chemical fume hood.*
 - 5.1 *Take shaker into fume hood for this purpose.*
 - 5.2 *Collect all silver stain waste in large beaker for hazardous waste disposal*
- 6 Wash sections in miliQ H₂O 2 times for 5 minutes each wash at room temp
 - 6.1 While sections are washing make up the A/B solution



- 6.2 Make 80ml per dish : 40ml sol A + 40ml sol B
- 7 Wash sections in A/B solution for 2 times for 10 min each wash at room temp
- 7.1 Use 25mL solution for each wash
- 8 Wash sections in A/B solution + E for 10 min at room temp
- 8.1 8mL sol A/B + 1 drop solution E
- 8.2 For 1 dish: 24ml sol A/B + 3 drops sol E
- 8.3 Wash sections in C/F solution for 2 min at room temp
- 8.4 50mL sol C + 2 drops sol F
- 8.5 Use 25mL per wash
- 9 Wash sections in C/F solution for 4 min at room temp
- 9.1 *Incubation longer than 4 min at this step decreases background and decreases signal. Shorter incubation increases signal but also increases background.*
- 10 Wash sections in D/F solution 1 time for 5 min at room temp
- 10.1 25mL sol D + 1 drop sol F



- 10.2 Make 25mL per dish
- 11 Wash sections in miliQ H₂O 2 times for 3 min each wash at room temp
- 12 Wash sections in 1X solution G 2 times for 5 min each wash
- 12.1 Dilute 10X sol G 1:10 in miliQ H₂O
- 12.2 Make 80mL per dish
- 12.3 Store sections in 1X sol G
- 13 Use a petri dish filled with 1X sol G to mount the sections on superfrost plus slides (Fisher Scientific 12-550-15)
- 14 This is easiest to do using a medium sized paint brush
- 15 Let sections dry at room temp after mounting
- 16 Protect slides from light
- 17 Clear sections in fresh xylene 3 times for 3 min each wash at room temp
- 17.1 **DO NOT** dehydrate in EtOH – it will cause loss of staining.
- 17.2 Use fresh xylene. Old xylene can contains residual EtOH from previous experiments. Residual EtOH will remove silver staining.
- 18 Use Entellan mounting medium (Electron Microscopy Sciences 14802) to coverslip slides



- 18.1 Apply one line of Entellan across the middle of the slide
- 18.2 Cover with glass coverslip and gently press down with forceps
- 18.3 Let slides dry in the fume hood for one hour – then dry on the bench overnight at room temp (protect slides from light) Gently remove excess Entellan with a razor blade