

Mar 31, 2023

Gcase co-immunoprecipitation

DOI

<https://dx.doi.org/10.17504/protocols.io.14egn2ompg5d/v1>

michela.deleidi¹, Federico Bertoli¹, Maria Jose Perez J.¹, Hariam Raji¹

¹German Center for Neurodegenerative Diseases (DZNE), Tübingen, 72076 Germany



Federico Bertoli

ASAP

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.14egn2ompg5d/v1>

Protocol Citation: michela.deleidi , Federico Bertoli, Maria Jose Perez J., Hariam Raji 2023. Gcase co-immunoprecipitation. protocols.io <https://dx.doi.org/10.17504/protocols.io.14egn2ompg5d/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: March 31, 2023

Last Modified: March 31, 2023

Protocol Integer ID: 79775

Keywords: protein interactions between the enzyme glucocerebrosidase, enzyme glucocerebrosidase, glucocerebrosidase, other proteins in human ipsc, neural precursor cell, protein, derived neural precursor cell, gcase, other protein, protein interaction, immunoprecipitation, enzyme, human ipsc

Abstract

We developed this protocol to identify protein-protein interactions between the enzyme glucocerebrosidase (GCase) and other proteins in human iPSC-derived Neural Precursor Cells.

Attachments



[676-1426.docx](#)

13KB

Materials

 Pierce & Warriner; IP Lysis Buffer **Thermo Fisher Catalog #87787**

 Pierce Protease and Phosphatase Inhibitor Mini Tablets **Thermo Fisher Catalog #A32959**

 Pierce & Warriner; Crosslink Magnetic IP/Co-IP Kit **Thermo Fisher Catalog #88805**

 Pierce Protein A/G Magnetic Beads **Thermo Fisher Scientific Catalog #88802**



Gcase co-immunoprecipitation

2h 5m

- 1 Wash cells 1X with phosphate-buffered saline (PBS, Sigma–Aldrich) and detach using Accutase.
- 2 Pellet the cell suspension at 280 rcf, 23°C, 00:05:00 . 5m
- 3 Lyse the pellets in IP/lysis buffer (Thermo Fisher, #87787) supplemented with a protease/phosphatase inhibitor cocktail (Pierce, #A32959).
- 4 Carry out coimmunoprecipitation using the Thermo Fisher Pierce Crosslink Magnetic IP/Co-IP Kit (#88805) according to the manufacturer's instructions:
- 5 Prewash 25 μL of Pierce protein A/G magnetic beads (Thermo Fisher, #88802-3) twice with 1X Modified Coupling Buffer and incubate with 10 μg of GBA MaxPab polyclonal rabbit antibody (Abnova) or normal rabbit IgG (Covalab, #pab01004-P) on a rotating wheel Overnight at 4 °C .
- 6 The following day, crosslink the antibody to the beads with a 0.25 millimolar (mM) DSS solution for 01:00:00 on a rotating wheel at Room temperature . 1h
- 7 Incubate crosslinked magnetic beads Overnight with a total of 7 mg of protein for each lysate. 1h
- 8 Elute Coimmunoprecipitated proteins from the beads with 60 μL of the kit-provided Elution buffer 2.0 (Pierce, #88805) and neutralize with 6 μL of Neutralization Buffer provided with the kit.
- 9 Prepare samples for Western blotting by adding 5x Lane buffer +10% DTT 1 Mass Percent to a final concentration of 1X.

Note

- Each western blot input sample loaded corresponds to a total of 50 μg of protein.
- Each western blot CoIP sample loaded corresponded to the total of each elution product (66 μL).

