



May 08, 2025

Gene expression analysis

 In 2 collections

DOI

<https://dx.doi.org/10.17504/protocols.io.yxmvmvmywebv3p/v1>

Lucas Blasco-Agell^{1,2}, Meritxell Pons-Espinal^{1,2}, Jara Montero-Muñoz^{1,2}, Angel Raya^{1,2}, Antonella Consiglio^{1,2}, Miquel Vila³

¹Department of Pathology and Experimental Therapeutics, Bellvitge University Hospital- IDIBELL, 08908 Hospitalet de Llobregat, Spain;

²Institute of Biomedicine of the University of Barcelona (IBUB), Carrer Baldori Reixac 15-21, Barcelona 08028, Spain;

³VHIR-CIBERNED-ASAP

Vilalab Public



Miquel Vila

VHIR-CIBERNED-ASAP

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.yxmvmvmywebv3p/v1>

Protocol Citation: Lucas Blasco-Agell, Meritxell Pons-Espinal, Jara Montero-Muñoz, Angel Raya, Antonella Consiglio, Miquel Vila 2025. Gene expression analysis. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.yxmvmvmywebv3p/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: April 17, 2025

Last Modified: May 08, 2025

Protocol Integer ID: 130936

Keywords: human microglia cell, neuromelanin, gene expression, neurotoxic factor, gene expression analysis

Funders Acknowledgements:

ASAP

Grant ID: ASAP-020505

Abstract

To evaluate whether human microglia cells (hMG) can be activated by neuromelanin (NM), hMG were cultured with and without the presence of NM for 24 h and then harvested to measure gene expression of selected pro-inflammatory and neurotoxic factors.

- 1 RNA is extracted using a RNeasy® Mini kit (Qiagen) according to manufacturer's instructions and quantified with a NanoDrop one spectrophotometer (Thermo Fisher Scientific).
- 2 Reverse transcription into complementary cDNA is performed with the Superscript III First-Strand Synthesis System (Invitrogen™) to a final concentration of 2 ng/μL.
- 3 Gene expression is measured by real time quantitative Reverse Transcription Polymerase Chain Reaction (qRT-PCR).
- 4 8 ng of cDNA are loaded in every well of a MicroAmp™ Optical 384-well Reaction Plate (Applied Biosystems™) with 10 nM of forward and reverse primers (Life Technologies) and PowerUp™ SYBR™ Green Master Mix (Applied Biosystems™).
- 5 Plate is run in a 7900HT Fast Real-Time PCR System with 384-well Block Module (Applied Biosystems™) following a standard cycling mode.
- 6 Cycling values are normalized to the housekeeping gene β-Actin and the relative fold gene expression of samples was calculated using the $2^{-\Delta\Delta C_t}$ method.