Dec 18, 2024

Version 1

General extraction of alkaloids in acidic medium V.1

DOI

<https://dx.doi.org/10.17504/protocols.io.5jyl8dek8g2w/v1>

Rômulo Pereira de Jesus¹, Ricardo Borges¹

¹UFRJ - Federal University of Rio de Janeiro, Brazil

LAABio-IPPN-UFRJ

Rômulo Pereira de Jesus

Federal University of Rio de Janeiro, Brazil.

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS

DOI: <https://dx.doi.org/10.17504/protocols.io.5jyl8dek8g2w/v1>

Protocol Citation: Rômulo Pereira de Jesus, Ricardo Borges 2024. General extraction of alkaloids in acidic medium.
protocols.io <https://dx.doi.org/10.17504/protocols.io.5jyl8dek8g2w/v1>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: In development

We are still developing and optimizing this protocol

Created: December 18, 2024

Last Modified: December 18, 2024

Protocol Integer ID: 116128

Keywords: general extraction of alkaloid, qualitative detection of alkaloid, separation of alkaloid, presence of alkaloid, alkaloid, phytochemical study, extraction in an acidic medium, suitable for phytochemical study, organic solvents in an acidic medium, plant sample, extraction, efficient method for the extraction, alkalization, acidic medium this protocol, ethyl acetate after alkalization, formation of orange, natural product, acidic medium, using organic solvent, solvent, ethyl acetate

Disclaimer

This protocol is for informational purposes only and does not substitute specific guidelines or applicable regulations. The use of the information contained herein is the responsibility of the user

Abstract

This protocol describes an efficient method for the extraction and qualitative detection of alkaloids in plant samples using organic solvents in an acidic medium. The procedure includes extraction in an acidic medium and partitioning with ethyl acetate after alkalization. The mixture is subjected to ultrasound and centrifugation processes, enabling the separation of alkaloids in the organic phase. The solvent is subsequently evaporated to obtain the compounds in their pure form. Qualitative identification is carried out using Dragendorff's reagent, which indicates the presence of alkaloids through the formation of orange to red precipitates. This method combines simplicity, reproducibility, and efficiency, making it suitable for phytochemical studies of natural products.

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