

Aug 24, 2023

Version 2

## General Fungal DNA Extraction V.2



Version 1 is forked from [Kasson Lab DNA Extraction](#)

DOI

[dx.doi.org/10.17504/protocols.io.3byl4q7rzvo5/v2](https://dx.doi.org/10.17504/protocols.io.3byl4q7rzvo5/v2)

Angie Macias<sup>1</sup>, Matthew T Kasson<sup>1</sup>, Brian Lovett<sup>1</sup>

<sup>1</sup>West Virginia University



Brian Lovett

USDA-ARS

### Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.3byl4q7rzvo5/v2>

**Protocol Citation:** Angie Macias, Matthew T Kasson, Brian Lovett 2023. General Fungal DNA Extraction. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.3byl4q7rzvo5/v2> Version created by **Brian Lovett**

**License:** This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working

**We use this protocol and it's working**

**Created:** August 24, 2023

**Last Modified:** August 24, 2023

**Protocol Integer ID:** 86920

**Keywords:** general fungal dna extraction, dna from various fungi, dna extraction, extracting dna, various fungi, fungi, molecular work such as pcr amplification, extraction method, pcr amplification, pcr, dna

## Abstract

This is a routine protocol for extracting DNA from various fungi. This extraction method is suitable for follow-up molecular work such as PCR amplification.

## Materials

Sterile micropestles, isopropyl alcohol, ethyl alcohol, cell lysis buffer, protein precipitation buffer, elution buffer, metal scraper.

## Protocol materials

✕ isopropyl alcohol **Merck MilliporeSigma (Sigma-Aldrich) Catalog #W292907**

✕ Cell Lysis Solution, 1000ml (for Wizard Genomic) **Promega Catalog #A7933**

✕ Nuclei Lysis Solution, 1000ml **Promega Catalog #A7943**

✕ isopropyl alcohol **Merck MilliporeSigma (Sigma-Aldrich) Catalog #W292907**

✕ Elution buffer pH 8.0 (250 mL) **Alfa Aesar Catalog #J61558**

✕ Ethyl Alcohol **Merck MilliporeSigma (Sigma-Aldrich) Catalog #E7023**

✕ Cell Lysis Solution, 1000ml (for Wizard Genomic) **Promega Catalog #A7933**

✕ Protein Precipitation Solution 350ml **Promega Catalog #A7953**

✕ isopropyl alcohol **Merck MilliporeSigma (Sigma-Aldrich) Catalog #W292907**

✕ Ethyl Alcohol **Merck MilliporeSigma (Sigma-Aldrich) Catalog #E7023**







✕ Elution buffer pH 8.0 (250 mL) **Alfa Aesar Catalog #J61558**

✕ Elution buffer pH 8.0 (250 mL) **Alfa Aesar Catalog #J61558**

✕ Cell Lysis Solution, 1000ml (for Wizard Genomic) **Promega Catalog #A7933**

## Troubleshooting

## Before you begin

- 1 Turn on hot water bath, set to  $65^{\circ}\text{C}$ .
  - 2 Pull two Eppendorf  $1.5\text{ mL}$  centrifuge tubes per sample.
    - 2.1 Label both sets of tubes with (short) sample names.
    - 2.2 Label one tube set for each sample with an "I" for  isopropyl alcohol **Merck MilliporeSigma (Sigma-Aldrich) Catalog #W292907**.
- 
- Sketch of "I"-labeled tubes (drawing from Angie Macias).
- 3 Add  $200\text{ }\mu\text{L}$  of  Cell Lysis Solution, 1000ml (for Wizard Genomic) **Promega Catalog #A7933** (or  Nuclei Lysis Solution, 1000ml **Promega Catalog #A7943**) to **tubes without "I"**.
  - 4 Add  $600\text{ }\mu\text{L}$  of  isopropyl alcohol **Merck MilliporeSigma (Sigma-Aldrich) Catalog #W292907** to **tubes labeled with "I"**.
  - 5 Place tube with  Elution buffer pH 8.0 (250 mL) **Alfa Aesar Catalog #J61558** into  $65^{\circ}\text{C}$  water bath.





## Extraction Protocol

1h 10m 3s

- 6 Sterilize some metal scrapers with flame and [M] 95 % (v/v)  
⊗ Ethyl Alcohol Merck MilliporeSigma (Sigma-Aldrich) Catalog #E7023 .
- 7 Add 1/2 pea-sized amount of fungal tissue (young hyphae) to each tube containing  
⊗ Cell Lysis Solution, 1000ml (for Wizard Genomic) Promega Catalog #A7933 .
- 7.1 Flame-sterilize and cool scrapers between samples.
- 7.2 Alternatively, pellet a pea-sized amount of mycelium grown in liquid culture and transfer to each tube.
- 8 Macerate each sample with a new, sterile micropestle until tissue is homogenous.
- 9 Add 🧪 400  $\mu$ L of  
⊗ Cell Lysis Solution, 1000ml (for Wizard Genomic) Promega Catalog #A7933 (to  
🧪 600  $\mu$ L total volume added).
- 10 Add tubes to a floating rack to allow samples to incubate directly in 🌡 65 °C water bath for ⌚ 00:30:00 . 30m
- 11 Remove samples and vortex for ⌚ 00:00:03 before returning to 🌡 65 °C water bath for ⌚ 00:30:00 . 30m 3s
- 11.1 Place a sufficient aliquot of  
⊗ Elution buffer pH 8.0 (250 mL) Alfa Aesar Catalog #J61558 in water bath to warm for Step 21.
- 12 Remove samples and allow them to cool on the bench for ⌚ 00:05:00 . 5m




13 Add  200  $\mu\text{L}$  of  Protein Precipitation Solution 350ml **Promega Catalog #A7953** to each tube and vortex for 10 seconds.

14 Centrifuge samples for  00:03:00 at  14.000 rpm .

3m



**Note**

Proteins will form a large pellet: unload samples carefully into rack.

15 Using a P1000 micropipette, transfer supernatant to each tube containing  isopropyl alcohol **Merck MilliporeSigma (Sigma-Aldrich) Catalog #W292907** and gently mix by inversion.

**Note**

It's better to leave some liquid than to carry bits of the protein pellet into the next step.

16 Centrifuge for  00:01:00 at  14.000 rpm .



1m

17 Carefully pour off the supernatant into waste container.

**Note**

Be careful to not lose your white DNA pellet!

18 Add  600  $\mu\text{L}$  of [M] 70 % (v/v)  Ethyl Alcohol **Merck MilliporeSigma (Sigma-Aldrich) Catalog #E7023** to each tube and mix gently by inversion.

19 Centrifuge for  00:01:00 at  14.000 rpm .

1m



20 Repeat Step 16.


21 Open and invert tubes onto a clean paper towel.

Note

A tube rack can be placed on the tube lids to secure inverted tubes onto the paper towel.

22 Add  100  $\mu\text{L}$  of warmed

 Elution buffer pH 8.0 (250 mL) **Alfa Aesar Catalog #J61558** to each tube.

23 Store fully-labeled tubes in a box (not a tube rack) in the   $-20\text{ }^{\circ}\text{C}$  freezer.