

Aug 13, 2023

## Generating regionally specified astrocytes

DOI

[dx.doi.org/10.17504/protocols.io.261ged3gov47/v1](https://dx.doi.org/10.17504/protocols.io.261ged3gov47/v1)

Tyra Fraser<sup>1</sup>, Lachlan Thompson<sup>1,2</sup>

<sup>1</sup>Florey Institute of Neuroscience and Mental Health; <sup>2</sup>University of Sydney



courtney.wright Wright

University of Sydney

### Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.261ged3gov47/v1>

**Protocol Citation:** Tyra Fraser, Lachlan Thompson 2023. Generating regionally specified astrocytes. **protocols.io**  
<https://dx.doi.org/10.17504/protocols.io.261ged3gov47/v1>

**License:** This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working

We use this protocol and it's working

**Created:** August 11, 2023



**Last Modified:** May 31, 2024

**Protocol Integer ID:** 86378

**Keywords:** ASAPCRN, astrocytes, stem cells, regional astrocytes, diseased astrocytes, astrocytes this protocol, specified astrocyte, cortical astrocyte, differentiation of stem cell, astrocyte, stem cell, specific ventral midbrain, midbrain, cell

**Funders Acknowledgements:**

**Michael J Fox Foundation**

**Grant ID:** ASAP-000497

## Abstract

This protocol outlines the differentiation of stem cells into regionally specific ventral midbrain and cortical astrocytes.

The protocol is a 4–5-month procedure and requires media changes 3 times a week.

## Guidelines

### Seeding density reference and expansion

- 3 spheres are plated into a 48 well plate > 24 well plate > 12 well plate > 6-well plate > T-25 > T75 > T175
- P0 is when you first seed them into 2D format
- At P10 VMDA astrocytes > MAX Passage
- At P12 Cortex astrocytes > MAX Passage
- At P10 Spinal cord astrocytes > MAX Passage



## Materials

### Material input


- PPMI stem cells
- D13 VMDA neurons
- D13 cortical neurons

### Key materials

-  Advanced DMEM/F-12 **Thermo Fisher Catalog #12634010**
-  Glutamax (100x) **Gibco - Thermo Fisher Scientific Catalog #35050-061**
-  Gibco™ N-2 Supplement (100X) **Thermo Fisher Scientific Catalog #17502048**
-  Antibiotic-Antimycotic 9100×0 [Anti-Anti] **Thermo Fisher Scientific Catalog #15240062**
-  B-27® Plus Supplement (50X) **Thermo Fisher Catalog #A3582801**
-  Basic FGF (FGF2), Human **Gold Biotechnology Catalog #1140-02**
-  EGF, Epidermal Growth Factor, human **Bio Basic Inc. Catalog #RC216-15.SIZE.100ug**
-  PDS Kit Papain Vial **Worthington Biochemical Corporation Catalog #LK003176**
-  DNase I, RNase free **Thermo Fisher Scientific Catalog #EN0525**

## Troubleshooting






## Safety warnings

 For hazard information and safety warnings, please refer to the SDS (Safety Data Sheet)










## Preparation

1 Prepare astrocyte base media by combining;

-  484.5 mL Advanced DMEM/F12
-  5 mL GlutaMAX
-  5 mL N2
-  5 mL ANTI-ANTI
-  500  $\mu$ L B-27+VITA

## Sphere formation

- 2
1. Disassociate D13 VMDA diff or D19 Cortex diff using Accutase for  00:03:00 to  00:06:00 1d 0h 15m
  2. Collect cells as small clumps and inhibit using respective diff media + Ri 1:1000
  3. Spin down cells at  1300 rpm for  00:03:00
  4. Aspirate media and resuspend in respective diff media (with factors required for diff) + Ri 1:1000.
  5. Seed 20 000 cells/ well in a 96-well round bottom plate
  6. Spin plate at  1300 rpm for  00:03:00
  7. After  24:00:00 , change the media with Astro base + EGF (20ug/ml stock; 20 ng/ml final 1:1000 dilution) and hLIF (10ug/ml stock; 20ng/ml final 1:500 dilution)








## Maintenance

- 3
- When spheres have been in hLIF and EGF for 4 weeks change the media to EGF (20ug/ml stock; 20 ng/ml final 1:1000 dilution) + FGF2 (100ug/ml stock; 20ng/ml final 1:5000 dilution)
  - EGF and FGF2 is maintained for a minimum of 2 months.

## Neurosphere slicing







- 4
1. When spheres become 0.5cm or larger it is time to chop them to axonotomise the neurons 1d 0h 5m



2. Aspirate media and slice spheres using an autoclaved blade or two needles. Slice spheres in one direction then rotate plate 90° and slice again. Continue to rotate until you've reached 360°
3. Add  5 mL of PBS and move sliced spheres into a 15ml Falcon tube
4. Let sphere settle and aspirate PBS
5. Resuspend in  0.5 mL DNase (1mL in large amount of starting material) and incubate at  37 °C for  00:05:00 . Shake tube every few minutes to distribute DNase through fragments
6. To DNase/cell mixture add  9 mL Astro base supplemented w/FGF2 and EGF and re-plate in a new low attachment 10cm plates
7. After  24:00:00 , media change spheres with  10 mL Astro base supplemented w/FGF2 and EGF to wash out DNase

## Plating down using Papain




53m

- 5
  1. When spheres have been in EGF+FGF2 for 2 months and have been sliced a minimum of 3 times they can be plated down in 2D format
  2. Coat plates using 1:80 MG
  3. Use the Worthington kit manual to reconstitute the Papain powder and make up all other reagents
  4. Transfer spheres to a 15ml falcon tube and wash with PBS-/- . Aspirate PBS -/- once spheres have settled
  5. Add Papain DNase mixture to the spheres (Usually double the amount of Papain to volume of spheres) and incubate at  37 °C for  00:15:00 to  00:30:00 . Perform triturates at every 5-minute interval  
\*\*NOTE: If your spheres are healthy and quite large you might need to manually pull them apart before adding Papain like what we do for cutting
  6. At the end of 15 minutes, triturate again and let large clumps settle to the bottom and take off all supernatant (this contains a lot of single cells)
  7. Add supernatant to new tube and inhibit with base media
  8. Add new Papain, and repeat step 3-5 until clumps are disassociated (you will not be able to dissociate the whole thing so some clumps will be left, no longer than 45 minutes).
  9. Spin cells that were inhibited at  1000 rpm for  00:04:00
  10. Aspirate supernatant and resuspend in  300 µL ovomucoid inhibitor +DNase (see Worthington manual)

53m





11. In a separate tube, add  1 mL of inhibitor albumin. Slowly transfer the cell solution onto the top of the protein gradient and leave to settle. \*\*Note: the Worthington manual says to spin down, this is not necessary
12. Transfer the top layer of cells to a fresh tube
13. Centrifuge at  1000 rpm for  00:04:00
14. Resuspend in EGF + FGF2 **(DO NOT ADD Ri)**

## Freezing astrocytes

- 6 To freeze astrocytes at the precursor stage, perform prep as normal however resuspend in Astrocyte base media + 10% DMSO

M

2w

- 7 To mature the astrocytes, maintain in CNTF 1:1000 + base media for 2 weeks