

Aug 10, 2018

Version 1

Generating Stable Transfection in Bodo saltans V.1

DOI

<https://dx.doi.org/10.17504/protocols.io.sh3eb8n>

Fatma Gomaa¹, Zuhong Li², Roberto Docampo², Peter Girguis¹, Virginia Edgcomb³

¹Harvard University; ²University of Georgia; ³Woods Hole Oceanographic Institution

Protist Research to Opti...



Fatma Gomaa

Harvard University

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.sh3eb8n>

Protocol Citation: Fatma Gomaa, Zuhong Li, Roberto Docampo, Peter Girguis, Virginia Edgcomb 2018. Generating Stable Transfection in Bodo saltans. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.sh3eb8n>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited



Protocol status: Working

We use this protocol and it's working

Created: August 10, 2018

Last Modified: August 10, 2018

Protocol Integer ID: 14619

Keywords: electroporated cell, stable transfection in bodo saltan, hours after electroporation, transfected cell, electroporation, gel electrophoresis, g418 resistant cell, agarose gel electrophoresis image, generating stable transfection, gel electrophoresis image, resistant cell, genotyping analysis, wave electroporator, wild type cell, wild cell, transfection, saltan, site of plasmid integration, dna, cells s1, cell, cells at different level, plasmid integration

Abstract

- *B. saltans* cells were electroporated using a square-wave electroporator (Nepa21, Bulldog Bio, Inc.) using one poring pulse of 200 volts with a pulse duration of 25 ms and five transfer pulses of 60 volts with a pulse duration of 99 ms, with plasmid targeting the 18S region (18S-GFP). A schematic representation of the plasmid, the target locus and the expected site of integration into the *B. saltans* genome is shown in Figure 1.
- Electroporated cells were selected with 1 µg/ml of G418, added 24 hours after electroporation. Cells were washed and subcultured into fresh selection medium every 3-4 days. G418 resistant cells started to emerge 7-9 days post-electroporation.
- Cells were processed for genotyping analysis to confirm plasmid integration 3 weeks post-electroporation. DNA was extracted from pools of transfected and wild cells using the Qiagen DNeasy Blood & Tissue kit.
- PCR analyses were used to characterize the 18S-GFP tagging using 6 sets of PCR primers, as shown in Figure 1 C.
- Gel electrophoresis image (Figure 2) showing the amplified PCR products at the expected sizes.
- Amplified PCR#1 with primer sets Ribo_tag_forward & GFP reverse (800 bp)
- Amplified PCR #2 with primer sets Neo_forward & Ribo_tag_reverse (1000 bp)
- Amplified PCR #3 with primer sets TubR & IG forward (3 bands)
- Amplified PCR #4 with primer sets Tub_forward & IG reverse (3 bands)
- Amplified PCR #5 with primer sets Ribo_tag_forward & Ribo_tag_reverse (Wild (C) cells band at 350 bp, transfected cells S1 and S2 two bands, 350 bp and 2800 bp)

Figure 1: Schematic representation of the (A); the 18S- GFP plasmid; (B) Ribosomal operon in *B. saltans* genome (C); and the expected site of plasmid integration in *B. saltans* genome through homologous regions 1 and 2 (HR1, HR2).

Figure 2: Agarose gel electrophoresis image of the amplified PCR products for *B. saltans* cells transfected with 18S-GFP Cassette (S1 and S2) and the wild type cells (C). The primers sets and the expected product sizes are mentioned above in the text.

Figure 3: *B. saltans* cells transfected with 18S-GFP tagging cassette. On the left: a light microscopy image showing the *B. saltans* cells (around 4 cells at different level). On the right: GFP fluorescent signal detected in the transfected cell. Scale bar 10 μ m.

Attachments



Protocol_stable tran...

1.9MB

