



May 31, 2024

# 🌐 Generation of Flp-In™ T-REx™ 293 cells stably expressing VPS13C<sup>mClover</sup> under a tetracycline inducible promoter

📖 [Nature Cell Biology](#)

DOI

[dx.doi.org/10.17504/protocols.io.36wgqnnnxgk5/v1](https://dx.doi.org/10.17504/protocols.io.36wgqnnnxgk5/v1)

Xinbo Wang<sup>1,2</sup>, Shujun Cai<sup>1,2</sup>, Will Hancock-Cerutti<sup>1,2,3</sup>, Pietro De Camilli<sup>1,2</sup>

<sup>1</sup>Departments of Neuroscience and of Cell Biology, Howard Hughes Medical Institute, Program in Cellular Neuroscience, Neurodegeneration and Repair, Yale University School of Medicine, New Haven, Connecticut 06510, USA;

<sup>2</sup>Aligning Science Across Parkinson's (ASAP) Collaborative Research Network, Chevy Chase, MD, 20815;

<sup>3</sup>Interdisciplinary Neuroscience Program and MD-PhD Program, Yale University School of Medicine, New Haven, Connecticut 06510, USA



**Xinbo Wang**

Yale University

## Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.36wgqnnnxgk5/v1>



**External link:** <https://doi.org/10.1038/s41556-025-01653-6>

**Protocol Citation:** Xinbo Wang, Shujun Cai, Will Hancock-Cerutti, Pietro De Camilli 2024. Generation of Flp-In™ T-REx™ 293 cells stably expressing VPS13C<sup>mClover</sup> under a tetracycline inducible promoter. **protocols.io**

<https://dx.doi.org/10.17504/protocols.io.36wgqnnnxgk5/v1>

**Manuscript citation:**

Wang, X., Xu, P., Bentley-DeSousa, A. *et al.* The bridge-like lipid transport protein VPS13C/PARK23 mediates ER-lysosome contacts following lysosome damage. *Nat Cell Biol* **27**, 776–789 (2025). <https://doi.org/10.1038/s41556-025-01653-6>

**License:** This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working

**We use this protocol and it's working**

**Created:** May 22, 2024

**Last Modified:** May 31, 2024

**Protocol Integer ID:** 100990

**Keywords:** ASAPCRN, tetracycline inducible promoter, generation of flp, rex, vps13c

## Abstract

This method describes the generation of Flp-In™ T-REx™ 293 cells stably expressing internally mClover tagged VPS13C under a tetracycline inducible promoter.

## Materials

### Specialized Reagents:

⊗ Flp-In™ T-REx™ 293 Cell Line Thermo Fisher Scientific Catalog #R78007

⊗ VPS13C<sup>m</sup>Clover3 addgene Catalog #118760

⊗ pcDNA<sup>+</sup>5/FRT/TO Vector Kit Thermo Fisher Catalog #V652020

⊗ pOG44 Flp-Recombinase Expression Vector Thermo Fisher Catalog #V600520

⊗ Blasticidin S HCl (10 mg/mL) Thermo Fisher Catalog #A1113903

⊗ Zeocin<sup>+</sup> Selection Reagent Thermo Fisher Catalog #R25001

Hygromycin B (ThermoFisher Catalog # J67371.8EQ)

Tetracycline (ThermoFisher Catalog # A39246)

### Reagents to prepare

### Culture media:

	A	B
	Dulbecco's modified Eagle's medium (DMEM)	
	fetal bovine serum (FBS)	10%
	zeocin	200 µg/ml
	blasticidin S.	15 µg/ml

### Post transfection media:

	A	B
	Dulbecco's modified Eagle's medium (DMEM)	
	fetal bovine serum (FBS)	10%
	blasticidin S.	15 µg/ml

### Selection media:

	A	B
	Dulbecco's modified Eagle's medium (DMEM)	



	A	B
	fetal bovine serum (FBS)	10%
	blasticidin S	15 µg/ml
	hygromycin B	200 µg/ml

**Expression media:**

	A	B
	Dulbecco's modified Eagle's medium (DMEM)	
	fetal bovine serum (FBS)	10%
	blasticidin S	15 µg/ml
	hygromycin B	200 µg/ml
	tetracycline	0.1 µg/ml

## Cloning:

- Takara Catalog # 638947
- Gibco manual (Catalog # R780-07)

**Troubleshooting****Safety warnings**

- ! All appropriate biosafety precautions should be observed when handling Cell culture and recombinant DNA.



## Cloning











- 1 Clone gene of interest into pcDNA5/FRT/TO. For VPS13C<sup>m</sup>Clover the In-Fusion system (Takara Catalog # 638947) was used.
- 2 Sequence plasmid.

## Cell culture

- 3 Thaw and culture Flp-In<sup>TM</sup> T-REx<sup>TM</sup> 293 cells according to Gibco manual (Catalog # R780-07) in a 10cm dish. Once cells are 90% Confluent, split.
- 4 Freeze some cells according to Gibco manual protocol for future use.

## Transfection and selection of stably integrated cells

15m

- 5 Plate cells in 6 well format to ~30% confluence in media containing no antibiotics.
- 6 The next day, transfect cells using Fugene HD and a ratio of 9:1 pOG44: pcDNA5/FRT/TO-VPS13C<sup>m</sup>Clover.
- 7 Mix  2 µg total DNA (  0.2 µg pcDNA5/FRT/TO-VPS13C<sup>m</sup>Clover,  1.8 µg pOG44) in  100 µL opti-MEM. 
- 8 Add  8 µL Fugene HD to DNA mixture. Let sit  00:15:00 . 
- 9 Add Fugene/DNA mixture to cells in a dropwise fashion.
- 10 At 24h post transfection, wash cells and add fresh post-transfection media containing Blasticidin (no Zeocin or Hygromycin). 
- 11 At 48h post-transfection, split cells and plate at <25% confluence in selection medium containing Blasticidin and  200 µg/ml Hygromycin B (no zeocin).



12 Feed the cells with Hygromycin B selection medium every 3-4 days.

13 When foci become apparent, pool Foci and allow cells to expand.

#### Note

- In our experience, at  $200\text{ }\mu\text{g/ml}$  Hygromycin B, it may appear that no colonies remain viable, however after 1-2 weeks a small amount of colonies were observed in which the gene of interest was stably integrated.
  - If all cells are nonviable under these conditions, it may be necessary to reduce the concentration of Hygromycin B in the selection media to  $100\text{ }\mu\text{g/ml}$  and/or optimize transfection conditions.
- Per manufacturers protocol, these surviving foci should be isogenic for the inserted gene.

14 Pooled cells can then be split, frozen, and/or used for experiments.

15 Freeze cells according to Gibco manual.

## Expression of VPS13C

16 Plate cells at desired confluence in format appropriate for experiment in expression media containing  $0.1\text{ }\mu\text{g/ml}$  -  $1\text{ }\mu\text{g/ml}$  tetracycline.

#### Note

In our hands, we observed no difference in expression levels of VPS13C<sup>mClover</sup> protein between  $0.1\text{ }\mu\text{g/ml}$  and  $1\text{ }\mu\text{g/ml}$  and thus we used the lower concentration, however the Gibco manual suggests using  $1\text{ }\mu\text{g/ml}$ .

17 After 24 hours, cells may be lysed for western blot or used for microscopy.