

Feb 03, 2025

Generation of human iPSC-derived tricultures

 In 1 collection

DOI

<https://dx.doi.org/10.17504/protocols.io.rm7vzknz8vx1/v1>

Maria Jose Perez J.¹

¹Institut Imagine

Team Deleidi



Federico Bertoli

ASAP

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.rm7vzknz8vx1/v1>

Protocol Citation: Maria Jose Perez J. 2025. Generation of human iPSC-derived tricultures. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.rm7vzknz8vx1/v1>

License: This is an open access protocol distributed under the terms of the **[Creative Commons Attribution License](#)**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: November 06, 2024

Last Modified: February 03, 2025

Protocol Integer ID: 111674

Keywords: tricultures generation of human ipsc, derived tricultures generation, generation of human ipsc, human ipsc, generation

Funders Acknowledgements:

ASAP

Abstract

Generation of human iPSC-derived tricultures

Materials

Final neuronal maturation medium: 1:1 DMEM-Ham's F-12 and neurobasal medium, 0.5% N2, 1% B27 without vitamin A, 1% P/S, and 1% GlutaMAX, supplemented with 200 μ M AA, 1 μ M purmorphamine, and 100 ng/mL FGF8 (PeproTech) + with 20 ng/mL BDNF (PeproTech) and 1 mM dibutyryl-cyclic AMP (dcAMP, Sigma–Aldrich).



Generation of human iPSC-derived tricultures

- 1 Differentiate iPSCs into cortical neurons as per the established protocol ("Generation of human iPSC-derived cortical neurons, astrocytes, and microglia") until they reach 14 days in vitro (DIV).
- 2 At 14 DIV, dissociate neurons using Accutase for 5' at 37C
- 3 Plate dissociated neurons onto glass coverslips pre-coated with Matrigel 1:30 in DMEM 1X for 30' at 37C
- 4 Maintain plated neurons in neuronal maturation medium supplemented with 20 ng/mL BDNF and 1 mM dcAMP
- 5 After four days in culture, seed iPSC-derived astrocytes onto the neurons at a 3:1 ratio (neurons:astrocytes)
- 6 Maintain neuron-astrocyte cocultures in the final neuronal maturation medium
- 7 Three days after adding astrocytes, add iPSC-derived microglial precursors to the coculture at a 3:1:1 ratio (neurons:microglia:astrocytes)
- 8 Maintain the tricultures in the final neuronal maturation medium for an additional seven days in final neuronal maturation medium supplemented with 100 ng/mL IL-34, 100 ng/mL M-CSF, and 10 ng/mL GM-CSF