

Jan 11, 2026

Generation of SLC45A1 Knockout iPSC Line (WTC11-Ngn2 Line)

DOI

<https://dx.doi.org/10.17504/protocols.io.e6nvw45zzlmk/v1>

Marc Gastou¹, Ali Ghoochani^{2,3,4,5}, Natalia Gomez-Ospina¹, Monther Abu-Remaileh^{2,3,4,5}

¹Department of Pediatrics, Stanford University School of Medicine, Stanford, CA, 94305, USA;

²Department of Chemical Engineering, Stanford University, Stanford, CA 94305, USA;

³Department of Genetics, Stanford University, Stanford, CA 94305, USA;

⁴The Institute for Chemistry, Engineering and Medicine for Human Health (Sarafan ChEM-H), Stanford University, Stanford, CA 94305, USA;

⁵Aligning Science Across Parkinson's (ASAP) Collaborative Research Network, Chevy Chase, MD, 20815, USA

asap



Ali Ghoochani

Stanford University

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.e6nvw45zzlmk/v1>



Protocol Citation: Marc Gastou, Ali Ghoochani, Natalia Gomez-Ospina, Monther Abu-Remaileh 2026. Generation of SLC45A1 Knockout iPSC Line (WTC11-Ngn2 Line). **protocols.io** <https://dx.doi.org/10.17504/protocols.io.e6nvw45zzlmk/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: July 24, 2025

Last Modified: January 11, 2026

Protocol Integer ID: 223162

Keywords: generation of slc45a1 knockout, slc45a1 knockout, pluripotent stem cell, ipsc line, slc45a1, spcas9 system, ipsc, gene, crispr, other any gene, wtc11, ngn2 line, induced pluripotent stem cell, stem cell

Abstract

This protocol describes the generation of an **SLC45A1** knockout (or other any gene) in WTC11-Ngn2 induced pluripotent stem cells (iPSCs) by targeting the coding sequence (CDS) using the CRISPR-spCas9 system.

Materials

Reagents:

- Thiazovivin (ROCK1 inhibitor), 10 μ M
- ReLeSRTM (Stemcell Technologies, for cell dissociation)
- Cas9 enzyme (IDT Alt-R S.p. Cas9 Nuclease V3), 10 μ g/ μ L
- sgRNA (Synthego, 2'-O-methyl modified), 3.2 μ g/ μ L
- P3 Primary Cell NucleofectorTM Solution (Lonza, V4XP-3032)
- 16-well NucleocuvetteTM Strips
- 4D-NucleofectorTM System (Lonza)
- iPSC growth medium (e.g., StemFlex, Stemcell Technologies)
- Sterile PBS
- 24-well tissue culture plates



Designing the Single Guide RNA (sgRNA)

- 1
 - Clone 1 sgRNA: **AGATGCCCTCTTCCCCAGCG**
 - Clone 2 sgRNA: **ACTTCAAACGACACCCCAAG**

Pre-treatment of iPSCs (Day -1)

- 2 **Objective:** Enhance iPSC viability prior to nucleofection.
 1. Plate iPSCs at ~50% confluency in iPSC growth medium.
 2. Add 10 μ M Thiazovivin.
 3. Incubate overnight at 37 °C with 5% CO₂.

Preparation of Cas9-sgRNA RNP Complex

- 3 Assemble the Cas9-sgRNA ribonucleoprotein (RNP) complex.
 - 3.1 In a PCR tube, mix 6 μ g Cas9 protein with 3.2 μ g sgRNA.
 - 3.2 Incubate at 37 °C for 15 minutes.
 - 3.3 Keep the RNP complex on ice until use.

Dissociation of iPSCs

- 4 Prepare a single-cell suspension for nucleofection.
 - 4.1 Aspirate the medium and add ReLeSR™ to the iPSCs.
 - 4.2 Incubate at 37 °C for 3–5 minutes.
 - 4.3 Gently tap the plate to detach cells and resuspend them in iPSC growth medium.



4.4 Count viable cells and adjust to 0.3×10^6 cells per nucleofection.

4.5 Centrifuge at **300 × g for 3 minutes**.

4.6 Wash the pellet with sterile PBS and centrifuge again at **300 × g for 3 minutes**.

4.7

Resuspend the cell pellet in 20 µL P3 Nucleofector™ Solution per nucleofection (16.4 µl of P3 solution was mixed with 3.6 µl of Supplement).

Post-Transfection Treatment (Day 1–3)

5 Facilitate iPSC recovery and expansion.

5.1 After 24 hours, replace the medium with fresh iPSC growth medium.

5.2 Continue treatment with 10 µM Thiazovivin.

5.3 Culture cells for an additional 72 hours.

Validation of SLC45A1 Knockout

6 Amplify the targeted genomic region via PCR.

6.1 Perform Sanger sequencing of the PCR amplicons.

6.2 Analyze chromatograms using ICE software (Synthego) to determine indel frequency and knockout efficiency by comparing edited cells to wild-type or mock-treated controls.