

Jul 09, 2024

Genetic expression suppressor screen

 [Science Advances](#)

 In 1 collection

DOI

<https://dx.doi.org/10.17504/protocols.io.ewov19z67lr2/v1>

Alexandros C Kokotos^{1,2}, Tim Ryan^{1,2}

¹Department of Biochemistry, Weill Cornell Medicine, New York, NY 10065, USA;

²Aligning Science Across Parkinson's (ASAP) Collaborative Research Network, Chevy Chase, Maryland 20815, USA



Alexandros C Kokotos

Weill Cornell Medicine

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.ewov19z67lr2/v1>

Protocol Citation: Alexandros C Kokotos, Tim Ryan 2024. Genetic expression suppressor screen. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.ewov19z67lr2/v1>

Manuscript citation:

Phosphoglycerate kinase is a central leverage point in Parkinson's Disease driven neuronal metabolic deficits

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: July 08, 2024

Last Modified: July 09, 2024

Protocol Integer ID: 103021

Keywords: ASAPCRN, Live-cell imaging, neurons, low-throughput screening, glycolytic enzyme, genetic expression suppressor screen, genetic expression suppressor screen this protocol, genetic screen, throughput genetic screen, glycolytic enzyme, neuronal function under hypometabolic condition, neuronal function

Funders Acknowledgements:

ASAP

Grant ID: 000580

Abstract

This protocol describes a low-throughput genetic screen to identify glycolytic enzymes that improve neuronal function under hypometabolic conditions.



Materials

Media Supplies

- 6-cyano-7-nitroquinoxaline-2,3-dione (CNQX; Alomone Labs, Cat. No. C-141)
- Calcium Chloride (Millipore-Sigma, Cat. No. C3881)
- DL-2-Amino-5-phosphonovaleric acid (APV; Millipore-Sigma A5282)
- Glucose (Millipore-Sigma, Cat. No. G7021)
- HEPES (Millipore-Sigma, Cat. No. H3375)
- Magnesium Chloride (Millipore-Sigma, Cat. No. M0250)
- Potassium Chloride (Millipore-Sigma, Cat. No. P9333)
- Sodium Chloride (Millipore-Sigma, Cat. No. S5761)

Stock solutions

APV (25mM Stock)

1. Add 20L ddH₂O to whole bottle of APV
2. Sonicate at ~55C for 30 minutes.
3. Make 1mL aliquots and store at -20C.

CNQX (10mM Stock)

1. Add 10 mL of ultrapure H₂O to the whole bottle of CNQX.
2. Vortex until dissolved.
3. Make 1 mL aliquots and store at -20 oC.

Media solutions

0.1mM glucose tyrodes (119mM NaCl, 2.5mM KCl, 0.1mM D-Glucose, 54.9mM HEPES, 2mM CaCl₂, 2mM MgCl₂, 50uM APV, 10uM CNQX)

NH₄Cl tyrodes (119mM NaCl, 2.5mM KCl, 5mM HEPES, 50mM NH₄Cl, 2mM CaCl₂, 2mM MgCl₂, 50uM APV, 10uM CNQX)

Equipment

- 40x, 1.3 N.A Objective (like Zeiss-Fluar)
- 488nm Laser (like Coherent-OBIS 488nm LX Laser)
- CO₂ Incubator (like VWR Symphony)
- EMCCD Camera (like Andor iXon Ultra 897)
- Epifluorescence Microscope (like Zeiss ObserverA1)
- Filter cube set for GFP excitation and emission (like Chroma filters)
- Laser speckle reducer (like Optotune LSR 3005)
- Laser/stage/stimulation-Master board controller (like a custom-designed **Arduino board**)
- Light transmission microscope (like Olympus CK40)
- Mercury arc lamp (like Zeiss HBO 100)
- Perfusion system (like Automate Scientific Valvelink 8.2)
- Microscope XYZ Stage (like Applied Scientific Instrumentation stage)



- Stimulus isolator (like World precision Instruments A385)
- Temperature controller for Objective (like Minco)

Plasmids

- vGLUT1-pHluorin
- mTagBFP2-N1
- hSyn HK1 (Addgene plasmid #220915; <http://n2t.net/addgene:220915>; RRID:Addgene_220915)
- hSyn AldoA (Addgene plasmid #220916; <http://n2t.net/addgene:220916>; RRID:Addgene_220916)
- hSyn AldoC (Addgene plasmid #220917; <http://n2t.net/addgene:220917>; RRID:Addgene_220917)
- hSyn GAPDH (Addgene plasmid #220918; <http://n2t.net/addgene:220918>; RRID:Addgene_220918)
- hSyn PGK1 (Addgene plasmid #220919; <http://n2t.net/addgene:220919>; RRID:Addgene_220919)
- hSyn PKM1 (Addgene plasmid #220920; <http://n2t.net/addgene:220920>; RRID:Addgene_220920)



Neuronal culture and transfection

- 1 Use cultured rodent primary hippocampal neurons, transfected with vGlut1-pH, mTagBFP2 and a construct expressing a single glycolytic enzyme, as explained in dx.doi.org/10.17504/protocols.io.ewov1qxr2gr2/v1.

Live-cell imaging

25m

- 2 Remove the cells from the incubator and bathe them in Tyrodes containing 0.1 mM glucose.
- 3 Mount the coverslips with the cells on the imaging chamber and image as explained in dx.doi.org/10.17504/protocols.io.q26g7pn4qgwz/v1.
- 4 Locate cells in an unbiased manner using the mTagBFP2 channel and center on their axon, judging by their morphology.
- 5 After a total of  00:05:00 that the cells have been in the 0.1 mM glucose tyrode buffer, start the live-cell imaging experiment in the vGlut1-pH channel. 5m
- 6 Stimulate the cells with 100 APs delivered at 10 Hz every  00:01:00 , for a total of 10 rounds of activity. 1m
- 7 At the end of the experiment, perfuse cells with Tyrodes containing 50 mM NH₄Cl.
- 8 Analyze the collected data accordingly.

Protocol references

dx.doi.org/10.17504/protocols.io.ewov1qxr2gr2/v1

dx.doi.org/10.17504/protocols.io.q26g7pn4qgwz/v1