

Dec 03, 2024

Genome Edited (KO) Cell Line Generation –

 [The Journal of Cell Biology](#)

DOI

<https://dx.doi.org/10.17504/protocols.io.dm6gp9bp5vzp/v1>

Amanda Bentley-DeSousa^{1,2,3,4,5}, Shawn Ferguson^{1,2,3,4,5,6}

¹Department of Cell Biology, Yale University School of Medicine, New Haven, Connecticut 06510, USA;

²Neuroscience, Yale University School of Medicine, New Haven, Connecticut 06510, USA;

³Program in Cellular Neuroscience, Neurodegeneration and Repair;

⁴Wu Tsai Institute Yale University School of Medicine, New Haven, Connecticut 06510, USA;

⁵Aligning Science Across Parkinson's (ASAP) Collaborative Research Network, Chevy Chase, MD, 20815, USA;

⁶Kavli Institute for Neuroscience, Yale University School of Medicine, New Haven, Connecticut 06510, USA



Amanda Bentley-DeSousa

Yale University

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.dm6gp9bp5vzp/v1>

Protocol Citation: Amanda Bentley-DeSousa, Shawn Ferguson 2024. Genome Edited (KO) Cell Line Generation –. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.dm6gp9bp5vzp/v1>

**Manuscript citation:**

<https://doi.org/10.1101/2023.10.31.564602>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: December 02, 2024

Last Modified: December 03, 2024

Protocol Integer ID: 113413

Keywords: RAW 264.7, KO, Synthego, ko cell lines in raw, ko cell line, cell line generation, cell lines in raw, genome edited, genome, cell, line, ko, raw

Funders Acknowledgements:

Aligning Science Across Parkinson's Disease

Grant ID: ASAP-000580

Abstract

This protocol details steps taken to create KO cell lines in RAW 264.7 cells.

Materials

DMEM (Thermo Fisher Scientific, 11965-092)

FBS (Thermo Fisher Scientific, 16140-071)

Penicillin/Streptomycin (Thermo Fisher Scientific, 15140122)

CellStripper (Corning, 356230)

Synthego CRISPR Gene Knockout V2 Mouse Kits

Lipofectamine CRISPRiMAX (ThermoFisher Scientific, CMAX00003)

Cas9 (Synthego CRISPR Gene Knockout V2)

Opti-Mem (Thermo Fisher Scientific, 31985062)



Day 0

- 1 Plate 2.5×10^5 RAW 264.7 cells per well in a 6-well dish.

Day 1

- 2 Transfect cells using Lipofectamine CRISPRiMAX, gene specific sgRNAs and recombinant Cas9 (Synthego CRISPR Gene Knockout V2).

2.1 CRISPRiMAX protocol

In Tube 1 combine  62.5 μL Opti-MEM,  3.25 μL sgRNA (from 3 uM stock),

 2.5 μL Cas9 (from 3 uM stock), and  2.5 μL Cas9+ then incubate for

 00:05:00 . In Tube 2 combine  62.5 μL Opti-MEM and  3.75 μL

CRISPRMAX then incubate for  00:05:00 . Combine tube 1 and tube 2. Incubate for

 00:20:00 . Add drop-wise to 1 well containing  1.8 mL of media.

Day 3

- 3 Change the media.

Day 4

- 4 Plate single cells into 96-well dishes to obtain clonal cell lines. After subculturing, confirm cell lines via immunoblotting