

Sep 17, 2019 Version 1

Genomic DNA extraction from diatom *P. multistriata* V.1

DOI

dx.doi.org/10.17504/protocols.io.7dghi3w

Francesco Manfellotto¹, Monia Teresa Russo¹, Rossella Annunziata², Mariella Ferrante²

¹SZN; ²Stazione Zoologica Anton Dohrn



Francesco Manfellotto

Stazione Zoologica Anton Dohrn

OPEN  ACCESS



DOI: dx.doi.org/10.17504/protocols.io.7dghi3w

Protocol Citation: Francesco Manfellotto, Monia Teresa Russo, Rossella Annunziata, Mariella Ferrante 2019. Genomic DNA extraction from diatom *P. multistriata*. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.7dghi3w>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: September 16, 2019

Last Modified: September 17, 2019

Protocol Integer ID: 27784

Abstract

Genomic DNA extraction from diatom *P. multistriata*



- 1 Grow cells in 250 ml
- 2 filter coltures with 1.2 μ m nitrocellulose membranes
- 3 Centrifuge at 6000 rpm for 5 minutes at 4 °C and remove medium
- 4 Resuspend cells with 500 μ l of TE buffer (10 mM TrisHCl pH 7.6 and 1 mM EDTA pH 8.0)
- 5 Add:
 - 400 mg of 0.2-0.3 mm zirconia/silica diameter beads,
 - 500 μ l phenol.
- 6 rinse filter with 1 ml f/2 medium in a falcon and then move cells to a 2 ml eppendorf
- 7 Mix with vortex 30 Hz 3 times for 85 seconds, each time put sample in ice for 60 seconds before vortex.
- 8 · Centrifuge at 10000 rpm for 5 minutes at 4°C.
- 9 Recover aqueous phase in new 1.5 ml eppendorf (about 600 μ l)
- 10 Add 500 μ l of PCI (Phenol:Chloroform: isoamyl alcohol 25:24:1 v/v) and mix by inversion.
- 11 Centrifuge at 10000 rpm for 5 minutes at 4°C
- 12 Move the aqueous phase in a new eppendorf and add 5 μ l of RNase-A 10 mg/ml
- 13 Incubate at 37 °C for 30 minutes.

3m



- 14 Add 500 µl di PCI (Phenol:Chloroform: isoamyl alcohol 25:24:1 v/v) mix by inversion.
- 15 Centrifuge at 10000 rpm for 5 minutes at 4°C.
- 16 Move the aqueous phase in a new 2 ml eppendorf and add:
 - 50 µl of 3 M NaAc (pH ± 5)
 - 1 ml of ethanol 96% (- 20 °C)
 - 2 µl glycogen (- 20 °C)
- 17 incubate over night at -20°C. 1h
- 18 Centrifuge the overnight samples at 13000 rpm for 30 minutes at 4°C 3m
- 19 discard aqueous phase
- 20 Wash the pellet by adding 1 ml ethanol 70% and mix gently by inversion
- 21 Centrifuge at 13000 rpm for 10 minutes at 4°C 1m
- 22 discard aqueous phase
- 23 Wash the pellet by adding 1 ml ethanol 70% and mix gently by inversion
- 24 Centrifuge at 13000 rpm for 10 minutes at 4°C 1m
- 25 discard aqueous phase
- 26 Remove aqueous phase and dry pellet at RT for at least 20 minutes 2m



27 add 50 μ l of Preheated TE 1X (pH 8) or sterile MilliQ water to pellet of DNA

Preheat the TE pH 8 or sterile MilliQ water at 55 °C

28 incubate at 55 °C for 20 minutes

2m

29 quantify DNA concentration by nanodrop or Qubit

30 in order to check DNA integrity, run a small amount of DNA with 1% agarose gel

31 The DNA is ready and store it at -20°C