

Jan 18, 2026

Genomic engineering of iPSCs

DOI

<https://dx.doi.org/10.17504/protocols.io.yxmvmdbdzb3p/v1>

Ellen Hertz¹, Yu Chen¹, Joel Blanchard², Ellen Sidransky¹

¹National Institutes of Health; ²Icahn School of Medicine at Mount Sinai



Ellen Hertz

National Institutes of Health

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.yxmvmdbdzb3p/v1>

Protocol Citation: Ellen Hertz, Yu Chen, Joel Blanchard, Ellen Sidransky 2026. Genomic engineering of iPSCs. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.yxmvmdbdzb3p/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: December 11, 2025

Last Modified: January 18, 2026

Protocol Integer ID: 234760

Keywords: ASAPCRN, genomic engineering of ipsc, genomic engineering, ipsc

Funders Acknowledgements:

ASAP

Grant ID: ASAP-000458

Abstract

Genomic engineering of iPSCs as in Hertz et al 2025



Prepare cells for transfection

- 1 Grow iPSCs in Essential 8 Medium (Cat.# A1517001 Thermo Fisher) until ca 70% confluency
- 2 Wash cells with PBS
- 3 Dissociate iPSC into single cells with 1 ml StemPro-Accutase (Cat.# A1110501; Gibco)
- 4 Incubate at 37 C for 5 min
- 5 Resuspend cells in DMEM and spin down at 300 x g for 3 min
- 6 Resuspend cells in Essential 8 Medium with Y-27632 and seed at a density of 200,000 cells/cm²
- 7 Let cells settle for 2 hours at 37 degrees

Transfection

- 8 Mix Lipofectamine Stem Transfection Reagent (Cat.# STEM00001; Invitrogen) and plasmids in OPTI-MEM (Cat # 31985062, Thermo) according to instructions
- 9 Incubate for 5 min at RT
- 10 Add drop wise to the cells
- 11 Incubate over night at 37 C and 5% CO

Selection of transfected cells



- 12 Reseed cells (as described above) and resuspend cells in StemFlex (Cat.# A3349401; Thermo Scientific), supplemented with Y-27632
- 13 Incubate 24 h
- 14 Start selection with 0.25 mg/ml puromycin (Cat.# A11138-03, Thermo Scientific)
- 15 Change media (including puromycin) every 2-3 days until the selected clones are almost confluent (days-weeks)

CRE-mediated excision of the selection cassette

- 16 Transfection of plasmid encoding GFP-tagged CRE-recombinase to remove the selection cassette as above. Incubate 24 h in cell incubator
- 17 Coat 96 well plate with laminin-521 overnight in incubator
- 18 FACS sort individual clones into laminin coated plate with Stemflex supplemented with 10 μ M Y-27632, 5 μ M Emricasan (Cat.# 7310, R&D), 1x penicillin-streptomycin, 0.7 μ M ISRIB (Cat.# 5284, R&D) and 1x polyamines (Cat.# P8483, Sigma)
- 19 Incubate cells 3 days
- 20 Change 1/2 media with Stemflex supplemented with penicillin-streptomycin, ISRIB and polyamines
- 21 Incubate 3 days
- 22 Change 1/2 media with Stemflex every 3 days until clones reach ca 50% confluency.
- 23 Expand clones