

Apr 27, 2022

Genomic RNA extraction by TRIzol method by using *Drosophila melanogaster* gut

DOI

dx.doi.org/10.17504/protocols.io.dm6gpb92plzp/v1



Mohamed Sabeelil Islam¹

¹Kuvempu University

Drosophila lab



Mohamed Sabeelil Islam

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.dm6gpb92plzp/v1>

Protocol Citation: Mohamed Sabeelil Islam 2022. Genomic RNA extraction by TRIzol method by using *Drosophila melanogaster* gut. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.dm6gpb92plzp/v1>

License: This is an open access protocol distributed under the terms of the **[Creative Commons Attribution License](#)**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited



Protocol status: Working

We use this protocol and it's working

Created: April 04, 2022

Last Modified: April 27, 2022

Protocol Integer ID: 60261

Keywords: genomic rna extraction by trizol method, genomic rna extraction, rna extraction, using drosophila melanogaster, drosophila melanogaster, rna, gut simple method, trizol method, extraction, genomic, gut

Disclaimer

DISCLAIMER – FOR INFORMATIONAL PURPOSES ONLY; USE AT YOUR OWN RISK

The protocol content here is for informational purposes only and does not constitute legal, medical, clinical, or safety advice, or otherwise; content added to protocols.io is not peer reviewed and may not have undergone a formal approval of any kind. Information presented in this protocol should not substitute for independent professional judgment, advice, diagnosis, or treatment. Any action you take or refrain from taking using or relying upon the information presented here is strictly at your own risk. You agree that neither the Company nor any of the authors, contributors, administrators, or anyone else associated with protocols.io, can be held responsible for your use of the information contained in or linked to this protocol or any of our Sites/Apps and Services.

Abstract

Simple method followed for RNA extraction of Drosophila melanogaster larval(3rd instar) gut for CDNA

Guidelines

ware gloves while extraction

Troubleshooting

- 1 Dissection was used to isolate the larval gut of *Drosophila* 3rd Instar.
(around 10-15 gut was taken)
- 2
The dissected gut was taken in a 1.5ml ependorf tube which contained 200µL of TRIzol.

(For a total of 800µL or 0.8 mL of TRIzol needed for the extraction)
To avoid spurt, 200µL were taken.
- 3 Homogenize by using a micropistle
- 4 Add 600µL or 0.6mL of TRIzol followed up by adding 200µL or 0.2mL of Chloroform
(after adding vortex it) and keep it for 5 minutes incubation
- 5 Centrifuge at 12000x for 10minutes at 4°C
- 6 Take supernatant into new tube and add equal amount of isopropanol
- 7 Centrifuge at 12000x for 10minutes at 4°C
- 8 discard the supernatant and retain the pellet
- 8.1 wash the pellet by using 1ml of 75% cold ethanol (twice)
(after adding ethanol vortex the solution and centrifuge at 7500x for 5 minutes)
- 9 Retain the pellet, add RNase free water and check the purity