



Jan 02, 2025

Genotyping *C. elegans*

 In 1 collection

DOI

<https://dx.doi.org/10.17504/protocols.io.yxmvmebbbg3p/v1>

Justin Donnelly¹

¹Arcadia Science

Arcadia Science



Arcadia Science

Arcadia Science

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.yxmvmebbbg3p/v1>

Protocol Citation: Justin Donnelly 2025. Genotyping *C. elegans*. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.yxmvmebbbg3p/v1>

Manuscript citation:

Avasthi P, Borges AL, Cheveralls K, Donnelly J, Mets DG, Reiter T. (2025). An experimental and computational workflow to characterize nematode motility behavior. <https://doi.org/10.57844/arcadia-b89a-7e76>



License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: August 23, 2024

Last Modified: March 13, 2025

Protocol Integer ID: 106384

Keywords: genomic extraction, pcr amplification, strain, pcr

Abstract

This protocol describes a procedure to genotype *C. elegans* strains via genomic extraction and PCR amplification. You should always confirm strain information upon receipt of the strain.

Materials

☒ Proteinase K, Molecular Biology Grade - 2 ml **New England Biolabs Catalog #P8107S**

☒ Magnesium Chloride **Fisher Scientific Catalog #AC223210010**

☒ Potassium chloride **P212121**

☒ Tris HCl **P212121**

PCR tubes

Thermocycler

Platinum wire

PCR master mix

PCR primers

1% agarose DNA gels

☒ New 6x Purple Loading Dye **New England Biolabs Catalog #B7024S**

Ultrapure water

Protocol materials

☒ New 6x Purple Loading Dye **New England Biolabs Catalog #B7024S**

☒ Proteinase K, Molecular Biology Grade - 2 ml **New England Biolabs Catalog #P8107S**

☒ Magnesium Chloride **Fisher Scientific Catalog #AC223210010**

☒ Potassium chloride **P212121**

☒ Tris HCl **P212121**



Genomic extraction from single worms

15s

- 1 Prepare lysis buffer. You can prepare lysis buffer in advance and store it at Room temperature before adding proteinase K.

Lysis buffer (without proteinase K):

[M] 10 millimolar (mM)

Tris HCl **P212121**

[M] 50 millimolar (mM)

Potassium chloride **P212121**

[M] 2 millimolar (mM)

Magnesium Chloride **Fisher Scientific Catalog #AC223210010**

8.3

- 2 Complete lysis buffer by adding

Proteinase K, Molecular Biology Grade - 2 ml **New England Biolabs Catalog #P8107S**

to a final concentration of [M] 0.5 mg/mL . Prepare $\geq$ 10 μ L complete lysis buffer per sample.

- 3 Transfer 10 μ L lysis buffer per tube to PCR tubes. Include an additional tube as a negative control.

- 4 **Pick a single worm from each strain and transfer to lysis buffer**

Flame-sterilize a platinum wire. Pick up one adult worm from plate, taking care to prevent puncturing the NGM. Transfer worm by gently submerging the platinum wire into the PCR tube under flame. Hold the wire in this position for several seconds to allow worms to detach. Flame-sterilize the wire before picking additional worms.

Note

Try to pick only a single worm per PCR tube. A couple of worms is fine, but too many can inhibit the reaction.

- 5 Immediately spin tubes in a benchtop microcentrifuge for 00:00:15 .

15s

- 6 Freeze tubes in LN2 or -80 $^{\circ}$ C for $\geq$ 00:10:00 .

10m



- 7 Lyse worms in a thermocycler according to the following protocol:

2h 45m

🔥 65 °C , ⌚ 01:00:00 – ⌚ 01:30:00 .

🔥 95 °C , ⌚ 00:15:00 (proteinase K inactivation step).

You can store extracted worm DNA indefinitely at 🔥 -80 °C .

PCR amplification

- 8 If using poison primer amplification, add 🧴 46 µL Q5 PCR master mix per tube to appropriate number of PCR tubes. To each, add 🧴 2 µL template (proteinase K extraction product) and 🧴 0.5 µL each primer.

If using standard amplification, add 🧴 47 µL Q5 PCR master mix per tube to appropriate number of PCR tubes. To each, add 🧴 2 µL template (proteinase K extraction product) and 🧴 0.5 µL each primer.

For negative control, add master mix and primers only.

- 9 Run PCR for 30–35 cycles according to Q5 manufacturer protocol. Determine annealing temperature based on primer sequences.

Analysis of PCR product via gel electrophoresis

- 10 Transfer 🧴 20 µL each PCR product to a new PCR tube and add 🧴 4 µL  New 6x Purple Loading Dye **New England Biolabs Catalog #B7024S** .
- 11 Load onto 1% agarose DNA gel. Electrophorese at 🧴 120 V until dye front reaches the end of the gel.
- 12 Visualize under blue light.



Note

For large-scale deletion strains with validated poison primer sequences, gel analysis is usually sufficient to confirm genotype. In other cases, further analysis via Sanger sequencing is required.

Analysis of PCR product via sequencing

- 13 In a clean PCR tube, add  9 μL ultrapure water,  1 μL PCR product, and  1 μL sequencing primer.
- 14 Submit sample for Sanger sequencing.

Protocol references

Protocol modified from:

<https://bpb-us-w2.wpmucdn.com/sites.wustl.edu/dist/0/1532/files/2018/07/Single-Worm-PCR-r8hwch.pdf>